

Qualification Documentation Alarm and Function Testing

Glaxowellcome Production France

CUC 2002

Cartoning machine

Serial Number	Overall order no.
100223	3100001673

© ROTZINGER PharmaPack GmbH

All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.

Table of contents 1/1

Title of Index Chapter Content **CUC 2002** Index no. **Cartoning machine**

	Alarm and Function Testing		
Alarm and	1	References	1
Function Testing Plan	2	Introduction	
	3	Scope	
	4	Test philosophy	
	5	Test scope	
	6	Glossary	
	7	Appendix	
Alarm and	1	Software version	2
Function Testing Execution	2	Check HMI menu, buttons and switches	
	3	Access protection	
	4	Audit trail	
	5	Version management for recipe	
	6	Batch Management	
	7	Reports and print outs	
	8	Basic function	
	9	Alarm Function Tests	
	10	Notes	
Alarm and Function Testing Report			3

Alarm and Function Testing Plan

CUC 2002 Cartoning machine

100223

Glaxowellcome Production France

Creation of this document supported by:



Valicare GmbH a Syntegon company Eschborner Landstraße 130-132

60489 Frankfurt a. Main Germany

www.valicare.com

Alarm and Function Testing Plan approval by Valicare GmbH and ROTZINGER PharmaPack GmbH (before execution)

Function	Name	Position/Company	Signature / Date
Created by	Lorena Fausten	Senior Validation Expert/ Valicare GmbH	28. Mar - 7073
Reviewed and approved by	Andreas Harzer	Validation Project Manager / Valicare GmbH	i. A. Marroi
Reviewed and approved	Florian Passow	Head of Documentation and Validation /	i.A.
by		ROTZINGER PharmaPack GmbH	28. Mar. 2023

Alarm and Function Testing Plan approval by Glaxowellcome Production (before execution)

Function	Name	Position/Company	Signature / Date
Approved by			
Approved by			
Approved by			

Document history:

Rev.	Date	Description	Author
1.0	05.Dec.2022	Version for approval	Lorena Fausten
2.0	19.Jan.2023	Version for approval Rev. 2.0, hardware and software changes after FAT 1	Lorena Fausten
3.0	17.Mar.2023	Version for approval Rev. 3.0, hardware and software changes after FAT 2	Lorena Fausten
4.0	28.Mar.2023	Version for approval Rev. 4.0, software changes after FAT 3	Lorena Fausten

Table of contents

1	References	4
2	Introduction	4
3	Scope	4
4	Test philosophy	5
4.1	Test approach	5
4.2	Test structure	5
4.2.1	Access protection	5
4.2.2	Version management	6
4.2.3	Audit trail	6
4.3	Test procedure	7
4.4	Test prerequisites	7
4.5	Data to be recorded	7
4.6	Test result	7
4.7	Handwritten entries and comments	7
4.8	Deviation handling	8
4.9	Change control	8
4.10	Test summary	8
4.11	Personnel involved in testing	8
5	Test scope	9
6	Glossary	9
7	Appendix	10
7.1	Identification of all Personnel involved	10
7.2	Test scope	11
7.3	Deviation sheet	19
7.4	Change control sheet	20
7.5	Forms	21

1 References

No.	Title
[1]	Order confirmation
[2]	SOP-AG-06.013e: Qualification deviation procedure
[3]	SOP-AG-23.007e: Handwritten documentation for qualification
[4]	SOP-AG-24.006e: Procedure for change control
[5]	Alarm and Function Testing Report

2 Introduction

This document forms the Alarm and Function Testing of the ROTZINGER PharmaPack GmbH equipment installed at Glaxowellcome Production, France. The execution of this document shows that the product administration as well as faults, messages, warnings and software switches function correctly.

3 Scope

This document contains the Alarm and Function Testing Plan and the references necessary to collect and record the evidence that shows the CUC 2002, 100223 operates correctly.

The test results will be approved after completion of the Alarm and Function Testing Report.

4 Test philosophy

The following chapter describes the general test philosophy for this Alarm and Function Testing.

4.1 Test approach

The objective of this Alarm and Function Testing is to gather and document evidence that verifies that the ROTZINGER PharmaPack GmbH equipment operates according to the manufacturer's requirements and updated design documents.

The Alarm and Function Testing plan has been produced according to the agreed scope of supply in [1], which is based on current GMP requirements.

It will be approved prior to execution by ROTZINGER PharmaPack GmbH and Glaxowellcome Production representatives. After testing, and closing of the nonconformances, the Alarm and Function Testing report can be issued and approved.

All information for an individual test as well as any additional documentation used or collected during the Alarm and Function Testing execution, will be attached to the individual test protocols and annotated with the test number, date and initials (see [2]).

The agreed acceptance criteria for each test are defined in the test scope (see chapter 7.2).

4.2 Test structure

All qualification tests of this document have the following structure:

- Test title
- Test number
- Test objective
- Test prerequisites
- Test procedure
- Acceptance criteria
- Data to be recorded (if applicable)
- Test result(s)
- Comments
- Final test assessment

4.2.1 Access protection

The machine software provides an access protection feature.

The access protection is connected to the customer network. In this special case, the user administration is not locally set, but as a domain, i.e. the Windows-domain, which is overlapping. This includes a multiple connection of different stations over one network.

An advantage of this administration system is that a new user can be added in the domain and all computers in the network dispose of this new user. User administration is simplified thereby.

User administration - domain

In the domain the user administration contains different levels:

- Groups: all users with the same characteristic parameters are summarized in a group and certain rights are assigned by the Administrator.
- User: A user is defined over its access data (user name, password)

Comment:

The configuration of user data is defined by the customers system and can not be set absolutely. As a reference the FDA 21 CFR Part 11/GMP Vol.4 Annex 11 provisions have to be maintained.

- The combination of user name and password is unique and allocated definitely to one access authorization.
 - (FDA 21 CFR Part 11 §11.300(a)/GMP Vol.4 Annex 11: 12.1)
- Only the Administrator can create or disable, but not delete, a user name. (FDA 21 CFR Part 11 §11.300(b)/GMP Vol.4 Annex 11: 12.1)
- The password becomes invalid automatically after x days, the user receives pertinent information y days previous to this date - x,y parameters can be set by the Administrator. (FDA 21 CFR Part 11 §11.300(b)/GMP Vol.4 Annex 11: 12.1)
- Further settings, concerning structure and length of the password and the definition of forbidden passwords ("No-go-List" for passwords) are defined by the Administrator.
- Consequences referring to invalid inputs and the corresponding inhibition of the system by Administrator

(FDA 21 CFR Part 11 §11.300(d)/GMP Vol.4 Annex 11: 12.3)

4.2.2 Version management

Version management - format - product - recipe data

- When creating a new format, a count up "number" is generated (e.g. "1"; "2"; "3"), the format name can be renamed, the text is open for configuration (e.g. "Aspirin 50 mg)
- One format can include different drafts or released versions with corresponding machine setting parameters. For a newly created draft version, the (version) number is increased by 1. When the draft is released to "released version" it changes its status from draft to released. (e.g. Format 1 and released version 4 will be named version 1.004)
- Preset format parameters can only be changed in a draft version. The draft version can be released, a released version can be activated for production.
- Versions can be disabled by a user (definition per matrix), however only the Admin can delete versions. All existing versions are shown in the HMI.

4.2.3 Audit trail

- All user actions related to process, quality and product as well as parameter changes must be recorded in the audit trail: e.g.:
 - Logon/logout of users, also invalid access attempts and system disable (FDA 21 CFR Part 11 §11.300(d)/GMP Vol.4 Annex 11: 12.1, 12.3)
 - Modifications of the access matrix (FDA 21 CFR Part 11 §11.10(d)/GMP Vol.4 Annex 11: 12.1, 12.3)
 - Faults, errors, message display
 - Creating, changing format/recipe data (filling/empty weights, ...)
 - Change of defaults of checking devices (code reader, camera, ...)
 - Checking devices On/Off
 - Options On/Off
 - Format change
 - Modifications of sizes / speeds / reject limits/ weights, ... relevant to product (machine speed, r.p.m, shift register ...)
 - Change of operating mode (Line, Automatic, Setup, ...)
- The audit trail is kept batch-related. A filter can be used to fade out events outside the batch (before and during interruptions of the batch).
- Example for batch-related: All events before starting the batch, during a batch interruption and during the batch are written to the audit trail, which is concluded at the end of the batch.
- If no batch management has been realized, the audit trails are managed by a different logical system, such as e.g. by days.
- Audit trails can be neither modified nor deleted.

- Audit trails are stored with a security system, any later entries/modifications are detected by the system (FDA 21 CFR Part 11 §11.10(c)/GMP Vol.4 Annex 11: 7.1, 12.1)
- All audit trail entries show the action, the complete user name, date and time (FDA 21 CFR Part 11 §11.10(e)/GMP Vol.4 Annex 11: 12.1, 12.3, 12.4)
- The audit trails can be according to criteria such as user, activity, date and time filtered and printed online.
- The user can enter a free comment to the log at any time, e.g. via the superimposed keyboard.
- The audit trails are part of the production report.

Example of an audit trail:

Time	User	Text
08/20/2020 3:11:44 pm	Maier	User logged on (name: group_xyyxx)
08/20/2020 3:12:44 pm	Maier	switch <30> Data from <off> to <on></on></off>
08/20/2020 3:15:44 pm	Maier	User logged out (name: group_xyyxx)

4.3 Test procedure

The procedure for executing each test is defined in the test scope table (see chapter 7.2) and/or in the test protocol.

4.4 Test prerequisites

If applicable the prerequisites for each test are listed in the test protocol.

4.5 Data to be recorded

A data specification that needs to be retained as documented evidence, will be stored in the relevant test procedure.

4.6 Test result

During test execution, individual test results have to be compared with the expected results, and an assessment has to be made for each acceptance criteria separately whether the result complies (yes) or does not comply (no).

A final assessment will close a test. "Results comply: Yes" (Passed) will be used if all acceptance criteria have been met. In case one or more acceptance criteria can not be met, the test will be assessed as failed ("Results comply: No").

The final assessment of a test (or retest) will then be recorded in the test scope (see chapter 7.2).

4.7 Handwritten entries and comments

Data entries,	comments,	corrections	or signatures	manually wri	tten onto the	prepared test	protocol wil	ll be
performed in	accordance	with [3] into	the rounded	boxes of the	document as	the test is ex	ecuted.	

Example:	
Example.	

Any correction to handwritten data will be made by the person entering the data. Each page used for executing the qualification is to be signed and dated.

Data that is gathered by instrumentation onto a recorder will be printed and attached to a test protocol. Applicable data calculations will be transferred to a computer and the calculation will be printed and attached to a test protocol. All attachments will be signed by the responsible author.

4.8 Deviation handling

All exceptions that are discovered during execution are to be entered into the deviations sheet according to [2] using the deviation sheet in chapter 7.3.

The deviation sheets are collected in the Alarm and Function Testing deviation chapter of the Alarm and Function Testing folder.

4.9 Change control

Changes during the qualification activities of this document shall be treated according to GDP.

4.10 Test summary

Status of the Alarm and Function Testing execution, and follow up actions after completion of the Alarm and Function Testing, are shown in [5].

4.11 Personnel involved in testing

All personnel involved in the qualification activities of this document will be listed using the table in chapter 7.1.

5 Test scope

A list of tests to be performed including the test procedures and acceptance criteria to be used during the Alarm and Function Testing can be found in the table in chapter 7.2: Test scope.

Following execution, the completed table will contain the final result for each test. It will then be placed in the "test scope chapter".

6 Glossary

AFT	Alarm and Function Testing	
CFR	Code of Federal Regulations	
FDA	Food and Drug Administration	
FLT	Fault (message)	
GDP	Good Documentation Practice	
GMP	Good Manufacturing Practice	
ME	Message	
PLC	Programmable Logic Controller	
QD	Qualification Documentation	
SOP	Standard Operating Procedure	
sws	Software Switch	
SWSOPM	Software Switch for Operation Mode	
URS	User Requirement Specification	
WA	Warning (message)	
n.a.	not applicable	

7 Appendix

7.1 Identification of all Personnel involved

All personnel involved in the qualification activities of this document will be listed in the table below and identified in the relevant test protocols.

The involved personnel must be authorized and familiar with the machine operation.

Name	Position/Company	Date	Signature	Initials

7.2 Test scope

A list of all tests to be performed during the Alarm and Function Testing (AFT), including information on the execution of the tests and acceptance criteria, are listed in the following table. After the execution of the tests, this table also contains all final results for each test.

The test protocols for the executed tests listed in the following table are filed as an appendix to the Alarm and Function Testing Plan in index 2 as "Alarm and Function Testing Execution".

Cartoning machine CUC 2002

No.	Test title	•	Deviation	Results comply		
			(No.)	Test (yes/no)	Retest (yes/no)	Date/Signature
1	Software version	The software version is documented.				
2	Check HMI menu, but	ttons and switches				
2.1	Spotcheck HMI	HMI menu and page tab spotcheck is OK				
3	Access protection					
3.1	USER CONFIGURATION FOR VALIDATION	User settings are adjusted according to data of table				
3.2	(Domain) Log-in	No function can be actuated until a valid user is logged in				
		Successful login, functions can be induced				
		Only one user can be logged in at a time				
3.3	(Domain) Automatic log-out	Automatic logout after x minutes inactivity x =				
		Entry to audit trail that user has been logged out automatically				

No.	Test title	title Acceptance criteria		Results co	Results comply		
			(No.)	Test (yes/no)	Retest (yes/no)	Date/Signature	
3.4	(Domain) Log-out	No operation can be actuated					
		Entry to audit trail about the implemented steps					
4	Audit trail						
4.1	Audit trail entries (parameter)	All entries show the activity, user name, date and time correctly					
4.2	Audit trail entries	All entries in the audit trail show the activity, user name, date and time					
		The entries can be filtered and printed online based on various criteria (user, activity, date and time)					
		Entries can not be modified or deleted					
5	Version management	for recipe					
5.1	Version management in case of new recipe (product version)	When creating a new recipe, it will be generated with the designation count up number "x"; a new version with count up number "y"					
		After enabling the draft version (only) the status changes					
5.2	Disable recipe	Versions can be deleted or disabled					
	version (product data)	When the corresponding filters are selected, all existing versions are open for viewing					
		Entry to audit trail is correct					
5.3	Recipe versions for production	Disabled recipe versions can not be activated for production					
		Only recipe versions with status "Draft" or "Enabled" can be activated for production					

12/21

No.	Test title	Acceptance criteria	Deviation	Results comply			
			(No.)	Test (yes/no)	Retest (yes/no)	Date/Signature	
5.4	Create new recipe on the basis of an existing recipe version	 A new recipe can be created on the basis of an existing recipe version A new draft version can be created on the basis of a current recipe version 					
5.5	Delete disabled recipe version	 Status of selected recipe version changes to "Deleted" Deleted recipe versions can no longer be enabled or activated 					
6	Batch Management			1			
6.1	Create new batch	A new batch can only be created if previous batch has been completed					
6.2	Start batch	Only the new created batch can be started					
6.3	Interrupt and resume batch	Batch can be interrupted and resumed					
6.4	Complete batch	After batch is completed, batch data can be saved or printed and batch cannot be resumed again					
6.5	Batch history	"Batch history" displays when and by whom status of batch has been changed					
6.6	Store batch data	Batches are stored as .zip-file in selected folder with correct name (+ date and time)					
6.7	Restore batch data	Batch is restored from .zip-file and listed in "Batch management"					

No.	Test title	Acceptance criteria	Deviation	Results co	Results comply			
			(No.)	Test	Retest	Date/Signature		
				(yes/no)	(yes/no)			
6.8	Delete batch data	Batches with status "Restored" can be deleted from "Batch management"						
		Batches with status "Stored" can be deleted from "Batch management"						
		Batches with status "Completed" can not be deleted from "Batch management"						
6.9	Print history of batch data	Recipe data and results can be viewed and printed						
6.10	Print batch data	Batch data can be printed when active batch is completed or after a batch has been selected for history						
		Batch data can be printed on connected printer						
		Batch data can be exported to a selected folder as .pdf-file						
7	Reports and print outs			,		,		
7.1	Batch print out	Batch data can be printed (as pdf-file) after a batch has been selected for history						
		Batch data can be exported						
7.2	Parameter: Machine settings (Recipe)	Print out from the machine settings (Recipe) is possible						
7.3	User right print out	It is possible to print or save the user rights						
8	Basic function	1	I	<u> </u>		I		

No.	Test title	est title Acceptance criteria I	Deviation	Results comply			
			(No.)	Test (yes/no)	Retest (yes/no)	Date/Signature	
8.1	Service: Clean Display	Display turns grey and no touch input can be done					
		After 30 seconds the HMI appears again					
8.2	Alarmlists	Lists of the Faults, Warnings and Messages can be selected					
		Lists of Faults, Warnings and Messages can be printed or saved as pdf					
8.3	Operating hours counter	Operating hours counter is shown					
8.4	Parameter limits	A value greater than the maximum limit cannot be confirmed					
		A value smaller than the minimum limit cannot be confirmed					
		 A value between the minimum and maximum limit can be confirmed 					
8.5	HMI Language	The HMI language switches to the chosen language					

9 A	larm Fun	ection Tests				
9.1		Test protocols	- Function tests software s	witches operating mode (SWSOPM)		
Chap./ SWS	SOPM No).	Deviation No.	Chap/ SWSOPM No.	Deviation No.	
.2		Test protocols	- Function tests software s	witches operating mode drive (SWSOPM_I	DRV)	
hap./ SWS	SOPM No).	Deviation No.	Chap/ SWSOPM No.	Deviation No.	
.3		Test protocols	- Function tests software s	witches (SWS)		
Chap./ SWS No. Deviation No.		Deviation No.	Chap./ SWS No.	Deviation No.		

9.4	Test protocols	- Messages (ME)			
Chap./ ME No		Deviation No.	Chap./ ME No.	Deviation No.	
	1				
9.5	Test protocols	- Warnings (WA)			
		T			
Chap./ WA No).	Deviation No.	Chap./ WA No.	Deviation No.	
9.6	Test protocols	Faulte (FLT)			
3.0	Test protocols	- i auto (i Li)			
Chap./ FLT No	D.	Deviation No.	Chap./ FLT No.	Deviation No.	
		<u> </u>	1	III.	

Results comply

yes/no

Date/Initials

18/21

7.3 **Deviation sheet**

After execution the filled in deviation sheets will be attached to this document and the final number of attached pages has to be documented.

Empty forr	ns are in	chapter	7.5 "Fo	rms", p	age 21.
------------	-----------	---------	---------	---------	---------

Number of attached pages: _____

7.4 Change control sheet

After execution the filled in change control sheets will be attached to this document and the final number of attached pages has to be documented.

Empty forms are in chapter 7.5 "Forms", page 21.

Number of attached pages: _____

7.5 **Forms**

INFORMATION



Please use the following forms as master copies when a deviation sheet or a change control sheet is required.

Test phase:	Deviation sheet / Formular für Abweic	hungen	Deviation Abweichu	
Refering to test (binder/ch dazugehöriger Test (Ord				
I) Description of the Devi Beschreibung der Abweid				
Name:		Date / Datum:		Initials / Handzeichen:
II) Proposed solution / Vorgeschlagene Lösung	(Syntegon):			
Name		D-1- /		Leitiele I
Name:		Date / Datum:		Initials / Handzeichen:
III) Decision – follow up a Entscheidung - weitere A				
Proposed solutio Vorgeschlagene	n accepted / Lősung akzeptiert	NOT accepted NICHT akzept Comments / Kommentar:	1 / iert =>	
Name:		Date / Datum:		Initials / Handzeichen:
IV) Result (customer) / Ergebnis (Kunde):				
Conform and acc		NOT conform NICHT konfo		accepted / CHT akzeptiert =>
NOT conform, bu NICHT konform a Rationale / Begründung:	ut accepted / aber akzeptiert	Comments / Kommentare:		
Name:		Date / Datum:		Initials / Handzeichen:
V) Deviation closed (cust Abweichung vollständig t				
Name:	The second secon	Date / Datum:		Initials / Handzeichen:
Name:		Date / Datum:		Initials / Handzeichen:

Test phase:	Deviation sheet / Formular für Abweic	hungen	Deviation Abweichu	
Refering to test (binder/ch dazugehöriger Test (Ord				
I) Description of the Devi Beschreibung der Abweid				
Name:		Date / Datum:		Initials / Handzeichen:
II) Proposed solution / Vorgeschlagene Lösung	(Syntegon):			
Name		D-1- /		Leitiele I
Name:		Date / Datum:		Initials / Handzeichen:
III) Decision – follow up a Entscheidung - weitere A				
Proposed solutio Vorgeschlagene	n accepted / Lősung akzeptiert	NOT accepted NICHT akzept Comments / Kommentar:	1 / iert =>	
Name:		Date / Datum:		Initials / Handzeichen:
IV) Result (customer) / Ergebnis (Kunde):				
Conform and acc		NOT conform NICHT konfo		accepted / CHT akzeptiert =>
NOT conform, bu NICHT konform a Rationale / Begründung:	ut accepted / aber akzeptiert	Comments / Kommentare:		
Name:		Date / Datum:		Initials / Handzeichen:
V) Deviation closed (cust Abweichung vollständig t				
Name:	The second secon	Date / Datum:		Initials / Handzeichen:
Name:		Date / Datum:		Initials / Handzeichen:

Change Formula	Change Control N°							
	Refering documentation (binder/chapter) / dazugehörige Dokumentation (Ordner/Kapitel):							
	I) Description / Beschreibung:							
	n for change / der Ånderung:							
	Re-Qualification necessary, if yes, Requalifizierung nötig, wenn ja, wie							
Nama:		Date /	Initials /					
Name:		Datum:	Handzeichen:					
	tomer's comments / entar des Kunden:							
	Change accepted / Change akzeptiert	Comments / Kommentar:						
•	NOT accepted, see comments / NICHT akzeptiert, siehe Kommentar							
Name:		Date / Datum:	Initials / Handzeichen:					
	low up action /		Trial National State of the Sta					
Aktion	sverfolgung:							
	Re-Qualification executed / Requalifizierung durchgeführt							
Name:		Date / Datum:	Initials / Handzeichen:					
IV) Rev	view:							
Name:		Date / Datum:	Initials / Handzeichen:					
Name:		Date / Datum:	Initials / Handzeichen:					

Change Formula	Change Control N°							
	Refering documentation (binder/chapter) / dazugehörige Dokumentation (Ordner/Kapitel):							
	I) Description / Beschreibung:							
	n for change / der Ånderung:							
	Re-Qualification necessary, if yes, Requalifizierung nötig, wenn ja, wie							
Nama:		Date /	Initials /					
Name:		Datum:	Handzeichen:					
	tomer's comments / entar des Kunden:							
	Change accepted / Change akzeptiert	Comments / Kommentar:						
•	NOT accepted, see comments / NICHT akzeptiert, siehe Kommentar							
Name:		Date / Datum:	Initials / Handzeichen:					
	low up action /		Trial National State of the Sta					
Aktion	sverfolgung:							
	Re-Qualification executed / Requalifizierung durchgeführt							
Name:		Date / Datum:	Initials / Handzeichen:					
IV) Rev	view:							
Name:		Date / Datum:	Initials / Handzeichen:					
Name:		Date / Datum:	Initials / Handzeichen:					

Alarm and Function Testing Execution

CUC 2002 Cartoning machine

100223

Glaxowellcome Production France

Table of contents

1	Software version	
2	Check HMI menu, buttons and switches	1
2.1	Spotcheck HMI	1
3	Access protection	1
3.1	USER CONFIGURATION FOR VALIDATION	1
3.2	(Domain) Log-in	1
3.3	(Domain) Automatic log-out	1
3.4	(Domain) Log-out	1
4	Audit trail	1
4.1	Audit trail entries (parameter)	1
4.2	Audit trail entries	1
5	Version management for recipe	2
5.1	Version management in case of new recipe (product version)	2
5.2	Disable recipe version (product data)	2
5.3	Recipe versions for production	2
5.4	Create new recipe on the basis of an existing recipe version	2
5.5	Delete disabled recipe version	2
6	Batch Management	2
6.1	Create new batch	2
6.2	Start batch	3
6.3	Interrupt and resume batch	3
6.4	Complete batch	3
6.5	Batch history	3
6.6	Store batch data	3
6.7	Restore batch data	3
6.8	Delete batch data	3
6.9	Print history of batch data	4
6.10	Print batch data	4
7	Reports and print outs	4
7.1	Batch print out	4
7.2	Parameter: Machine settings (Recipe)	4
7.3	User right print out	4
8	Basic function	4
8.1	Service: Clean Display	4

QD 2 Alarm and Function	Testing	Execution
-----------------------------	---------	-----------

8.2	Alarmlists
8.3	Operating hours counter
8.4	Parameter limits
8.5	HMI Language
9	Alarm Function Tests
9.1	Test protocols - Function tests software switches operation mode (SWSOPM)
9.1.1	SWSOPM 1: AUTOMATIC
9.1.2	SWSOPM 2: SETUP
9.1.3	SWSOPM 3: INCHING PANEL
9.1.4	SWSOPM 4: INCHING CABLE
9.1.5	SWSOPM 6: OPEN STOP BRAKE
9.1.6	SWSOPM 7: ZERO DRIVE
9.1.7	SWSOPM 8: INCHING SINGLE AXIS
9.1.8	SWSOPM 11: INCHING INFEED
9.1.9	SWSOPM 60: LINE MODE DOWNSTREAM EQUIPMENT
9.1.10	SWSOPM 62: MACHINE: AUTOMATIC RESTART
9.1.11	SWSOPM 63: LINE MODE UPSTREAM EQUIPMENT
9.2	Test protocols - Function tests software switches operation mode drive (SWSOPM_DRV)
9.2.1	SWSOPM_DRV 1002: INSERTION
9.2.2	SWSOPM_DRV 1003: PREINSERTION HORIZONTAL
9.2.3	SWSOPM_DRV 1004: PREINSERTION VERTICAL
9.2.4	SWSOPM_DRV 1007: CARTON SETUP
9.2.5	SWSOPM_DRV 1008: LEAFLET DEVICE
9.2.6	SWSOPM_DRV 1015: MAIN DRIVE
9.2.7	SWSOPM_DRV 1016: INSERTION TABLE
9.2.8	SWSOPM_DRV 1017: CARTON ALIGNMENT
9.2.9	SWSOPM_DRV 1054: ROBOT 1: DRILL AXIS
9.2.10	SWSOPM_DRV 1055: ROBOT 1
9.2.11	SWSOPM_DRV 1059: ROBOT 2: DRILL AXIS
9.2.12	SWSOPM_DRV 1060: ROBOT 2
9.2.13	SWSOPM_DRV 1061: ROBOT 1: CYCLE BELT 1
9.2.14	SWSOPM_DRV 1062: ROBOT 1: CYCLE BELT 2
9.3	Test protocols - Function tests software switches (SWS)
9.3.1	SWS 52: COMPRESSED AIR
9.3.2	SWS 53: PERMANENT READING CODE
9.3.3	SWS 54: VACUUM
9.3.4	SWS 56: MACHINE INSIDE LIGHT OFF
9.3.5	SWS 100: INSERT DEFECTIVE PRODUCT
9.3.6	SWS 106: CONTINUOUS INSERTION
9.3.7	SWS 119: CONTINUOUS PREINSERTION
9.3.8	SWS 120: PREINSERTION UPPER FRONT POSITION
9.3.9	SWS 121: PREINSERTION PLATE 1 DOWN
9.3.10	SWS 122: PREINSERTION PLATE 2 DOWN

QD 2	Alarm and Function Testing Execution	4/317
9.3.11	SWS 123: PREINSERTION PLATE 3 DOWN	110
9.3.12	SWS 124: PREINSERTION PLATE 4 DOWN	112
9.3.13	SWS 127: PREINSERTION	114
9.3.14	SWS 200: CODE READER CARTON	116
9.3.15	SWS 201: CARTON LASER PRINTING	118
9.3.16	SWS 203: BLOWING AIR CARTON	120
9.3.17	SWS 208: CONTINUOUS CALL CARTON	122
9.3.18	SWS 256: CHECK CARTON OPEN COVER FLAP	124
9.3.19	SWS 257: CHECK CARTON OPEN SIDE FLAP	126
9.3.20	SWS 300: LEAFLET DEVICE	128
9.3.21	SWS 302: CODE READER LEAFLET REAR	129
9.3.22	SWS 303: CHECK LEAFLET PRESENT	131
9.3.23	SWS 304: INSERT DEFECTIVE LEAFLET	133
9.3.24	SWS 307: CONTINUOUS CALL LEAFLET	135
9.3.25	SWS 400: ROBOT 1	137
9.3.26	SWS 401: ROBOT 1: RUN EMPTY	138
9.3.27	SWS 405: ROBOT 1: CALL PRODUCT	139
9.3.28	SWS 500: ROBOT 2	140
9.3.29	SWS 505: ROBOT 2: CALL PRODUCT	141
9.3.30	SWS 900: BROCHURE INFEED	142
9.3.31	SWS 901: BROCHURE INFEED: CODE READER	143
9.3.32	SWS 902: BROCHURE INFEED: CONTINUOUS CALL	145
9.3.33	SWS 2001: ROBOT 2: RUN WITH FLOW-PACK	147
9.3.34	SWS 2002: CODE READER GLUED-IN BROCHURE	149
9.4	Test protocols - Messages (ME)	151
9.4.1	ME 25: BATCH CONTROL: BATCH IS INTERRUPTED	151
9.4.2	ME 26: BATCH CONTROL: BATCH IS NOT ACTIVE	152
9.4.3	ME 27: BATCH CONTROL: BATCH IS ACTIVE	153
9.4.4	ME 40: OPERATING MODE: AUTOMATIC	154
9.4.5	ME 41: OPERATING MODE: SETUP	155
9.4.6	ME 42: OPERATING MODE: INCHING MODE PANEL	156
9.4.7	ME 43: OPERATING MODE: INCHING MODE CABLE FRONT	157
9.4.8	ME 47: OPERATING MODE: INCHING MODE SINGLE AXES	158
9.4.9	ME 49: OPERATING MODE: INCHING MODE INFEED	159
9.4.10	ME 160: MACHINE: OPERATIONAL	160
9.4.11	ME 448: ROBOT 1: REFERENCE RUN	161
9.4.12	ME 548: ROBOT 2: REFERENCE RUN	162
9.5	Test protocols - Warnings (WA)	163
9.5.1	WA 27: CONTROL-PC: M-GUARD ENABLED	163
9.5.2	WA 28: CONTROL-PC: REMOTE MAINTENANCE ENABLED	164
9.5.3	WA 29: OPERATOR: NO USER LOGGED IN	165
9.5.4	WA 45: CONTROL CABINET: UPS NOT READY	166
9.5.5	WA 220: CARTON: PREWARNING LOW STOCK - IB32_B01	167

9.5.6	WA 380: LEAFLET INSERTER TURRET MAGAZINE: PREWARNIN LOW LEVEL	168
9.5.7	WA 381: LEAFLET INSERTER TURRET MAGAZINE: NOT READY	169
9.5.8	WA 800: LEAFLET INSERTER : MINIMAL SUPPLY	170
9.5.9	WA 900: BROCHURE INFEED: MINIMAL SUPPLY	171
9.5.10	WA 980: BROCHURE INFEED: TURRET MAGAZINE: MINIMAL SUPPLY	172
9.5.11	WA 981: BROCHURE INFEED: TURRET MAGAZINE: NOT IN OPERATION MODE	173
9.5.12	WA 2001: ROBOT 1: UPLINE MACHINE IS NOT RUNNING	174
9.5.13	WA 2002: ROBOT 2: UPLINE MACHINE IS NOT RUNNING	175
9.5.14	WA 2003: ROBOT 2: DOWNLINE MACHINE IS NOT READY	176
9.6	Test protocols - Faults (FLT)	177
9.6.1	FLT 3: OPERATOR: NORMAL STOP ACTIVATED	177
9.6.2	FLT 4: OPERATOR: INCORRECT OPERATING MODE SELECTED	178
9.6.3	FLT 5: COMPRESSED AIR: MONITORING INLET PRESSURE	179
9.6.4	FLT 6: ETHERCAT: FAULT	180
9.6.5	FLT 7: COMPRESSED AIR: MONITORING SAFE SHUTDOWN	181
9.6.6	FLT 11: MOTION DRIVES: OVERLOAD POWER SUPPLY - IK00_Q10	182
9.6.7	FLT 17: OPERATOR: OPERATING MODE CHANGE	183
9.6.8	FLT 22: VACUUM: OVERLOAD PUMP	184
9.6.9	FLT 23: VACUUM: MONITORING NEGATIVE PRESSURE	185
9.6.10	FLT 34: MACHINE: REPAIR SWITCH DRIVES SWITCHED OFF - IW00_Q02	186
9.6.11	FLT 35: CONTROL CABINET: MONITORING MAIN SWITCH - IW00_Q01	187
9.6.12	FLT 42: CONTROL CABINET S1: OVERTEMPERATURE - IW00_B13	188
9.6.13	FLT 79: EMERGENCY STOP ACTIVATED: LINE	189
9.6.14	FLT 81: EMERGENCY STOP ACTIVATED: OPERATING PANEL - IF91_S90	190
9.6.15	FLT 86: GUARD OPEN: BOTTOM [7]	191
9.6.16	FLT 87: GUARD OPEN: CARTON BELT [5]	192
9.6.17	FLT 88: GUARD OPEN: INSERTION REAR [20]	193
9.6.18	FLT 111: PRODUCT SENSING: CONSECUTIVE FAULT - IB10_B01	194
9.6.19	FLT 112: PRODUCT SENSING: SENSOR MONITORING - IB10-B01	195
9.6.20	FLT 120: INSERTION: SAFETY SENSOR - IB36_B01	196
9.6.21	FLT 122: PRODUCT CHAIN: OVERFILL COLLECTING BOX - IB80_B02	197
9.6.22	FLT 126: INSERTION: OVERLOAD PREINSERTION	198
9.6.23	FLT 129: INSERTION: MONITORING SENSOR SAFETY SENSOR - IB36_B01	199
9.6.24	FLT 200: MACHINE: JAM AT DISCHARGE	200
9.6.25	FLT 201: DOWNSTREAM: STOP FROM MACHINE 1	201
9.6.26	FLT 203: DOWNSTREAM: STOP FROM MACHINE 2	202
9.6.27	FLT 210: GUARD OPEN: CARTON MAGAZINE [6]	203
9.6.28	FLT 212: GUARD OPEN: DISCHARGE TOP [2]	204
9.6.29	FLT 213: CARTON: OVERLOAD SIDE FLAP FOLDER - IB37_B01	205
9.6.30	FLT 220: CARTON: STOCK AT THE END - IA30_B02	206
9.6.31	FLT 221: MACHINE: COVER RAIL OPEN - IA30-B03	207
9632	FLT 222: CARTON: CONSECUTIVE FAULT CROSSCHECK - IA30, B01	208

QD 2 Alarm and Function Te	esting Execution
--------------------------------	------------------

9.6.33	FLT 223: CARTON: PICKUP CHECK - IA30_B04	209
9.6.34	FLT 224: CARTON BELT: OVERLOAD	210
9.6.35	FLT 225: CARTON: CONSECUTIVE FAULT CODE READER - IB75_B01	211
9.6.36	FLT 226: CARTON:\FAULT CODE READER - IB75_B01	212
9.6.37	FLT 227: CARTON: SENSOR MONITORING CROSSCHECK - IA30_B01	213
9.6.38	FLT 237: CARTON: SENSOR MONITORING PICKUP CHECK - IA30_B04	214
9.6.39	FLT 243: LASER PRINTING: DEVICE NOT READY	21
9.6.40	FLT 250: CARTON EJECTION 1: OVERFILL COLLECTING CONTAINER - IQ82_B03	216
9.6.41	FLT 251: CARTON EJECTION 1: CROSS CHECK - IQ82_B01	217
9.6.42	FLT 254: CARTON REJECT 1: SENSOR MONITORING CROSS CHECK - IQ82_B01	218
9.6.43	FLT 270: CARTON: OPEN FLAP: SENSOR MONITORING COVER FLAP FRONT - IB86_B01	219
9.6.44	FLT 271: CARTON: OPEN FLAP: SENSOR MONITORING COVER FLAP REAR - IB86_B02	220
9.6.45	FLT 272: CARTON: OPEN FLAP: SENSOR MONITORING SIDE FLAP FRONT - IB86_B03	22
9.6.46	FLT 273: CARTON: OPEN FLAP: SENSOR MONITORING SIDE FLAP REAR - IB86_B04	222
9.6.47	FLT 289: CARTON: OPEN FLAP: CONSECUTIVE FAULT COVER FLAP -IB86_B01 / IB86_B02	223
9.6.48	FLT 294: CARTON: OPEN FLAP: CONSECUTIVE FAULT SIDE FLAP -IB86_B03 / IB86_B04	22
9.6.49	FLT 301: LEAFLET: CONSECUTIVE FAULT CROSS CHECK IN PRODUCT CHAIN - IX161_B10	227
9.6.50	FLT 302: LEAFLET: CROSS CHECK EJECTION - IX161_B12	229
9.6.51	FLT 304: LEAFLET: CONSECUTIVE FAULT CROSS CHECK IN UNIT - IX161_B11	230
9.6.52	FLT 321: LEAFLET: CONSECUTIVE FAULT CODE READER REAR - IB75_B05	232
9.6.53	FLT 322: LEAFLET: CONSECUTIVE FAULT CROSS CHECK IN CARTON - IB75_B04	234
9.6.54	FLT 324: LEAFLET: FAULT CODE READER REAR - IB75_B05	236
9.6.55	FLT 325: LEAFLET: SENSOR MONITORING CROSS CHECK IN CARTON - IB75-B04	237
9.6.56	FLT 326: LEAFLET: SENSOR MONITORING CROSS CHECK IN PRODUCT CHAIN - IX161_B10	238
9.6.57	FLT 408: ROBOT 1: SENSOR MONITORING PRODUCT SENSING EXCEED HEIGHT	239
9.6.58	FLT 412: ROBOT 1: EMERGENCY STOP FRONT ACTIVATED	240
9.6.59	FLT 421: ROBOT 1: GUARD OPEN: FRONT LEFT [17]	241
9.6.60	FLT 422: ROBOT 1: GUARD OPEN: FRONT RIGHT [16]	242
9.6.61	FLT 430: ROBOT 1: BELT 1: OVERLOAD	243
9.6.62	FLT 432: ROBOT 1: BELT 2: OVERLOAD	244
9.6.63	FLT 450: ROBOT 1: PRODUCT SENSING: EXCEED HEIGHT	245
9.6.64	FLT 451: ROBOT 1: BELT 1: PRODUCT JAM	246
9.6.65	FLT 457: ROBOT 1: MONITORING COMPRESSED AIR	247
9.6.66	FLT 477: ROBOT 1: CONTROL CABINET: OVERTEMPERATURE	248
9.6.67	FLT 512: ROBOT 2: EMERGENCY STOP FRONT ACTIVATED	249

QD 2	Alarm and Function Testing Execution	7/317
9.6.68	FLT 521: ROBOT 2: GUARD OPEN: FRONT LEFT [19]	250
9.6.69	FLT 522: ROBOT 2: GUARD OPEN: FRONT RIGHT [18]	251
9.6.70	FLT 547: ROBOT 2: EMPTYING	252
9.6.71	FLT 577: ROBOT 2: CONTROL CABINET: OVERTEMPERATURE	253
9.6.72	FLT 807: LEAFLET INSERTER: NOT IN OPERATIONAL	254
9.6.73	FLT 808: LEAFLET INSERTER: GUARD OPEN: REAR [23]	255
9.6.74	FLT 822: LEAFLET INSERTER: EMERGENCY STOP ACTIVATED	256
9.6.75	FLT 843: LEAFLET INSERTER: CONTROL CABINET: VOLTAGE MONITORING	257
9.6.76	FLT 902: BROCHURE INFEED: CONSECUTIVE FAULT CODE READER	258
9.6.77	FLT 903: BROCHURE INFEED: CONSECUTIVE FAULT CROSS CHECK	259
9.6.78	FLT 907: BROCHURE INFEED: NOT IN OPERATIONAL	261
9.6.79	FLT 916: BROCHURE INFEED: GUARD OPEN: PRODUCT CHAIN INFEED [14]	262
9.6.80	FLT 920: BROCHURE INFEED: NOT IN CORRECT POSITION [15]	263
9.6.81	FLT 922: BROCHURE INFEED: EMERGENCY STOP ACTIVATED	264
9.6.82	FLT 943: BROCHURE INFEED: CONTROL CABINET: VOLTAGE MONITORING	265
9.6.83	FLT 1212: GUARD OPEN: COLLECTING CONTAINER 1 [3]	266
9.6.84	FLT 1213: CARTON EJECTION 1: OVERFILL GUARD COLLECTING CONTAINER	267
9.6.85	FLT 1230: DOWNSTREAM: GUARD OPEN: MACHINE 1	268
9.6.86	FLT 1413: ROBOT 1: REMOVE PRODUCT BEFORE REFERENCE/ SYNCHRON RUN	269
9.6.87	FLT 1416: ROBOT 1: EMERGENCY STOP REAR ACTIVATED	
9.6.88	FLT 1417: ROBOT 1: GUARD OPEN: REAR LEFT [25]	
9.6.89	FLT 1418: ROBOT 1: GUARD OPEN: REAR RIGHT [24]	
9.6.90	FLT 1419: ROBOT 1: RANGE LIMIT REACHED	273
9.6.91	FLT 1516: ROBOT 2: EMERGENCY STOP REAR ACTIVATED	274
9.6.92	FLT 1517: ROBOT 2: GUARD OPEN: REAR LEFT [27]	275
9.6.93	FLT 1518: ROBOT 2: GUARD OPEN: REAR RIGHT [26]	
9.6.94	FLT 1519: ROBOT 2: RANGE LIMIT REACHED	
9.6.95	FLT 2001: GUARD OPEN: BACK SIDE BELOW LOADING [21]	

FLT 2002: GUARD OPEN: BACK SIDE VACUUM PUMP [22]

FLT 2004: GUARD OPEN: PRODUCT EJECTION COLLECTING BOX [4]

FLT 2009: LEAFLET INSERTER: POSITION CHECK LEAFLET INSERTER

9.6.102 FLT 2012: LEAFLET INSERTER: GUARD OPEN: TOP LEFT [9]

9.6.103 FLT 2013: LEAFLET INSERTER: GUARD OPEN: TOP RIGHT [8]

9.6.104 FLT 2014: LEAFLET INSERTER: POSITION CHECKINFEED BELT [10]

9.6.96

9.6.97

9.6.98

9.6.99

9.6.105

FLT 2015: LEAFLET INSERTER: EMERGENCY STOP

9.6.106 FLT 2016: BROCHURE INFEED: EMERGENCY STOP

ACTIVATED: TURRET MAGAZINE

ACTIVATED: TURRET MAGAZINE

FLT 2005: LASER: SENSOR MONITORING PRODUCT TRIGGER

279

280

281

282

283

284

285

286

287

288

289

QD 2 Alarm and Function T	esting Execution
-------------------------------	------------------

<u>QD Z </u>	Alaim and Function Testing Execution	0/01/
9.6.107	FLT 2017: BROCHURE INFEED: COVER OPEN: TRANSFER SECTION [13]	290
9.6.108	FLT 2019: GUARD OPEN: LASER DOOR [1]	291
	FLT 2021: EMERGENCY STOP ACTIVATED: LASER	
9.6.110	FLT 2022: LEAFLET: MONITORING PROPER SIGNAL CODE READING LEAFLET REAR - IB75_B05	293
9.6.111	FLT 2023: BROCHURE INSERTER: MONITORING PROPER SIGNAL CODE READING BROCHURE - IB75_B06	294
9.6.112	FLT 2024: BROCHURE: MONITORING PROPER SIGNAL CODE READING BROCHURE IN CARTON - IB75_B07	295
9.6.113	FLT 2025: CODE READER GLUED-IN BROCHURE NOT READY	296
9.6.114	FLT 2026: GLUED IN BROCHURE: CONSECUTIVE FAULT CODE READER - IB75_B07	297
9.6.115	FLT 2029: LASER: PRODUCT JAM	298
9.6.116	FLT 2036: ROBOT 1: CROSS CHECK EJECTOR	299
9.6.117	FLT 2037: ROBOT 1: EXCESS FILL COLLECTING BOX ACCELERATION BELT	300
9.6.118	FLT 2038: ROBOT 1: EXCESS FILL COLLECTING BOX CYCLE BELT	301
9.6.119	FLT 2040: ROBOT 2: EMERGENCY STOP ACTIVATED: PUCK OUTFEED BELT	302
9.6.120	FLT 2041: ROBOT 2: EMERGENCY STOP ACTIVATED: DISKUS BELT	303
9.6.121	FLT 2043: ROBOT 2: PUCK OUTFEED BELT: OVERLOAD	304
9.6.122	FLT 2044: ROBOT 2: PUCK INFEED BELT: OVERLOAD	305
9.6.123	FLT 2045: ROBOT 2: DISKUS BELT: OVERLOAD	306
9.6.124	FLT 2046: ROBOT 2: PUCK NOT EMPTY AT OUTFEED BELT	307
9.6.125	FLT 2055: ROBOT 1: MONITORING SENSOR PRODUCT CYCLE BELT -W166_B06	308
9.6.126	FLT 2056: ROBOT 1: MONITORING SENSOR FLOWPACK CYCLE BELT -W166_B07	309
9.6.127	FLT 2057: CARTON REJECT 1: CROSS CHECK CARTON REJECT IQ82_B10	310
9.6.128	FLT 2058: CARTON REJECT 1: SENSOR MONITORING CROSS CHECKCARTON REJECTED - IQ82_B10	312
9.6.129	FLT 2059: ROBOT 1: MONITORING SENSOR DISTANCE CONTROL 1 INFEED BELT -W166-B01	313
9.6.130	FLT 2060: ROBOT 1: MONITORING SENSOR DISTANCE CONTROL 2 INFEED BELT -W166-B02	314
9.6.131	FLT 2061: INSERTION: MONITORING SENSOR: OVERLOAD PREINSERTION -W40-B02	315
9.6.132	FLT 2062: ROBOT 1: MONITORING SENSOR: EXCESS FILL COLLECTION BOX CYCLE BELT -W166_B10	316
10	Notes	247

1 Software version

Test objective	Documentation of software version used for test execution
Test procedure	
Required operations	Write down in table below the installed software version (Baseline) before execution of alarm and function testing.
Software/Firmware	Version
Baseline	
Test result	yes/no
Acceptance criteria	The software version is documented.
Comments	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

2 Check HMI menu, buttons and switches

2.1 Spotcheck HMI

Test objective	Menu check of HMI menu tree	
Test procedure		
Test prerequisites	Machine is ready in automatic modeUser Admin is logged in	
Required operations	 Copy the menu tree from the operation manual chap. 5 Select longest tree as example (or one of the longest trees) Go step by step through the selected tree Control the menu and the names of the page tabs Attach the copy to this test 	
Comments	This is a spotcheck test. Only the longest menu tree needs to be check	ked
Acknowledgement	None	
Test result		yes/no
Acceptance criteria	HMI menu and page tab spotcheck is OK	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

3 Access protection

USER CONFIGURATION FOR VALIDATION 3.1

Test objective	This configuration is neccessary to create an user for the following varieties procedure	validation
Test procedure		
Required operations	Log in with user "adminlocal"	
	Select the menu "System configuration"	
	Select the submenu "User management"	
	 In submenu "User Management", press "General settings" Create following configuration: 	
	GENERAL SETTINGS	
	User	
	Display last user (min.) 0	
Test result		yes/no
Acceptance criteria	User settings are adjusted according to data of table	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

3.2 (Domain) Log-in

Test objective	No action without logged in user
,	Login requires the individual user name and password
	Checking the number of user-login actions
	Chooking the hamber of according to
Test procedure	
Test prerequisites	 Adminlocal has created a new group to execute all the following tests Name of the group Machine is operational No user is logged in
Required operations	 Press "Start" Warning 29 "Operator: No user logged in" is indicated Touch any function key (no function can be actuated) Call up function "Log in user" Log in with user name and password Operate any function key (function is actuated) Try to log in a second user with name and password (not possible)
Comments	Activation of machine stop (Normal Stop) and Emergency-Stop possible without logged in user
Acknowledgement	• None

Test result		yes/no
Acceptance criteria	No function can be actuated until a valid user is logged in	
	Successful login, functions can be induced	
	Only one user can be logged in at a time	
Comments		

Results comply	yes/no Date/Initials
D	
Results approved	Date/Initials

3.3 (Domain) Automatic log-out

Test objective	Check of access protection: Time log-out	
Test procedure		
Test prerequisites	Machine is operationalNo user is logged in	
Required operations	 Call up function "Log in user" Enter user name and password No action for x minutes Time of last action: Enter user name and password Select submenu "Diagnostics" → press "Audit view" Check in audit trail the correct entry with userID, date, time and actiog-out" 	ition "Auto
Comments	Print out audit trail after having finished all tests	
Acknowledgement	• None	
Test result		yes/no
Acceptance criteria	 Automatic logout after x minutes inactivity x = Entry to audit trail that user has been logged out automatically 	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

3.4 (Domain) Log-out

Test objective	Check of User log-out	
Test procedure		
Test prerequisites	Machine is operationalUser is logged in	
Required operations	Log out user	
	Time of log-out:	
	Press any function key - no function can be activated	
	Log in user > enter user name and password	
	 Select submenu "Diagnostics" → "Audit view" 	
	 Check in audit trail the correct entry with userID, date, time and action "Auto log-out" 	
Comments	Activation of machine stop (Normal Stop) and Emergency-Stop possible without logged in user	
Acknowledgement	• None	
Test result	yes/ne	0
Acceptance criteria	No operation can be actuated	
	Entry to audit trail about the implemented steps	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

4 Audit trail

4.1 Audit trail entries (parameter)

Test objective	All user actions related to process, quality and product as well as all parameter changes must be recorded in the audit trail	
Test procedure		
Test prerequisites	 User "adminlocal" is logged in Machine is ready in automatic mode Date and time are identical to external clock 	
Required operations	Time code Confirmation in audit trail	
	Press Emergency Stop	
	Release Emergency Stop	
	Complete the batch	
	 Select submenu "Diagnostics" > press "Audittrail viewer" > press "Display/ Viewer" 	
	 Check entries for correctness and confirm with "\u03c4" 	
	Print out the "Audit trail" and attach the document to the test protocol	
Comments	 21 CFR Part 11: §11.10(d,g), §11.200(a), §11.300(a,b,d) GMP Vol.4 Annex 11: 12.1; 12.3 	
Acknowledgement	• None	

Test result		yes/no
Acceptance criteria	All entries show the activity, user name, date and time correctly	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

4.2 Audit trail entries

Test objective	Structure of entries
Test procedure	
Test prerequisites	Machine is operational
Required operations	 Log in user "adminlocal" Log out user "adminlocal" Open guard door Close guard door: Press "Reset" Select "Diagnostics" Call up function "Audittrail viewer" Press "Filter setting" Select filter "User name" by activating the checkmark Press "Setting" within the filter setting "User name" Move one of the given users from field "All" to the field "Selected" (e.g. validation) Select filter "Time period" Press "Setting" within the filter setting "Time period" Adjust the time period to the last 10 minutes Confirm with the checkmark Close the menu "Filter setting" by using of the checkmark Press "Display" within the "Audittrail viewer" The corresponding audit view of user "adminlocal" is displayed Call up function "Print" The audit trail printout displays the selected date and time Select filter "Faults" Confirm with the checkmark Press "Refresh" The corresponding audit view is displayed Try to modify or delete last entry in audit trail Print out audit trail after completion of the tests
Consequence	 All entries in the audit trail show the activity, user name, date and time The entries can be filtered and printed online based on various criteria (user, activity, date and time, batch) Entries can not be modified or deleted

Test procedure		
Comments	 Select an audit trail in "Audittrail viewer" within a time frame showin 21 CFR Part 11: §11.10(d,g), §11.200(a), §11.300(a,b,d) GMP Vol.4 Annex 11: 12.1; 12.3 	g actions
Acknowledgement	 Remove "Filter setting" Close "Audit trail viewer" Several filter settings can be activated for the next start of the "Audiviewer" 	ittrail
Test result		yes/no
Acceptance criteria	All entries in the audit trail show the activity, user name, date and time	
	The entries can be filtered and printed online based on various criteria (user, activity, date and time)	
	Entries can not be modified or deleted	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

5 Version management for recipe

5.1 Version management in case of new recipe (product version)

Test objective	All recipe data are managed in menu "Recipe" > submenu "Recipe management"
Test procedure	
Test prerequisites	Machine is operationalUser "adminlocal" is logged in
Required operations	 Select menu "Recipe" > select submenu "Recipe management" Select one recipe and press "Versions" Select the current version and press "Create new recipe" Recipe version number "x" = is created Change the recipe name of recipe version "x" (e.g. "Validation A 500 mg") Call up the new generated recipe in the window "recipe version" (e.g. "Validation A 500 mg") Select a version and press "New draft version" A new version draft with count up number "y" = Select draft "y" and press "Enable" Confirm message Select the current recipe version in the window Press "Activate" In the overview "Current recipe" the recipe is shown (e.g. Validation A 500 mg) Version x.y = (e.g. 0002.0001) released
Comments	 21 CFR Part 11: §11.10 (b) GMP Vol.4 Annex 11: 8.1

Test result		yes/no
Acceptance criteria	When creating a new recipe, it will be generated with the designation count up number "x"; a new version with count up number "y"	
	After enabling the draft version (only) the status changes	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

5.2 Disable recipe version (product data)

Test objective	Check whether recipe version can be deleted
Test procedure	
Test prerequisites	Machine is operationalUser "adminlocal" is logged in
Test 1	Delete a version
Required operations	 Select menu "Recipe" > selet submenu "Select recipe management" Select the new recipe (e.g. Validation A 500 mg) Press "Versions" Create a new draft version Select the new draft version Press "Delete" and confirm (note time of action) Leave menu "Recipe versions product" Leave menu "Recipe list product" Select menu "System configuration" > select submenu "Diagnostic" Select "Audittrail viewer" Select "Filter settings" Activate checkmarks for "Function" Select "Function setting" Select "Size management" and change to selected Select "Size management" and press "Events" Select "Recipe deleted" and change to selected Confirm tree times with "Checkmark" and press button "Display" Check the entry in the audit trail Print out audit trail after completion of the tests
Consequence	 Versions can be deleted When the corresponding filters are selected, all existing versions are open for viewing Entry to audit trail is correct
Comments	 If no active version is available create a "New draft version" Inactive versions can be deleted by the Admin only Enabled recipes can not be deleted
Acknowledgement	None

Test procedure	
Test 2	Block a recipe version
Required operations	 Select menu "Recipe" > selet submenu "Recipe management" Select the new recipe (e.g. Validation A 500 mg) Press "Versions" Create a new draft version of the draft version Select the new draft version Press "Disable" (note time of action) Press filter "Disable" Leave menu "Recipe versions product" Leave menu "Recipe list product" Select menu "System configuration" > select submenu "Diagnostic" Select "Audittrail viewer" Select "Filter settings" Activate checkmarks for "Function" Select "Function setting" Select "Size management" and change to selected Select "Recipe disabled" and change to selected Confirm tree times with "Checkmark" and press button "Display" Check the entry in the audit trail Print out audit trail after completion of the tests
Consequence	 Versions can be blocked When the corresponding filters are selected, all existing versions are open for viewing Entry to audit trail is correct
Comments	If necessary create a "New draft version" to perform the test
Acknowledgement	• None

Test result		yes/no
Acceptance criteria	Versions can be deleted or disabled	
	When the corresponding filters are selected, all existing versions are open for viewing	
	Entry to audit trail is correct	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

5.3 Recipe versions for production

Test objective	Check that only recipe versions with status "Draft" or "Enabled" can be activated for production
Test procedure	
Test prerequisites	 User "adminlocal" is logged in Machine is ready in automatic mode Recipe e.g. "Validation A 500mg" is created "Recipe Version" is created as disabled "Recipe Version" is created as deleted "Recipe Version" is created as enabled version "Recipe Version" is created as draft version
Required operations	 Select menu "recipe management" > select "recipe management" > SWS "recipe management" An overview of all recipes is displayed In overview window of recipes select recipe e.g. "Validation A 500mg" and press SWS "Versions" All filters are switched on Select a recipe version with status "Disabled" SWS "Activate" is not active
Consequence	Selected recipe version cannot be activated
Required operations	Select a recipe version with status "Deleted"SWS "Activate" is not active
Consequence	Selected recipe version cannot be activated
Required operations	 Select a recipe version with status "Enabled" Press SWS "Activate"
Consequence	 Activated version is highlighted in green All settings for this recipe are taken over by program
Required operations	 Select a recipe version with status "Draft" Press SWS "Activate"
Consequence	 Activated version is highlighted in green All settings for this recipe are taken over by program
Required operations	Press SWS "Print" and attach printout to the test *)
Comments	*) If no printer is available create pdf-file and print on an external printer

Test result		yes/no
Acceptance criteria	Disabled recipe versions can not be activated for production	
	Only recipe versions with status "Draft" or "Enabled" can be activated for production	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

Create new recipe on the basis of an existing recipe version 5.4

Test objective	Create new recipe on the basis of an existing recipe version
Test procedure	
Test prerequisites	 Machine is operational User "adminlocal" is logged in
Required operations	 Select menu "Recipe" > select submenu "Select recipe management" Select one recipe and press "Versions" Select the current version and press "Create new recipe"
Consequence	 Recipe version number "x" = is created A new recipe can be created on the basis of an existing recipe version
Required operations	 Press "Change Recipe name", change the name of recipe version "x" (e.g. "Validation A 500 mg") Call up the new generated recipe in the window "recipe version" (e.g. "Validation A 500 mg") Select a version and press "New draft version"
Consequence	 A new version draft with count up number "y" = A new draft version can be created on the basis of a current recipe version
Test result	yes/no
Acceptance criteria	 A new recipe can be created on the basis of an existing recipe version A new draft version can be created on the basis of a current recipe version
Comments	
Results comply	yes/no Date/Initials
Results approved	Date/Initials

5.5 Delete disabled recipe version

Test objective	
rest objective	Check whether disabled recipe version can be deleted
Test procedure	
Test prerequisites	 User "adminlocal" is logged in Machine is ready in automatic mode Recipe version to be deleted is not active (e.g. "Validation A 500 mg")
Required operations	 Select menu "recipe management" > select "recipe management" > SWS "recipe management" An overview of all recipes is displayed In overview window of recipes select recipe e.g. "Validation A 500 mg" and press SWS "Versions" All filters are switched on Select a recipe version with status "Disabled" Press SWS "Delete" and confirm indicated message
Consequence	Status of selected recipe version changes to "Deleted" and can not be enabled or activated
Comments	• None
Acknowledgement	• None
Test result	yes/no
Acceptance criteria	 Status of selected recipe version changes to "Deleted" Deleted recipe versions can no longer be enabled or activated
	2 Bolotod recipe versione dan no longer be enabled of delivated
Comments	
Comments	
Comments Results comply	yes/no Date/Initials

6 Batch Management

6.1 Create new batch

Test objective	Check whether new batch can be created		
Test procedure			
Test prerequisites	Recipe is enabled and active		
Required operations	 Select menu "production" > select submenu "production" Enter in window "Batch" number "123" and confirm New batch is created Press SWS "Start" to start new created batch and confirm indicated message Try to insert name and number for a new batch in window "Batch" 		
Consequence	 Creation of new batch is not possible (Number and name are deactivated fields while the previous batch is activated) 		
Required operations	 Press SWS "Complete" to complete the active batch and confirm indicated message Press SWS "X" to close window "Print batch data" and confirm indicated message Enter in window "Batch" number "1234" and confirm New batch is created 		
Consequence	A new batch can only be created if previous batch has been completed		
Comments	• None		
Acknowledgement	• None		
Test result	yes/no		
Acceptance criteria	A new batch can only be created if previous batch has been completed		
Comments			
Results comply	yes/no Date/Initials		
Results approved	Date/Initials		

6.2 Start batch

Test objective	Check whether batch can be started		
Test procedure			
Test prerequisites	A new batch has been created but not yet started		
Required operations	Press SWS "Start" and confirm indicated message		
Consequence	 Status of created and started batches changes from "Completed" to "Started" Batch information is displayed in submenu "Production", SWS "Show batch data" in window "Batch" 		
Comments	It is not possible to create or start a new batch		
Acknowledgement	• None		
Test result	yes/no		
Acceptance criteria	Only the new created batch can be started		
Comments			

Results comply	yes/no Date/Initials
Results approved	Date/Initials

6.3 Interrupt and resume batch

Test objective	Check whether batch can be interrupted and resumed again
Test procedure	
Test prerequisites	A batch has been started
Required operations	 Select menu "production" > select submenu "production" In window "Batch" press SWS "Interrupt" and confirm indicated message Press SWS "X" to close window "Print batch data"
Consequence	In overview Batch management status of active batch changes from "Started" to "Interrupted"
Required operations	 Select menu "production" > select submenu "production" In window "Batch" press SWS "Resume" and confirm indicated message
Consequence	In overview Batch management status of active batch changes from "Interrupted" to "Resumed"
Comments	• None
Acknowledgement	None
Test result	yes/no
Acceptance criteria	Batch can be interrupted and resumed
Comments	

Results comply	yes/no	Date/Initials
Results approved	Date/Initials	

6.4 Complete batch

Test objective	Check whether batch can be completed		
Test procedure			
Test prerequisites	A batch has been started or interrupted		
Required operations	 Select menu "production" > select submenu "production" In window "Batch" press SWS "Complete" and confirm indicated message Press SWS "X" to close window "Print batch data" and confirm indicated message 		
Consequence	In overview Batch management status of active batch changes to "Completed"		
Required operations	Try to resume batch		
Consequence	After batch is completed it cannot be resumed again		
Comments	• None		
Acknowledgement	None		
Test result	yes/no		
Acceptance criteria	After batch is completed, batch data can be saved or printed and batch cannot be resumed again		
Comments			

Results comply	yes/no Date/Initials
Results approved	Date/Initials

6.5 Batch history

Test objective	Check whether batch history can be displayed
Test procedure	
Test prerequisites	Batch is completed
Required operations	 Select menu "Batch history" > select SWS "Batch management" Select a batch with status "Completed" and press SWS "Show" Window "Batch view" opens and batch information such as "Batch-ID.", "Batch name", "Batch size" and "Batch part size" as well as active recipe version are displayed Press SWS "History"
Consequence	"Batch history" displays when and by whom the status of the batch has been changed
Comments	• None
Acknowledgement	• None
Test result	yes/no
Acceptance criteria	"Batch history" displays when and by whom status of batch has been changed
Comments	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

6.6 Store batch data

Test objective	Check whether batch data can be stored
Test procedure	
Test prerequisites	A batch with status "Completed" is created
Required operations	 Select menu "Batch History" > select SWS "Batch management" Select a batch with status "Completed" Batch-ID.:
	 Press SWS "Store" Select a folder and press SWS ">"" "Message ("Store") " Do you want to delete the sources after batch storage is complete?" is displayed" is displayed Press SWS "No" Confirm indicated message ("Store") Close window In "Batch History" press SWS "Batch management" Choose the batch and press SWS "Show" Press SWS "History" Note time when .zip-file has been stored
Consequence	 Batch is still listed in overview "Batch management" Status has changed from "Completed" to "Stored" Batch is stored in selected folder with file name: batch number +_Date +_Time?? + ".zip" Example: batch no.: 12345 file name: 12345_20140128_122935.zip

Test procedure	
Required operations	 Select menu "Batch management" > select submenu "Batch history" > select SWS "Batch management"
	Select a batch with status "Completed"
	Batch-ID.:
	Press SWS "Store"
	 Select a folder and press SWS " "" Message ("Store") "The storage result is good. Do you want to delete the sources?" is displayed
	Press SWS "Yes"
	Confirm indicated message ("Store")
	Close window
	Press SWS "Batch management"
Consequence	Batch is no longer listed in overview
Comments	• None
Acknowledgement	• None
Test result	yes/no
Acceptance criteria	Batches are stored as .zip-file in selected folder with correct name (+ date and time)
Comments	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

6.7 Restore batch data

Test objective	Check whether stored batch data can be restored
Test procedure	
Test prerequisites	A batch with status "Stored" is created
Required operations	 Select menu "Batch management" > select submenu "Batch history" > select SWS "Batch management" Press SWS "Restore" Select a stored file Batch-ID.: Select a folder and press SWS "✓" Message ("Store") "The target already exists. Do you really want to overwrite?" is displayed Press SWS "Yes" Select restored file in window "Batch management" and press SWS "Show" Press SWS "History"
Consequence	 .zip-file is not deleted Batch is listed in overview with status "Restored" Existing data is replaced in overview
Required operations	 Select menu "Batch management" > select submenu "Batch history" > select SWS "Batch management" Press SWS "Restore" Select same file as before and press SWS "✓" Message ("Store") "The target already exists. Do you really want to overwrite?" is displayed Press SWS "No"
Consequence	Action is interrupted
Comments	• None
Acknowledgement	• None

Test result		yes/no
Acceptance criteria	Batch is restored from .zip-file and listed in "Batch management"	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

6.8 Delete batch data

Test objective	Check whether only stored batch data can be deleted
Test procedure	
Test prerequisites	Batch with status "Stored", "Restored" and "Completed" are created
Required operations	Select menu "Batch management" > select "Batch history" > Select SWS "Batch management"
	Select batch with status "Stored"
	 Press SWS "Delete" Message ("Store") "Do you really want to delete 1 batches?" is displayed
	Press SWS "Yes"
Consequence	Batch is no longer listed in overview
Required operations	Select batch with status "Restored"
	Press SWS "Delete"
	Message ("Store") "Do you really want to delete 1 batches?" is displayed
	Press SWS "Yes"
Consequence	Batch is no longer listed in overview
Required operations	Select batch with status "Completed"
Consequence	SWS "Delete" is not active
	Batch cannot be deleted
Comments	• None
Acknowledgement	• None

Test result		yes/no
Acceptance criteria	Batches with status "Restored" can be deleted from "Batch management"	
	Batches with status "Stored" can be deleted from "Batch management"	
	Batches with status "Completed" can not be deleted from "Batch management"	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

6.9 Print history of batch data

Test objective	Check whether history of batch data can be printed
Test procedure	
Test prerequisites	Batch with status "Completed" is created
Required operations	Select menu "Batch management" > select submenu "Batch history" > select SWS "Batch management"
	Select batch with status "Completed"
	 Press SWS "History" Message ("History") "The selected batch was chosen for the history." is displayed
	Press SWS "Ok"
	Select menu "Batch history"
	Recipe data and results can be viewed for the selected batch
	Press SWS "Print batch"
Consequence	Selected batch data can be printed as hard copy on connected printer or as pdf-file to be saved on a selected drive
Comments	• None
Acknowledgement	• None
Test result	yes/no
Acceptance criteria	Recipe data and results can be viewed and printed
Comments	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

6.10 Print batch data

Test objective	Print batch data
Test procedure	
Test prerequisites	Machine is ready in automatic mode
	User "validation" is logged in
Required operations	Select menu "Batch management" > select submenu "Batch History" > select SWS "Batch management"
	Select batch with status "Completed"
	 Press SWS "History" Message ("History") "The selected batch was chosen for the history" is displayed
	Press SWS "Ok"
	Change to submenu "History"
	Result and Recipe data can be viewed for the selected batch
	Press SWS "Print batch"
	Select pdf for "Result" and "Recipe data"
	Select printer for "Result" and "Recipe data"
	In "Configuration" select target File path for storage and printer
	• Press SWS ""
Consequence	 Result and recipe data are exported to selected folder as .pdf-file with file name "batch number + _Result.pdf" resp. "batch number + _RecipeData.pdf" *)cted folder as .pdf-file with file name "batch number + _Result.pdf" *)
	Recipe data are printed on connected printer
Required operations	Select menu "Production" > select submenu "Production"
	Create a new batch
	Start the new batch and confirm indicated message
	 Complete the new batch Message "Do you really want to complete the batch?" is displayed
	Press SWS "OK"Message "Print batch data" is displayed
	Select pdf for "Alarms" and "Parameter change"
	Select printer for "Parameter change"
	In "Configuration" select target File path for storage and printer
	• Press SWS ""
Consequence	Batch alarm list and parameter change are exported to selected folder as .pdf-files with file names "batch number + _Alarm.pdf" and + "batch number + _Parameter.pdf" *)
	Batch parameter change is printed on connected printer.

C_print-Bdat(1)

Test procedure		
Comments	*) Date and time can be added in file name if needed	
Acknowledgement	• None	
Test result		yes/no
Acceptance criteria	Batch data can be printed when active batch is completed or after a batch has been selected for history	
	Batch data can be printed on connected printer	
	Batch data can be exported to a selected folder as .pdf-file	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

7 Reports and print outs

7.1 Batch print out

Test objective	Check whether batch data can be printed
Test procedure	
Test prerequisites	 Machine is ready in automatic mode User "adminlocal" is logged in
Required operations	 Select "Archive" Select menu "Batch management" Select batch with status "Completed" Press SWS "History" Message ("History") "The selected batch was chosen for the history" is displayed Press SWS "Ok" Press "Show > Press "History" Change to submenu "Batch history" > Press "Show batch data" > Press "Print" Result and Recipe data can be viewed for the selected batch Change to "Batch history" Select pdf for "Result" and "Recipe data" In "Configuration" select target File path for storage Press SWS "\nu"
Consequence	Result and recipe data are exported to selected folder as .pdf-file with file name "batch number + _Result.pdf" resp. "batch number + _RecipeData.pdf" *) * * * * * * * * * * * *
Comments	 Attach the printout to this test (can be done at the end of testing). *) Date and time can be added in file name if needed
Acknowledgement	• None

Test result		yes/no
Acceptance criteria	Batch data can be printed (as pdf-file) after a batch has been selected for history	
	Batch data can be exported	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

7.2 Parameter: Machine settings (Recipe)

Test objective	Parameter: Mechanical settings	
Test procedure		
Test prerequisites	 Machine is ready in automatic mode User "adminlocal" is logged in *) 	
Required operations	Call up "Recipe"Press "Print active recipe"	
Comments	 Attach the printout to this test (can be done at the end of testing). *)Test can be done by all users of the "Group Management" 	
Acknowledgement	• None	
Test result		yes/no
Acceptance criteria	Print out from the machine settings (Recipe) is possible	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

7.3 User right print out

Test objective	Check whether list of user rights can be printed
Test procedure	
Test prerequisites	Machine is ready in automatic modeUser "adminlocal" is logged in
Required operations	 Call up "Settings" -> page tab "User management" -> select "Group Management" Choose the user group "Admin customer" (or "Administrator") and press "Configure" -> select "all rights" Press "Print", save as a PDF-File and print it later.
Comments	Attach the printout to this test (can be done at the end of testing).
Acknowledgement	• None
Test result	yes/no
Acceptance criteria	It is possible to print or save the user rights
Comments	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

8 Basic function

8.1 Service: Clean Display

Test objective	Check of the "Clean Display" function	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
	User adminlocal is logged in	
Required operations	Call up "Settings" -> "Service" and press "Clean Display"	
	 If you press the button, a dialog window appears, where you can choose "O or "Cancel" 	k"
	• With "Ok" the display turns grey and no touch inputs can be done, and the time of 30 seconds runs down and is shown on the display	
	After 30 seconds the HMI appears again	
Comments	• None	
Acknowledgement	• None	
Test result	yes/n	0
Acceptance criteria	Display turns grey and no touch input can be done	
	After 30 seconds the HMI appears again	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

8.2 Alarmlists

Test objective	Availability of Alarmlists
Test procedure	
Test prerequisites	Machine is ready in automatic mode
Required operations	 Call up "Settings" -> "Documentation" Lists of Faults, Warnings or Messages can be selected The complete lists of Faults, Warnings and Messages can be printed or saved as pdf
Comments	• None
Acknowledgement	• None
Test result	yes/no
Test result Acceptance criteria	 Lists of the Faults, Warnings and Messages can be selected Lists of Faults, Warnings and Messages can be printed or saved as pdf

Results comply	yes/no Date/Initials
Results approved	Date/Initials

8.3 Operating hours counter

Test objective	Display of Operating hours counter
Test procedure	
Test prerequisites	Machine is ready in automatic mode
Required operations	 Go to "Settings" -> "Information" submenu "Machine" to see "Operating hours" On the display the operating hours counter is shown
Comments	• None
Acknowledgement	• None
Test result	yes/no
Acceptance criteria	Operating hours counter is shown
Comments	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

8.4 Parameter limits

Test objective	Test of limitation of Parameter limits	
Test procedure		
Test prerequisites	Machine is ready in automatic modeUser adminlocal is logged in	
Required operations	 Call up "Recipe Management" -> "Recipe Data" Select Parameter [1] in "Speeds": [1] Automatic Click on value Try to enter a value greater than the maximum limit of the machine selected a value smaller than the minimum limit of the machine selected a value within the limits 	
Comments	Limit values are shown in window	
Acknowledgement	Enter original value	
Test result		yes/no
Acceptance criteria	A value greater than the maximum limit cannot be confirmed	
	A value smaller than the minimum limit cannot be confirmed	
	A value between the minimum and maximum limit can be confirmed	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

HMI Language 8.5

Test objective	Availability and change of HMI languages
Test procedure	
Test prerequisites	Machine is ready in automatic mode
Required operations	 Call up "System configuration" -> "Language" and press "German" After pressing the button "German", the complete HMI texts change into German After pressing the button "English", the complete HMI texts change into English
	 Call up "System configuration" -> "Language" and press "French" After pressing the button "French", the complete HMI texts change into French
Comments	• None
Acknowledgement	• None
Test result	yes/no
Acceptance criteria	The HMI language switches to the chosen language
Comments	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9 Alarm Function Tests

Important Information

The following table lists the fault messages, warning messages, messages and software switches which are displayed on the control panel but can be simulated only by program changes and are therefore not tested.





WARNING

To perform a test procedure, it may be necessary to start the machine while the doors are open and the door contact switches are bridged.

Safety functions are not-active.

Severe injuries may result.

- Only authorized/trained persons are allowed to run the machine.
- Special attention must be given and care taken when working in the danger zone without protection.
- Please note that some fault messages, warnings or messages may contain a variable which is generated e.g. from the servo drive or profibus and gives further information about this device. The variable depends on the actual fault and therefore may vary.

Message number	Message text
162	Motion control: Data transmission to motion drives
164	Motion control: Synchronizing
2003	Codereader leaflet: read error
2004	Codereader carton: read error
2005	Codereader brochure: read error
2006	Codereader glued-in brochure: read error

Warning number	Warning text
	Machine: Compressed air switched off
	Control-PC: Used disk space exceeds 80%

Fault number	Fault text
8	HMI: Fault communication to control
12	Motion controller: 0 : No error code active
13	Machine: Reference not set to zero

Fault number	Fault text
15	Control-PC: Used disk space exceeds 90%
20	Batch control: Partial batch obtained
21	Batch control: Complete batch obtained
26	Control cabinet: Com error safety plc
27	Servo drives: Incorrect cam data
28	Control cabinet: Monitoring error safety plc
43	Control cabinet: Voltage monitoring
45	Control cabinet: UPS not ready
77	Machine: Fault monitoring contactors dropped away
78	Emegency stop: Safety circuit not active
82	Machine: Fault speed monitoring Inching mode - IF91_F50
83	Guard basic machine: Safety circuit not active
178	Motion drive slide in table: Fault
179	Motion drive slide in table: 0 : No error code active
186	Motion drive horizontal preinsertion: Fault
187	Motion drive vertical preinsertion: Fault
188	Motion drive insertion: Fault
191	Motion drive insertion: 0 : No error code active
192	Motion drive horizontal preinsertion: 0 : No error code active
193	Motion drive vertical preinsertion: 0 : No error code active
196	Motion drive main drive: Fault
198	Motion drive main drive: 0 : No error code active

Fault number	Fault text
285	Motion drive carton setup: Fault
288	Motion carton setup: 0 : No error code active
401	Robot 1: Motion drive: Fault
455	Robot 1: Collision range
478	Robot 1: Motion control: Fault
480	Robot 1: Motion drive axis A: Fault
481	Robot 1: Motion drive axis B: Fault
482	Robot 1: Motion drive axis C: Fault
483	Robot 1: Motion drive rotary axis: Fault
486	Robot 1: Motion drive axis A: 0 : No error code active
487	Robot 1: Motion drive axis B: 0 : No error code active
488	Robot 1: Motion drive axis C: 0 : No error code active
489	Robot 1: Motion drive rotary axis: 0 : No error code active
555	Robot 2: Collision range
580	Robot 2: Motion drive axis A: Fault
581	Robot 2: Motion drive axis B: Fault
582	Robot 2: Motion drive axis C: Fault
583	Robot 2: Motion drive rotary axis: Fault
586	Robot 2: Motion drive axis A: 0 : No error code active

Fault number	Fault text
587	Robot 1: Motion drive axis B: 0 : No error code active
588	Robot 2: Motion drive axis C: 0 : No error code active
589	Robot 2: Motion drive rotaty axis: 0 : No error code active
809	Leaflet inserter: Motion drive: Overload power supply
812	Leaflet inserter: Motion drive: Fault
819	Leaflet inserter: Motion drive: 0 : No error code active
901	Brochure infeed: Fault code reader
1238	Motion drive outfeed belt, slide in: Fault
1239	Motion drive outfeed belt, slide in: 0 : No error code active
1240	Motion drive outfeed belt, opp. slide in: Fault
1241	Motion drive outfeed belt, opp. slide in: 0 : No error code active
1433	Robot 1: Vacuum monitor: Suction cup 1
1434	Robot 1: Vacuum monitor: Suction cup 2
1435	Robot 1: Vacuum monitor: Suction cup 3
1443	Robot 1: Control cabinet: Voltage monitoring
1533	Robot 2: Vacuum monitor: Suction cup 1
1534	Robot 2: Vacuum monitor: Suction cup 2
1535	Robot 2: Vacuum monitor: Suction cup 3

Fault number	Fault text
1543	Robot 2: Control cabinet: Voltage monitoring
2006	Robot 2: Frequency conveyer puck outfeed belt
2007	Robot 2: Frequency conveyer puck infeed belt
2008	Robot 2: Frequency conveyer diskus belt
2028	Leaflet inserter: Motion drive leaflet reject wheel: 0 : No error activated
2030	Motion drive discharge belt front: @Fault
2031	Motion drive discharge belt front: 0 : No error code active
2032	Motion drive discharge belt back: @Fault
2033	Motion drive discharge belt back: 0 : No error code active
2034	Motion drive discharge belt top: @Fault
2035	Motion drive discharge belt top: 0 : No error code active
2047	Code reader carton: no Code loaded
2048	Code reader leaflet: no Code loaded
2049	Code reader brochure: no Code loaded
2050	Code reader glued-in brochure: no Code loaded
2051	Robot 1: Motion drive infeed belt master: Fault
2052	Robot 1: Motion drive infeed belt master: 0 : No error code active
2053	Robot 1: Motion drive infeed belt ahead: Fault
2054	Robot 1: Motion drive infeed belt ahead: 0 : No error code active
2063	Infeed 2: Monitoring sensor: Reducer accumulation pressure 1 -W166_B64
2064	Infeed 2: Monitoring sensor: Reducer accumulation pressure 1 -W166_B65
2065	Infeed 2: Monitoring sensor: Reducer accumulation pressure 2 -W166_B67
2066	Infeed 2: Monitoring sensor: Reducer accumulation pressure 2 -W166_B68
2067	Number of requested tara parts reached

Test protocols - Function tests software switches operation mode (SWSOPM) 9.1

SWSOPM 1: AUTOMATIC 9.1.1

Test objective	Function test of operating mode switch (SWSOPM)		
Test procedure			
	TEST 1		
Test prerequisites	Machine is ready in setup mode		
Required operations	 Activate SWSOPM 1 "Automatic" Press "Reset" Press "Start" 		
Consequence	Message is displayed on control panel (ME40 "Operating mode: Automatic")		
Comments	Test can be done together with test ME 40		
Acknowledgement	Press "Stop"Press "Reset"		
	TEST 2		
Test prerequisites	Machine is ready in automatic mode		
Required operations	 Activate SWSOPM 2 "Setup" Press "Reset" Press "Start" 		
Consequence	Message is displayed on control panel (ME41 "Operating mode: Setup")		
Comments	Test can be done together with test ME 41		
Acknowledgement	Press "Stop"Press "Reset"		

Test result		yes/no
Acceptance criteria	Test 1: SWSOPM 1 "Automatic" is activated, Operating mode "Automatic" is active	
	 Test 1: SWSOPM 1 "Automatic" is activated, Message is displayed on control panel (ME 40) 	
	Test 2: SWSOPM 2 "Setup" is activated, Operating mode "Automatic" is not active	
	Test 2: SWSOPM 2 "Setup" is activated, Message is displayed on control panel (ME41)	
Comments		
)

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWSOPM 2: SETUP 9.1.2

Test objective	Function test of operating mode switch (SWSOPM)
Test procedure	
	TEST 1
Test prerequisites	Machine is ready in automatic mode
Required operations	Activate SWSOPM 2 "Setup"
	Press "Reset"
	Press "Start"
Consequence	Message is displayed on control panel (ME41 "Operating mode: Setup")
Comments	Test can be done together with test of ME 41
Acknowledgement	Press "Stop"
	Press "Reset"
	TEST 2
Test prerequisites	Machine is ready in setup mode
Required operations	Activate SWSOPM 1 "Automatic"
	Press "Reset"
	Press "Start"
Consequence	Message is displayed on control panel (ME40 "Operating mode: Automatic")
Comments	Test can be done together with test of ME 41
Acknowledgement	Press "Stop"
	Press "Reset"

Test result		yes/no
Acceptance criteria	Test 1: SWSOPM 2 "Setup" is activated, Operating mode "Setup" is active	
	Test 1: SWSOPM 2 "Setup" is activated, Message is displayed on control panel (ME 41)	
	Test 2: SWSOPM 1 "Automatic" is activated, Operating mode "Setup" is not active	
	 Test 2: SWSOPM 1 "Automatic" is activated, Message is displayed on control panel (ME 40) 	
Comments		
)

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWSOPM 3: INCHING PANEL 9.1.3

Test objective	Function test of operating mode switch (SWSOPM)
Test procedure	
	TEST 1
Test prerequisites	Machine is ready in automatic mode
Required operations	 Activate SWSOPM 3 "Inching panel" Press "Reset" Press "Start" and hold the button "Start"
Consequence	 Message is displayed on control panel (ME42 "Operating mode: Inching mode panel") Operating mode "Inching panel" is active
Comments	Test can be done together with test of ME 42
Acknowledgement	Press "Stop"Press "Reset"
	TEST 2
Test prerequisites	Machine is ready in inching mode panel
Required operations	 Activate SWSOPM 1 "Automatic" Press "Reset" Press "Start"
Consequence	 Message is displayed on control panel (ME40 "Operating mode: Automatic") Operating mode "Inching panel" is not active
Comments	Test can be done together with test of ME 40
Acknowledgement	Press "Stop"Press "Reset"

Test result		yes/no
Acceptance criteria	Test 1: SWSOPM 3 "Inching panel" is activated, Operating mode "Inching panel" is active	
	 Test 1: SWSOPM 3 "Inching panel" is activated, Message is displayed on control panel (ME 42) 	
	Test 2: SWSOPM 1 "Automatic" is activated, Operating mode "Inching panel" is not active	
	Test2: SWSOPM 1 "Automatic" is activated, Message is displayed on control panel (ME 40)	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWSOPM 4: INCHING CABLE 9.1.4

Test objective	Function test of operating mode switch (SWSOPM)
Test procedure	
	TEST 1
Test prerequisites	Machine is ready in automatic mode
Required operations	 Plug in inching cable Activate SWSOPM 4 "Inching cable" Run the machine with inching cable
Consequence	 Message is displayed on control panel (ME 43 "Operating mode: Inching mode cable front") Operating mode "Inching cable" is active
Comments	Test can be done together with test of ME 43
Acknowledgement	• None
	TEST 2
Test prerequisites	Machine is ready in inching mode cable
Required operations	 Remove inching cable Press "Reset" Activate SWSOPM 1 "Automatic" Press "Reset" Press "Start"
Consequence	 Message is displayed on control panel (ME 40 "Operating mode: Automatic") Operating mode "Inching cable" is not active
Comments	Test can be done together with test of ME 43
Acknowledgement	Press "Stop"Press "Reset"

Test result		yes/no
Acceptance criteria	Test 1: SWSOPM 4 "Inching cable" is activated, Operating mode "Inching cable" is active	
	 Test 1: SWSOPM 4 "Inching cable" is activated, Message is displayed on control panel (ME 43) 	
	Test 2: SWSOPM 1 "Automatic" is activated, Operating mode "Inching cable" is not active	
	 Test 2: SWSOPM 1 "Automatic" is activated, Message is displayed on control panel (ME 40) 	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWSOPM 6: OPEN STOP BRAKE 9.1.5

Test objective	Function test of operating mode switch (SWSOPM)	
Test procedure		
Test prerequisites	Machine is ready in inching mode	
Required operations	Activate SWSOPM 6 "Open stop brake"Select one of the drives	
Consequence	Chosen drive can be moved manuallyWarning is displayed on control panel	
Comments	• None	
Acknowledgement	• None	
Test result	yes/no)
Acceptance criteria	SWSOPM 6 "Open stop brake" is activated, chosen drive can be moved manually	
	SWSOPM 6 "Open stop brake" is activated, Warning is displayed on control panel	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWSOPM 7: ZERO DRIVE 9.1.6

Test objective	Function test of operating mode switch (SWSOPM)	
Test procedure		
Test prerequisites	Machine is ready in setup modeServo drives are set to zero position	
Required operations	 Activate SWSOPM 7 "Zero drive" Choose one of the drives Confirm with the checkmark 	
Consequence	Drive is zeroed	
Comments	• None	
Acknowledgement	None	
Test result	yes	/no
Acceptance criteria	SWSOPM 7 "Zero drive" is activated, Drive is zeroed	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWSOPM0007_AFT_CUC_CUK(2)

SWSOPM 8: INCHING SINGLE AXIS 9.1.7

Test objective	Function test of operating mode switch (SWSOPM)
Test procedure	
	TEST 1
Test prerequisites	Machine is ready in inching mode
Required operations	 Activate SWSOPM 8 "Inching single axis" Select one of the robot axis Press "Start" and hold the button "Start"
Consequence	 Message is displayed on control panel (ME47: "Operating mode: Inching mode single axis") Operating mode "Inching single axis" is active
	TEST 2
Test prerequisites	Machine is ready in inching mode single axis
Required operations	 Activate SWSOPM 1 "Automatic" Press "Reset" Press "Start"
Consequence	 Message is displayed on control panel (ME40 "Operating mode: Automatic") Operating mode "Inching single axis" is not active
Comments	• None
Acknowledgement	Press "Stop"Press "Reset"

	yes/no
Test 1: Operating mode "Inching single axis" is active	
Test 1: Message is displayed on control panel (ME 47)	
Test 2: Operating mode "Inching single axis" is not active	
Test 2: Message is displayed on control panel (ME 40)	
	 Test 1: Message is displayed on control panel (ME 47) Test 2: Operating mode "Inching single axis" is not active

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWSOPM 11: INCHING INFEED 9.1.8

Test objective	Function test of operating mode switch (SWSOPM)
Test procedure	
Test prerequisites	Machine is ready in automatic mode
Required operations	 Activate SWSOPM 11 "Inching infeed" Press "Reset" Press "Start" at the operating panel of the infeed system
Consequence	 Message is displayed on control panel (ME49 "Operating mode: Inching mode infeed") Operating mode "Inching infeed" is active
Comments	Test can be done together with test of ME 49
Acknowledgement	Press "Stop"Press "Reset"
	TEST 2
Test prerequisites	Machine is ready in inching mode panel
Required operations	 Activate SWSOPM 1 "Automatic" Press "Reset" Press "Start"
Consequence	 Message is displayed on control panel (ME40 "Operating mode: Automatic") Operating mode "Inching infeed" is not active
Comments	Test can be done together with test of ME 40
Acknowledgement	Press "Stop"Press "Reset"

Test result		yes/no
Acceptance criteria	Test 1: SWSOPM 11 "Inching infeed" is activated, Operating mode "Inching panel" is active	
	 Test 1: SWSOPM 11 "Inching infeed" is activated, Message is displayed on control panel (ME 49) 	
	Test 2: SWSOPM 1 "Automatic" is activated, Operating mode "Inching infeed" is not active	
	Test2: SWSOPM 1 "Automatic" is activated, Message is displayed on control panel (ME 40)	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWSOPM 60: LINE MODE DOWNSTREAM EQUIPMENT 9.1.9

Test objective	Function test of operating mode switch (SWSOPM)
Test procedure	
	TEST 1
Test prerequisites	Machine is ready in automatic mode
Required operations	 Activate SWSOPM 60 "Line mode downstream equipment" Press "Start" Create fault at downstream machine
Consequence	Machine stopsFault message is displayed on control panel
Comments	• None
Acknowledgement	Press "Stop"Press "Reset"
	TEST 2
Test prerequisites	Machine is ready in automatic mode
Required operations	 Deactivate SWSOPM 60 "Line mode downstream equipment" Press "Start" Try to create fault at downstream machine
Consequence	 Machine is running in automatic mode Fault message is not displayed in control panel
Comments	• None
Acknowledgement	Press "Stop"Press "Reset"

Test result		yes/no
Acceptance criteria	Test 1: SWSOPM 60 "Line mode downstream equipment" is activated: Machine is in line mode. In case of faults fault message will be displayed on control panel and machine stops.	
	 Test 2: SWSOPM 60 "Line mode downstream equipment" is deactivated: Machine is not in line mode. In case of faults fault message is not displayed on control panel, machine continues run in automatic mode. 	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWSOPM 62: MACHINE: AUTOMATIC RESTART 9.1.10

Test objective	Function test of operating mode switch (SWSOPM)
Test procedure	
	TEST 1
Test prerequisites	 SWSOPM 60 "Line mode downstream equipment" is activated SWSOPM 62 "Machine: Automatic restart" is activated Line is running at reduced speed
Required operations	Create FLT 2041 "Robot 2:Emergency stopp activated: Diskus belt"
Consequence	Fault message is displayed on control panel
Required operations	Press "Reset"
Consequence	Machine restarts automatically after acknowledgement of FLT 2041
Comments	Test can be done together with test of FLT 2041
Acknowledgement	Press "Stop"
	TEST 2
Test prerequisites	 Machine is ready in automatic mode SWSOPM 62 "Machine: Automatic restart" is deactivated
Required operations	 Press "Start" Create FLT 2041 "Robot 2:Emergency stopp activated: Diskus belt"
Consequence	Machine stopsFault message is displayed on control panel
Required operations	Acknowledge FLT 2041Press "Reset"
Consequence	Machine does not restarts automatically after acknowledgement of FLT 2041
Comments	• None
Acknowledgement	• None

Test result		yes/no
Acceptance criteria	Test 1: SWSOPM 62 "Machine: Automatic restart" is activated, Machine restarts automatically after acknowledgement of the fault message	
	 Test 2: SWSOPM 62 "Machine: Automatic restart" is deactivated, Machine does not restart automatically after acknowledgement of the fault message 	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWSOPM 63: LINE MODE UPSTREAM EQUIPMENT 9.1.11

Test objective	Function test of operating mode switch (SWSOPM)
Test procedure	
	TEST 1
Test prerequisites	Machine is ready in automatic mode
Required operations	Activate SWSOPM 63 "Line mode upstream equipment"
	Activate SWS 400 "Robot 1"
	Activate SWS 500 "Robot 2"
	Press "Start"
Consequence	Machine is running in line mode
Comments	• None
Acknowledgement	• None
	TEST 2
Test prerequisites	Machine is ready in automatic mode
Required operations	Deactivate SWSOPM 63 "Line mode upstream equipment"
	Deactivate SWS 400 "Robot 1"
	Deactivate SWS 500 "Robot 2"
	Press "Start"
Consequence	Machine is running without product infeed
Comments	• None
Acknowledgement	• None

Test result		yes/no
Acceptance criteria	Test 1: SWSOPM 63 "Line mode upstream equipment" is activated, Machine is running with product infeed	
	Test 2: SWSOPM 63 "Line mode upstream equipment" is deactivated, Machine is running without product infeed	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

Test protocols - Function tests software switches operation mode drive 9.2 (SWSOPM_DRV)

SWSOPM_DRV 1002: INSERTION 9.2.1

Test objective	Function test of operating mode drive (SWSOPM_DRV)
Test procedure	
Test prerequisites	Machine is ready in automatic mode
Required operations	 Set drive "Insertion" mechanically to position "zero" Close guard door: Press "Reset" Activate SWSOPM 7 "Zero drive" Activate SWSOPM_DRV 1002
Consequence	Message is displayed on control panel "Do you really want to set drive to zero"
Comments	The position of the drive should not be adjusted, as it has already been adjusted.
Acknowledgement	 Message "Do you really want to set drive to zero" - Cancel input Press "Reset"
Test result	yes/no
Acceptance criteria	SWSOPM_DRV 1002 "Insertion" is activated, Message is displayed on control panel "Do you really want to set drive to zero"
Comments	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWSOPM_DRV 1003: PREINSERTION HORIZONTAL 9.2.2

Test objective	Function test of operating mode drive (SWSOPM_DRV)	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Activate SWSOPM 6 "Open stop brake"Activate SWSOPM_DRV 1003	
Consequence	Brake can be moved manually	
Comments	• None	
Acknowledgement	None	
Test result	у	es/no
Acceptance criteria	Brake is released/ open	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWSOPM_DRV 1004: PREINSERTION VERTICAL 9.2.3

Test objective	Function test of operating mode drive (SWSOPM_DRV)	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Activate SWSOPM 6 "Open stop brake"Activate SWSOPM_DRV 1004	
Consequence	Brake can be moved manually	
Comments	• None	
Acknowledgement	None	
Test result		yes/no
Acceptance criteria	Brake is released/ open	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWSOPM_DRV 1007: CARTON SETUP 9.2.4

Test objective	Function test of operating mode drive (SWSOPM_DRV)	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Activate SWSOPM 6 "Open stop brake"Activate SWSOPM_DRV 1007	
Consequence	Brake can be moved manually	
Comments	• None	
Acknowledgement	• None	
Test result		yes/no
Acceptance criteria	Brake is released/ open	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWSOPM_DRV 1008: LEAFLET DEVICE 9.2.5

Test objective	Function test of operating mode drive (SWSOPM_DRV)
Test procedure	
Test prerequisites	Machine is ready in automatic mode
Required operations	 Set drive "Insertion" mechanically to position "zero" Close guard door: Press "Reset" Activate SWSOPM 7 "Zero drive" Activate SWSOPM_DRV 1008
Consequence	Message is displayed on control panel "Do you really want to set drive to zero"
Comments	The position of the drive should not be adjusted, as it has already been adjusted.
Acknowledgement	 Message "Do you really want to set drive to zero" - Cancel input Press "Reset"
Test result	yes/no
Acceptance criteria Comments	SWSOPM_DRV 1008 "Leaflet device" is activated, Message is displayed on control panel "Do you really want to set drive to zero"

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWSOPM_DRV 1015: MAIN DRIVE 9.2.6

Test objective	Function test of operating mode drive (SWSOPM_DRV)
Test procedure	
Test prerequisites	Machine is ready in setup mode
Required operations	 Set drive mechanically to position "zero" Close guard door: Press "Reset" Activate SWSOPM 7 "Zero drive" Activate SWSOPM_DRV 1015 "Main drive"
Consequence	Message is displayed on control panel "Do you really want to set drive to zero"
Comments	The position of the drive should not be adjusted, as it has already been adjusted.
Acknowledgement	 Message "Do you really want to set drive to zero" - Cancel input Press "Reset"
Test result	yes/no
Acceptance criteria	SWSOPM_DRV 1015 "Main drive" is activated, Message is displayed on control panel "Do you really want to set drive to zero"
Comments	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWSOPM_DRV 1016: INSERTION TABLE 9.2.7

Test objective	Function test of operating mode drive (SWSOPM_DRV)
Test procedure	
Test prerequisites	Machine is ready in automatic mode
Required operations	 Set drive "Insertion" mechanically to position "zero" Close guard door: Press "Reset" Activate SWSOPM 7 "Zero drive" Activate SWSOPM_DRV 1016
Consequence	Message is displayed on control panel "Do you really want to set drive to zero"
Comments	The position of the drive should not be adjusted, as it has already been adjusted.
Acknowledgement	 Message "Do you really want to set drive to zero" - Cancel input Press "Reset"
Test result	yes/no
Acceptance criteria Comments	SWSOPM_DRV 1016 "Insertion table" is activated, Message is displayed on control panel "Do you really want to set drive to zero"

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWSOPM_DRV 1017: CARTON ALIGNMENT 9.2.8

Test objective	Function test of operating mode drive (SWSOPM_DRV)
Test procedure	
Test prerequisites	Machine is ready in automatic mode
Required operations	 Set drive "Insertion" mechanically to position "zero" Close guard door: Press "Reset" Activate SWSOPM 7 "Zero drive" Activate SWSOPM_DRV 1017
Consequence	Message is displayed on control panel "Do you really want to set drive to zero"
Comments	The position of the drive should not be adjusted, as it has already been adjusted.
Acknowledgement	 Message "Do you really want to set drive to zero" - Cancel input Press "Reset"
Test result	yes/no
Acceptance criteria Comments	SWSOPM_DRV 1017 "Carton alignment" is activated, Message is displayed on control panel "Do you really want to set drive to zero"

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWSOPM_DRV 1054: ROBOT 1: DRILL AXIS 9.2.9

Test objective	Function test of operating mode drive (SWSOPM_DRV)	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Activate SWSOPM 6 "Open stop brake"Activate SWSOPM_DRV 1054	
Consequence	Brake can be moved manually	
Comments	• None	
Acknowledgement	None	
Test result		yes/no
Acceptance criteria	Brake is released/ open	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.2.10 SWSOPM_DRV 1055: ROBOT 1

Test objective	Function test of operating mode drive (SWSOPM_DRV)	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Press button "CAR1-W150-S20P"	
Consequence	Brake can be moved manually	
Comments	• None	
Acknowledgement	None	
Test result		yes/no
Acceptance criteria	Brake is released/ open	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWSOPM_DRV 1059: ROBOT 2: DRILL AXIS 9.2.11

Test objective	Function test of operating mode drive (SWSOPM_DRV)	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Activate SWSOPM 6 "Open stop brake"Activate SWSOPM_DRV 1059	
Consequence	Brake can be moved manually	
Comments	None	
Acknowledgement	None	
Test result		yes/no
Acceptance criteria	Brake is released/ open	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.2.12 SWSOPM_DRV 1060: ROBOT 2

Test objective	Function test of operating mode drive (SWSOPM_DRV)	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Press button "CAR1-W250-S20P	
Consequence	Brake can be moved manually	
Comments	• None	
Acknowledgement	None	
Test result		yes/no
Acceptance criteria	Brake is released/ open	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWSOPM_DRV 1061: ROBOT 1: CYCLE BELT 1 9.2.13

Test objective	Function test of operating mode drive (SWSOPM_DRV)	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Activate SWSOPM 6 "Open stop brake"Activate SWSOPM_DRV 1061	
Consequence	Brake can be moved manually	
Comments	• None	
Acknowledgement	None	
Test result		yes/no
Acceptance criteria	Brake is released/ open	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWSOPM_DRV 1062: ROBOT 1: CYCLE BELT 2 9.2.14

Test objective	Function test of operating mode drive (SWSOPM_DRV)
Test procedure	
Test prerequisites	Machine is ready in automatic mode
Required operations	Activate SWSOPM 6 "Open stop brake"Activate SWSOPM_DRV 1062
Consequence	Brake can be moved manually
Comments	• None
Acknowledgement	• None
Test result	yes/no
Acceptance criteria	Brake is released/ open
Comments	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

Test protocols - Function tests software switches (SWS) 9.3

9.3.1 SWS 52: COMPRESSED AIR

Test objective	Function test of software switch (SWS)
Test procedure	
	TEST 1
Test prerequisites	Machine is ready in setup modeSWS 400 "Robot 1" is activated
Required operations	 Activate SWS 405 "Robot 1: Call product" Activate SWS 52
Consequence	Robot 1 is running
Comments	• None
Acknowledgement	• None
	TEST 2
Test prerequisites	Machine is ready in setup modeRobot 1 is running
Required operations	Deactivate SWS 52
Consequence	Robot 1 stops running
Comments	• None
Acknowledgement	• None

Test result		yes/no
Acceptance criteria	Test 1: If SWS 52 is activated, the robot 1 is running	
	Test 2: If SWS 52 is deactivated, the robot 1 stops running	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWS 53: PERMANENT READING CODE 9.3.2

Test objective	Function test of software switch (SWS)
Test procedure	
	TEST 1
Test prerequisites	Machine is ready in setup mode
Required operations	 Activate SWS 53 "Permanent reading code" Activate button for network navigator keyence and log in Within the menu select the active format Open "Production" Select "code reader, sensor carton" In navigator keyence setup select "diagnose" Place carton in front of reading head of carton
Consequence	 Code is displayed on control panel within the "Keyence navigator-menu" Code reading is active at standstill TEST 2
Toot proroquiaitos	
Required operations	 Code reading is active at standstill Deactivate SWS 53 "Permanent reading code"
Consequence	 Code is not displayed on control panel within the "Keyence navigator-menu" Code reading is not active at standstill
Comments	During continuous code reading a lamp at the code reader carton is blinking
Acknowledgement	Deactivate SWS 53 "Permanent reading code"

Test result Acceptance criteria Test 1: Code is displayed on control panel within the "Keyence navigator-menu"	0
Test 1: Code reading is active at standstill	
Test 2: Code is not displayed on control panel within the "Keyence navigator-menu"	
Test 2: Code reading is not active at standstill	
Comments	

Results comply	yes/no	Date/Initials
Results approved	Date/Initials	

SWS 54: VACUUM 9.3.3

Test objective	Function test of software switch (SWS)
Test procedure	
	TEST 1
Test prerequisites	Machine is ready in setup mode
Required operations	 Activate SWS 208 "Continuous call carton" Activate SWS 54 "Vacuum"
Consequence	Vacuum pump is running
Comments	• None
Acknowledgement	• None
	TEST 2
Test prerequisites	Machine is ready in setup modeVacuum pump is running
Required operations	Deactivate SWS 54 "Vacuum"
Consequence	Vacuum pump is not running
Comments	• None
Acknowledgement	• None

Test result		yes/no
Acceptance criteria	Test 1: If SWS 54 "Vacuum" is activated, the vacuum pump is running	
	Test 2: If SWS 54 "Vacuum" is deactivated, the vacuum pump is not running	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWS 56: MACHINE INSIDE LIGHT OFF 9.3.4

Test objective	Function test of software switch (SWS)	
Test procedure		
	TEST 1	
Test prerequisites	Machine is ready in automatic mode	
Required operations	Activate SWS 56 "Machine inside light off"	
Consequence	Lamps in production room don't shine	
Comments	• None	
Acknowledgement	• None	
	TEST 2	
Test prerequisites	Machine is ready in automatic mode	
Required operations	Deactivate SWS 56 "Machine inside light off"	
Consequence	Lamps in production room shine	
Comments	• None	
Acknowledgement	• None	
Test result		yes/no
Test result Acceptance criteria	Test 1: SWS 56 "Machine inside light off" is activated: Lamps in production room don't shine	yes/no
	Test 1: SWS 56 "Machine inside light off" is activated: Lamps in	yes/no
	 Test 1: SWS 56 "Machine inside light off" is activated: Lamps in production room don't shine Test 2: SWS 56 "Machine inside light off" is deactivated: Lamps in 	yes/no
Acceptance criteria	 Test 1: SWS 56 "Machine inside light off" is activated: Lamps in production room don't shine Test 2: SWS 56 "Machine inside light off" is deactivated: Lamps in 	yes/no
Acceptance criteria	 Test 1: SWS 56 "Machine inside light off" is activated: Lamps in production room don't shine Test 2: SWS 56 "Machine inside light off" is deactivated: Lamps in 	yes/no

9.3.5 SWS 100: INSERT DEFECTIVE PRODUCT

Test objective	Function test of software switch (SWS)
Test procedure	
	TEST 1
Test prerequisites	Machine is ready in automatic mode
Required operations	 Activate SWS 100 "Insert defective product" Insert leaflet with faulty code in stack Insert some product in the infeed Press "Start" Call some products
Consequence	Carton with faulty leaflet is loaded but ejected on discharge belt
Comments	• None
Acknowledgement	Press "Stop"Press "Reset"
	TEST 2
Test prerequisites	Machine is ready in automatic mode
Required operations	 Deactivate SWS 100 "Insert defective product" Insert leaflet with faulty code in stack Insert some product in the infeed Press "Start" Call some products
Consequence	 Leaflet is ejected No carton call Product is ejected at insertion
Comments	• None
Acknowledgement	Press "Stop"Press "Reset"

Test result		yes/no
Acceptance criteria	Test 1: SWS 100 is activated: In case of defective product or defective leaflet, carton will be loaded but rejected on discharge belt	
	Test 2: SWS 100 is deactivated: In case of defective product or defective leaflet, product will be rejected without loading	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWS 106: CONTINUOUS INSERTION 9.3.6

Test objective	Function test of software switch (SWS)
Test procedure	
	TEST 1
Test prerequisites	Machine is ready in setup mode
Required operations	Activate SWS 106 "Continuous insertion"Press "Start"
Consequence	Insertion pushers are activated
Acknowledgement	Press "Stop"Press "Reset"
	TEST 2
Test prerequisites	Machine is ready in setup mode
Required operations	Deactivate SWS 106 "Continuous insertion"Press "Start"
Consequence	Insertion pushers are deactivated
Comments	• None
Acknowledgement	Press "Stop"Press "Reset"

Test result		yes/no
Acceptance criteria	Test 1: Insertion pushers are activated	
	Test 2: Insertion pushers are deactivated	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWS 119: CONTINUOUS PREINSERTION 9.3.7

Test objective	Function test of software switch (SWS)	
Test procedure		
	TEST 1	
Test prerequisites	Machine is ready in setup mode	
Required operations	Activate SWS 119 "Continuous preinsertion"Press "Start"	
Consequence	Preinsertion moves in continuous operation	
Comments	• None	
Acknowledgement	• None	
	TEST 2	
Test prerequisites	Machine is ready in setup mode	
Required operations	Deactivate SWS 119 "Continuous preinsertion"Press "Start"	
Consequence	Preinsertion does not move	
Comments	• None	
Acknowledgement	None	
Test result		yes/no
Acceptance criteria	Test 1: Preinsertion moves in continuous operation	
	Test 2: Preinsertion does not move	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.3.8 SWS 120: PREINSERTION UPPER FRONT POSITION

Test objective	Function test of software switch (SWS)	
Testdurchführung		
	TEST 1	
Test prerequisites	Machine is ready in setup modeGuard doors are closed	
Required operations	Activate SWS 120 "Preinsertion upper front position"Press "Reset"	
Consequence	Preinserter moves into an upper front position (assembly position)	
Comments	• None	
Acknowledgement	• None	
	TEST 2	
Test prerequisites	 Machine is ready in setup mode Guard doors are closed Preinserter is in upper front position 	
Required operations	Deactivate SWS 120 "Preinsertion upper front position"Press "Reset"	
Consequence	Preinserter moves back to original position	
Comments	• None	
Acknowledgement	• None	
Test result		yes/no
Acceptance criteria	Test 1: Preinserter moves into an upper front position (assembly position)	
	Test 2: Preinserter moves back to original position	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.3.9 SWS 121: PREINSERTION PLATE 1 DOWN

Test objective	Function test of software switch (SWS)
Test procedure	
	TEST 1
Test prerequisites	Machine is ready in setup modePreinsertion is located outside of folding carton
Required operations	 Activate SWS 121 "Preinsertion plate 1 down" Press "Reset"
Consequence	Plate 1 at preinsertion lowers
Comments	Plate 1 at preinsertion lowers only if it is located outside folding carton
Acknowledgement	• None
	TEST 2
Test prerequisites	Machine is ready in setup modePreinsertion is located outside of folding carton
Required operations	Deactivate SWS 121 "Preinsertion plate 1 down"
Consequence	Plate 1 at preinsertion raises
Comments	Plate 1 at preinsertion lowers only if it is located outside folding carton
Acknowledgement	• None
Test result	yes/no
Acceptance criteria	Test 1: If SWS 121 is activated plate 1 at preinsertion lowers
	Test 2: If SWS 121 is deactivated plate 1 at preinsertion raises
Comments	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWS 122: PREINSERTION PLATE 2 DOWN 9.3.10

Test objective	Function test of software switch (SWS)
Test procedure	
	TEST 1
Test prerequisites	 Machine is ready in setup mode Preinsertion is located outside of folding carton
Required operations	Activate SWS 122 "Preinsertion plate 2 down"Press "Reset"
Consequence	Plate 2 at preinsertion lowers
Comments	Plate 2 at preinsertion lowers only if it is located outside folding carton
Acknowledgement	• None
	TEST 2
Test prerequisites	 Machine is ready in setup mode Preinsertion is located outside of folding carton
Required operations	Deactivate SWS 122 "Preinsertion plate 2 down"
Consequence	Plate 2 at preinsertion raises
Comments	Plate 2 at preinsertion lowers only if it is located outside folding carton
Acknowledgement	• None
Test result	yes/no
Acceptance criteria	Test 1: If SWS 122 is activated plate 2 at preinsertion lowers
	Test 2: If SWS 122 is deactivated plate 2 at preinsertion raises
Comments	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.3.11 SWS 123: PREINSERTION PLATE 3 DOWN

Test objective	Function test of software switch (SWS)
Test procedure	
	TEST 1
Test prerequisites	 Machine is ready in setup mode Preinsertion is located outside of folding carton
Required operations	Activate SWS 123 "Preinsertion plate 3 down"Press "Reset"
Consequence	Plate 3 at preinsertion lowers
Comments	Plate 3 at preinsertion lowers only if it is located outside folding carton
Acknowledgement	• None
	TEST 2
Test prerequisites	 Machine is ready in setup mode Preinsertion is located outside of folding carton
Required operations	Deactivate SWS 123 "Preinsertion plate 3 down"
Consequence	Plate 3 at preinsertion raises
Comments	Plate 3 at preinsertion lowers only if it is located outside folding carton
Acknowledgement	• None
Test result	yes/no
Acceptance criteria	Test 1: If SWS 123 is activated plate 3 at preinsertion lowers
	Test 2: If SWS 123 is deactivated plate 3 at preinsertion raises
Comments	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWS 124: PREINSERTION PLATE 4 DOWN 9.3.12

Test objective	Function test of software switch (SWS)
Test procedure	
	TEST 1
Test prerequisites	 Machine is ready in setup mode Preinsertion is located outside of folding carton
Required operations	Activate SWS 124 "Preinsertion plate 4 down"Press "Reset"
Consequence	Plate 4 at preinsertion lowers
Comments	Plate 4 at preinsertion lowers only if it is located outside folding carton
Acknowledgement	• None
	TEST 2
Test prerequisites	 Machine is ready in setup mode Preinsertion is located outside of folding carton
Required operations	Deactivate SWS 124 "Preinsertion plate 4 down"
Consequence	Plate 4 at preinsertion raises
Comments	Plate 4 at preinsertion lowers only if it is located outside folding carton
Acknowledgement	• None
Test result	yes/no
Acceptance criteria	Test 1: If SWS 124 is activated plate 4 at preinsertion lowers
	Test 2: If SWS 124 is deactivated plate 4 at preinsertion raises
Comments	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.3.13 SWS 127: PREINSERTION

Test objective	Function test of software switch (SWS)	
Test procedure		
	TEST 1	
Test prerequisites	Machine is ready in setup mode	
Required operations	Activate SWS 127 "Preinsertion"Press "Start"	
Consequence	Preinsertion is active	
Comments	• None	
Acknowledgement	Press "Stop"	
	TEST 2	
Test prerequisites	Machine is ready in setup mode	
Required operations	Deactivate SWS 127 "Preinsertion"Press "Start"	
Consequence	Preinsertion is not active	
Comments	• None	
Acknowledgement	Press "Stop"Press "Reset"	
Test result	yes/r	10
Acceptance criteria	Test 1: If SWS 127 "Preinsertion" is activated, preinsertion is active	
	Test 2: If SWS 127 "Preinsertion" is deactivated, preinsertion is not active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWS 200: CODE READER CARTON 9.3.14

Test objective	Function test of software switch (SWS)
Test procedure	
	TEST 1
Test prerequisites	 Machine is ready in automatic mode Prepare two (2) cartons with wrong / manipulated carton code Product is available
Required operations	 Insert cartons with wrong / manipulated code into magazine Activate SWS 200 "Code reader carton" Press "Start"
Consequence	Wrong or non-coded cartons are ejected
	TEST 2
Required operations	 Insert cartons with wrong / manipulated code into magazine Deactivate SWS 200 "Code reader carton" Press "Start"
Consequence	Wrong or non-coded cartons are not recognized and ejected
Comments	• None
Acknowledgement	 Activate SWS 200 "Code reader carton" Press "Stop" Press "Reset"

Test result		yes/no
Acceptance criteria	 Test 1: SWS200 is activated: wrong or non-coded cartons are ejected Test 2: SWS200 is deactivated: wrong or non-coded cartons are not recognized and ejected 	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWS 201: CARTON LASER PRINTING 9.3.15

Test objective	Function test of software switch (SWS)
Test procedure	
	TEST 1
Test prerequisites	Machine is ready in automatic mode
Required operations	 Activate SWSOPM 63 "Line mode upstream equipment" Activate SWS 201 "Carton laser printer" Press "Start"
Consequence	Machine is running in line modeLaser printer unit is running
Comments	• None
Acknowledgement	Press "Stop"Press "Reset"
	TEST 2
Test prerequisites	 Machine is ready in automatic mode SWSOPM 63 "Line mode upstream equipment" is activated SWS 201 "Carton laser printer" is activated
Required operations	Press "Start"Deactivate SWS 201 "Carton laser printer"
Consequence	 Carton laser unit stops running Machine stops Fault message is displayed on control panel
Comments	• None
Acknowledgement	Press "Stop"Press "Reset"

Test result		yes/no
Acceptance criteria	Test 1: If SWS 201 "Carton laser printer" is activated: Machine is in line mode and laser printer unit is running.	
	 Test 2: If SWS 201 "Carton laser printer" is deactivated: Machine is in line mode but carton laser unit is not running. In case of faults fault message will be displayed on control panel and machine stops. 	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWS 203: BLOWING AIR CARTON 9.3.16

Test objective	Function test of software switch (SWS)
Test procedure	
	TEST 1
Test prerequisites	 Machine is ready in automatic mode Folded cartons are available
Required operations	Activate SWS 203 "Blowing air carton"Press "Start"
Consequence	Folding of the side flap is supported by blowing air
Comments	• None
Acknowledgement	Press "Stop"Press "Reset"
	TEST 2
Required operations	Deactivate SWS 203 "Blowing air carton"
Consequence	Folding of the side flap is not supported by blowing air
Comments	• None
Acknowledgement	Press "Stop"Press "Reset"

Acceptance criteria Test 1: Folding of the side flap is supported by blowing air Test 2: Folding of the side flap is not supported by blowing air Comments	s/no
Comments	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWS 208: CONTINUOUS CALL CARTON 9.3.17

Test objective	Function test of software switch (SWS)
Test procedure	
	TEST 1
Test prerequisites	Machine is ready in setup modeCartons are available
Required operations	 Activate SWS 54 "Vacuum" (only with vacuum pump) Activate SWS 208 "Continuous call carton" Press "Start"
Consequence	Folding cartons are drawn off and erected without product with each cycle
Comments	No packaged goods required for the function "Continuous call carton"
Acknowledgement	Press "Stop"Press "Reset"
	TEST 2
Test prerequisites	Machine is ready in setup modeCartons are available
Required operations	Deactivate SWS 208 "Continuous call carton"
Consequence	Cartons are not drawn off and not erected without product with each cycle
Comments	No packaged goods required for the function "Continuous call carton"
Acknowledgement	Press "Stop"Press "Reset"

Test result		yes/no
Acceptance criteria	 Test 1: Folding cartons are drawn off and erected without product with each cycle Test 2: Folding cartons are not drawn off and not erected without product with each cycle 	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWS 256: CHECK CARTON OPEN COVER FLAP 9.3.18

Test objective	Function test of software switch (SWS)
Test procedure	
	TEST 1
Test prerequisites	 Machine is ready in automatic mode Product is available SWS 256 "Check carton open cover flap" is activated Counter 256 "Carton: Consecutive fault open flap" is set to 3
Required operations	 Press "Start" Press "Stop" until some cartons has been closed Put out one cover flap of the carton before the cartons reach the sensors "CAR1.B86-B01"/ "CAR1.B86-B02" Press "Start"
Consequence	 Carton with open cover flap is detected Carton with open cover flap is carried out as bad
Acknowledgement	• None
Comments	• None
	TEST 2
Required operations	 Deactivate SWS 256 "Check carton open cover flap" Press "Start" Press "Stop" until some cartons has been closed Put out one cover flap of the carton before the cartons reach the sensors "CAR1.B86-B01"/ "CAR1.B86-B02" Press "Start"
Consequence	 Carton with open cover flap is not detected Carton with open cover flap is carried out as good
Acknowledgement	• None
Comments	• None

QD 2 Alarm and Fur	nction Testing Execution	125/317
Test result		yes/no
Acceptance criteria	 Test 1: If SWS 256 "Check carton open cover flap" is activated: Carton with open cover flap is detected Test 2: If SWS 256 "Check carton open cover flap" is deactivated: Carton with open cover flap is not detected 	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWS 257: CHECK CARTON OPEN SIDE FLAP 9.3.19

Test objective	Function test of software switch (SWS)
Test procedure	
	TEST 1
Test prerequisites	 Machine is ready in automatic mode Product is available SWS 257 "Check carton open side flap" is activated Counter 257 "Carton: Consecutive fault open flap" is set to 3
Required operations	 Press "Start" Press "Stop" until some cartons has been closed Put out one side flap of the carton before the cartons reach the sensors "CAR1.B86-B03"/ "CAR1.B86-B04" Press "Start"
Consequence	 Carton with open side flap is detected Carton with open side flap is carried out as bad
Acknowledgement	• None
Comments	• None
	TEST 2
Required operations	 Deactivate SWS 257 "Check carton open side flap" Press "Start" Press "Stop" until some cartons has been closed Put out one side flap of the carton before the cartons reach the sensors "CAR1.B86-B03"/ "CAR1.B86-B04" Press "Start"
Consequence	 Carton with open side flap is not detected Carton with open side flap is carried out as good
Acknowledgement	• None
Comments	None

QD 2 Alarm and Fu	nction Testing Execution	127/317
Test result		yes/no
Acceptance criteria	 Test 1: If SWS 257 "Check carton open side flap" is activated: Carton with open side flap is detected Test 2: If SWS 257 "Check carton open side flap" is deactivated: Carton with open side flap is not detected 	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.3.20 SWS 300: LEAFLET DEVICE

Test objective	Function test of software switch (SWS)	
Test procedure		
	TEST 1	
Test prerequisites	Machine is ready in automatic mode with product	
Required operations	Activate SWS 300 "Leaflet device"Press "Start"	
Consequence	Leaflet device is activeLeaflets are called	
	TEST 2	
Required operations	Deactivate SWS 300 "Leaflet device"	
Consequence	Leaflet device is not activeLeaflets are not called	
Comments	• None	
Acknowledgement	Press "Stop"Press "Reset"	
Test result		yes/no
Test result Acceptance criteria	 Test 1: SWS 300 is activated: Leaflet device is active, leaflets are called Test 2: SWS 300 is deactivated: Leaflet device is not active, leaflets are not called 	yes/no
	 Test 2: SWS 300 is deactivated: Leaflet device is not active, 	yes/no
Acceptance criteria	 Test 2: SWS 300 is deactivated: Leaflet device is not active, 	yes/no
Acceptance criteria	 Test 2: SWS 300 is deactivated: Leaflet device is not active, 	yes/no

SWS 302: CODE READER LEAFLET REAR 9.3.21

Test objective	Function test of software switch (SWS)
Test procedure	
	TEST 1
Test prerequisites	 Machine is ready in automatic mode Prepare the leaflet with the manipulated code (e.g. use a black pen to draw an additional line in the code or combine two narrow lines to form a thick line).
Required operations	 Activate SWS 300 "Leaflet device" Activate SWS 302 "Code reader leaflet rear" Insert leaflets with wrong or manipulated code (rear) Press "Start"
Consequence	 All wrong or non-coded leaflets are ejected Counter "Code reader rear" increases by one per defective leaflet
Comments	 Consecutive fault after repeated defective leaflets Leaflet ejection in leaflet unit
Acknowledgement	• None
	TEST 2
Test prerequisites	 Machine is ready in automatic mode Prepare the leaflet with the manipulated code (e.g. use a black pen to draw an additional line in the code or combine two narrow lines to form a thick line).
Required operations	 Activate SWS 300 "Leaflet device" Deactivate SWS 302 "Code reader leaflet rear" Insert leaflets with wrong or manipulated code (rear) Press "Start"
Consequence	All wrong or non-coded leaflets are not ejected
Comments	Leaflet ejection in leaflet unit
Acknowledgement	• None

Test result		yes/no
Acceptance criteria	Test 1: If SWS 302 "Code reader leaflet rear" is activated, all wrong or non-coded leaflets are ejected	
	Test 1: If SWS 302 "Code reader leaflet rear" is activated, counter "Code reader rear" increases by one per defective leaflet	
	Test 2: If SWS 302 "Code reader leaflet rear" is deactivated, leaflets with wrong or missing code on rear side are not ejected	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWS 303: CHECK LEAFLET PRESENT 9.3.22

Test objective	Function test of software switch (SWS)
Test procedure	
	TEST 1
Test prerequisites	Machine is ready in automatic mode
Required operations	 Run the machine with product Activate SWS 300 "Leaflet device" Deactivate SWS 303 "Check Leaflet present" Press "Stop" Open guard and remove leaflet out of carton before leaflet sensor "=CAR1.B75-B04" Press "Reset" Press "Start"
Consequence	Carton without leaflet will not be ejected
Acknowledgement	• None
	TEST 2
Test prerequisites	Machine is ready in automatic mode
Required operations	 Run the machine with product Activate SWS 300 "Leaflet device" Activate SWS 303 "Check Leaflet present" Press "Stop" Open guard and remove leaflet out of carton before leaflet sensor "=CAR1.B75-B04" Press "Reset" Press "Start"
Consequence	 Carton without leaflet will be ejected Counter 315 "Leaflet: Fault cross check in carton" increases by one per defective leaflet
Acknowledgement	Press "Stop"Press "Reset"
Comments	• None

Test result		yes/no
Acceptance criteria	Test 1: SWS 303 "Check Leaflet present" is deactivated: Carton without leaflet will not be ejected	
	Test 2: SWS 303 "Check Leaflet present" is activated: Carton without leaflet will be ejected	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.3.23 SWS 304: INSERT DEFECTIVE LEAFLET

Test objective	Function test of software switch (SWS)
Test procedure	
	TEST 1
Test prerequisites	Machine is ready in automatic mode
Required operations	 Activate SWS 304 "Insert defective leaflet" Insert leaflet with faulty code in stack Insert some product in the infeed Press "Start" Call some products
Consequence	Carton with faulty leaflet is loaded but ejected on discharge belt
Comments	• None
Acknowledgement	Press "Stop"Press "Reset"
	TEST 2
Test prerequisites	Machine is ready in automatic mode
Required operations	 Deactivate SWS 304 "Insert defective leaflet" Insert leaflet with faulty code in stack Insert some product in the infeed Press "Start" Call some products
Consequence	 Leaflet is ejected No carton call Product is ejected at insertion
Comments	• None
Acknowledgement	Press "Stop"Press "Reset"

QD 2 Alaitii aliu Fui	iction resting execution	134/317
Test result		yes/no
Acceptance criteria	 Test 1: SWS 304 is activated: In case of defective leaflet, leaflet will be loaded but rejected on discharge belt Test 2: SWS 304 is deactivated: In case of defective leaflet, leaflet 	
Comments	will be rejected without loading	

Results comply	yes/no	Date/Initials
Results approved	Date/Initials	

SWS 307: CONTINUOUS CALL LEAFLET 9.3.24

Test objective	Function test of software switch (SWS)
Test procedure	
	TEST 1
Test prerequisites	Machine is ready in setup mode
Required operations	 Activate SWS 300 "Leaflet device" Activate SWS 54 "Vacuum" (only with vacuum pump) Activate SWS 307 "Continuous call leaflet" Press "Start"
Consequence	Leaflets are drawn off and folded without product with each cycle
Comments	• None
Acknowledgement	Press "Stop"Press "Reset"
	TEST2
Test prerequisites	Machine is ready in setup mode
Required operations	 Activate SWS 300 "Leaflet device" Activate SWS 54 "Vacuum" (only with vacuum pump) Deactivate SWS 307 "Continuous call leaflet" Press "Start"
Consequence	Leaflets are not drawn off and folded without product with each cycle
Comments	• None
Acknowledgement	Press "Stop"Press "Reset"

yes/no
are
ets

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.3.25 SWS 400: ROBOT 1

Test objective	Function test of software switch (SWS)
Test procedure	
	TEST 1
Test prerequisites	Machine is ready in automatic mode
Required operations	Activate SWS 400 "Robot 1"Press "Start"
Consequence	"Robot 1" is active
	TEST 2
Test prerequisites	Cartoning machine is running and "Robot 1" is active
Required operations	 Press "Stop" Deactivate SWS 400 "Robot 1" Press "Start"
Consequence	"Robot 1" stops and is not active
Comments	• None
Acknowledgement	Press "Reset"
Test result	yes/no
Acceptance criteria	 Test 1: If SWS 400 "Robot 1" is activated, "Robot 1" is active Test 2: If SWS 400 "Robot 1" is deactivated, "Robot 1" stops and is not active
Comments	
Results comply	yes/no Date/Initials

9.3.26 SWS 401: ROBOT 1: RUN EMPTY

Test objective	Function test of software switch (SWS)
Test procedure	
	TEST 1
Test prerequisites	Machine is ready in automatic modeSWS 400 "Robot 1" is activated
Required operations	Activate SWS 401 "Robot 1: Run empty"Press "Start"
Consequence	 Cartoning machine starts operation and product on infeed belt run empty Machine is running with predefined cycles
	TEST 2
Required operations	Deactivate SWS 401 "Robot 1: Run empty"Press "Stop"
Consequence	Cartoning machine and product infeed stop
Comments	• None
Acknowledgement	• None
Test result	yes/no
Test result Acceptance criteria	 Test 1: If SWS 401 "Robot 1: Run empty" is activated, infeed belt runs empty and stops automatically after predefined cycles Test 2: If SWS 401 "Robot 1: Run empty" is deactivated, cartoning
Acceptance criteria	Test 1: If SWS 401 "Robot 1: Run empty" is activated, infeed belt runs empty and stops automatically after predefined cycles
	 Test 1: If SWS 401 "Robot 1: Run empty" is activated, infeed belt runs empty and stops automatically after predefined cycles Test 2: If SWS 401 "Robot 1: Run empty" is deactivated, cartoning
Acceptance criteria	 Test 1: If SWS 401 "Robot 1: Run empty" is activated, infeed belt runs empty and stops automatically after predefined cycles Test 2: If SWS 401 "Robot 1: Run empty" is deactivated, cartoning
Acceptance criteria	 Test 1: If SWS 401 "Robot 1: Run empty" is activated, infeed belt runs empty and stops automatically after predefined cycles Test 2: If SWS 401 "Robot 1: Run empty" is deactivated, cartoning
Acceptance criteria Comments	Test 1: If SWS 401 "Robot 1: Run empty" is activated, infeed belt runs empty and stops automatically after predefined cycles Test 2: If SWS 401 "Robot 1: Run empty" is deactivated, cartoning machine and product infeed stop
Acceptance criteria Comments	Test 1: If SWS 401 "Robot 1: Run empty" is activated, infeed belt runs empty and stops automatically after predefined cycles Test 2: If SWS 401 "Robot 1: Run empty" is deactivated, cartoning machine and product infeed stop

SWS 405: ROBOT 1: CALL PRODUCT 9.3.27

Test objective	Function test of software switch (SWS)
Test procedure	
	TEST 1
Test prerequisites	Machine is ready in automatic modeSWS 400 "Robot 1" is activated
Required operations	 Activate SWS 405 "Robot 1: Call product" Press "Start"
Consequence	Cartoning machine is running and products are called at the infeed
	TEST 2
Required operations	Deactivate SWS 405 "Robot 1: Call product"
Consequence	No products are called
Comments	• None
Acknowledgement	• None
Test result	yes/no
Acceptance criteria	 Test 1: If SWS 405 "Robot 1: Call product" is activated, products are called at the infeed belt Test 2: If SWS 405 "Robot 1: Call product" is deactivated, no products are called at infeed belt
Comments	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.3.28 SWS 500: ROBOT 2

Test objective	Function test of software switch (SWS)
Test procedure	
	TEST 1
Test prerequisites	Machine is ready in automatic mode
Required operations	Activate SWS 500 "Robot 2"Press "Start"
Consequence	Robot 2 is active
	TEST 2
Test prerequisites	Cartoning machine is running and Robot 2 is active
Required operations	 Press "Stop" Deactivate SWS 500 "Robot 2" Press "Start"
Consequence	"Robot 2" stops and is not active
Comments	• None
Acknowledgement	Press "Reset"
Test result	yes/no
Acceptance criteria	 Test 1: If SWS 500 "Robot 2" is activated, "Robot 2" is active Test 2: If SWS 500 "Robot 2" is deactivated, "Robot 2" stops and is not active
Comments	
Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWS 505: ROBOT 2: CALL PRODUCT 9.3.29

Test objective	Function test of software switch (SWS)
Test procedure	
	TEST 1
Test prerequisites	 Machine is ready in automatic mode SWS 500 "Robot 2" is activated
Required operations	Activate SWS 505 "Robot 2: Call product"Press "Start"
Consequence	Cartoning machine is running and products are called at the infeed
	TEST 2
Required operations	Deactivate SWS 505 "Robot 2: Call product"
Consequence	No products are called
Comments	• None
Acknowledgement	• None
Test result	yes/no
Acceptance criteria	 Test 1: If SWS 505 "Robot 2: Call product" is activated, products are called at the product infeed Test 2: If SWS 505 "Robot 2: Call product" is deactivated, no products are called at infeed belt
Comments	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWS 900: BROCHURE INFEED 9.3.30

Test objective	Function test of software switch (SWS)
Test procedure	
	TEST 1
Test prerequisites	Machine is ready in automatic mode
Required operations	Activate SWS 900 "Brochure infeed"Press "Start"
Consequence	Infeed of brochures is running
	TEST 2
Required operations	Deactivate SWS 900 "Brochure infeed"
Consequence	Cartoning line is in operation without brochures infeed
Comments	• None
Acknowledgement	Press "Stop"Press "Reset"
Test result	yes/no
Acceptance criteria	 Test 1: If SWS 900 "Brochure infeed" is activated, infeed of brochures is running Test 2: If SWS 900 "Brochure infeed" is deactivated, cartoning line is in operation without brochures infeed
Comments	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWS 901: BROCHURE INFEED: CODE READER 9.3.31

Test objective •	Function test of software switch (SWS)
Test procedure	
	TEST 1
Test prerequisites	 Machine is ready in automatic mode Prepare brochures with manipulated code (e.g. put a line between the code with a black pen or combine two small lines to one bold line)
Required operations	 Activate SWS 900 "Brochure infeed" Activate SWS 901 "Brochure infeed: Code reader" Insert the brochures with wrong or manipulated code Press "Start"
Consequence	 Manipulated brochures and product are rejected, no carton call Counter 901 "Brochure infeed: Consecutive fault code reading" increases by one per manipulated brochure code
	TEST 2
Test prerequisites	Machine is ready in automatic mode
Required operations	 Deactivate SWS 901 "Brochure infeed: Code reader" Insert brochures with wrong or manipulated code Press "Start"
Consequence	Manipulated brochure is packed together with the product into the carton
Comments	• None
Acknowledgement	Press "Stop"Press "Reset"

		yes/no
	 Test 1: If SWS 901 "Brochure infeed: Code reader" is activated, brochure will be rejected in case of wrong code Test 1: If SWS 901 "Brochure infeed: Code reader" is activated, counter 901 "Brochure infeed: Consecutive fault code reading" increases by one per manipulated brochure code Test 2: If SWS 901 "Inserter: Code reader" is deactivated, insert is inserted into carton together with product 	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWS 902: BROCHURE INFEED: CONTINUOUS CALL 9.3.32

Test objective	Function test of software switch (SWS)
Test procedure	
	TEST 1
Test prerequisites	 Machine is ready in setup mode SWS 900 "Brochure infeed" is activated SWS 902 "Brochure infeed: Continuous call" is activated The magazine belt of the brochure infeed is empty
Required operations	Press "Start"
Consequence	 The brochure infeed can be started without product (brochures) The brochures are called continuously
Acknowledgement	Press "Stop"Press "Reset"
TEST 2	
Test prerequisites	Machine is ready in setup mode
Required operations	 Deactivate SWS 902 "Brochure infeed: Continuous call" Press "Start"
Consequence	 The brochure infeed cannot be started No brochures are called
Comments	• None
Acknowledgement	• None

Test result		yes/no
Acceptance criteria	Test 1: If SWS 902 "Brochure infeed: Continuous call" is activated, brochures are called continuously	
	Test 2: If SWS 902 "Brochure infeed: Continuous call" is deactivated, no brochures are called	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWS 2001: ROBOT 2: RUN WITH FLOW-PACK 9.3.33

Test objective	Function test of software switch (SWS)	
	r direction took or continue officer (evro)	
Test procedure		
	TEST 1	
Test prerequisites	 Machine is ready in automatic mode SWSOPM 63 "Line mode upstream equipment" is activated SWS 400 "Robot 1" is activated SWS 500 "Robot 2" is activated 	
Required operations	Activate SWS 2001 "Robot 2: Run with Flow-Pack"Press "Start"	
Consequence	Machine is running with flow-packs	
	TEST 2	
Test prerequisites	Cartoning machine is running with flow-packs	
Required operations	 Press "Stop" Deactivate SWS 2001 "Robot 2: Run with Flow-Pack" Press "Start" 	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Press "Reset"	
Test result		yes/no
Acceptance criteria	 Test 1: SWS 2001 "Robot 2: Run with Flow-Pack" is activated: Machine is running with flow-packs Test 2: SWS 2001 "Robot 2: Run with Flow-Pack" is deactivated: Fault message is displayed on control panel Test 2: SWS 2001 "Robot 2: Run with Flow-Pack" is deactivated: Machine cannot be started as long as fault is active 	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

SWS 2002: CODE READER GLUED-IN BROCHURE 9.3.34

Test objective	Function test of software switch (SWS)
Test procedure	
	TEST 1
Test prerequisites	 Machine is ready in automatic mode Prepare two (2) booklets with wrong / manipulated code Product is available
Required operations	 Insert booklets with wrong / manipulated code into magazine Activate SWS 2002 "Code reader glued-in brochure" Press "Start"
Consequence	Wrong or non-coded booklets are ejected
	TEST 2
Required operations	 Insert booklets with wrong / manipulated code into magazine Deactivate SWS 2002 "Code reader glued-in brochure" Press "Start"
Consequence	Wrong or non-coded booklets are not recognized and ejected
Comments	• None
Acknowledgement	 Activate SWS 2002 "Code reader glued-in brochure" Press "Stop" Press "Reset"

Test result		yes/no
Acceptance criteria	 Test 1: SWS2002 is activated: wrong or non-coded booklets are ejected Test 2: SWS200 is deactivated: wrong or non-coded booklets are not recognized and ejected 	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.4 Test protocols - Messages (ME)

9.4.1 ME 25: BATCH CONTROL: BATCH IS INTERRUPTED

Test objective	Test whether correct message is displayed on control panel
Test procedure	
Test prerequisites	 Machine is ready in automatic mode A released recipe is active A new batch has been created but not started
Required operations	Start batchPress "Interrupt" on "Batch control"
Consequence	Message is displayed on control panel
Comments	Test can be done together with test of ME 27
Acknowledgement	Press "Interrupt" button and confirm indicated message to activate the batch again
Test result	yes/no
Acceptance criteria	Message is displayed on control panel
Comments	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

ME 26: BATCH CONTROL: BATCH IS NOT ACTIVE 9.4.2

Test objective	Test whether correct message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Release a recipe	
Consequence	Message is displayed on control panel	
Comments	• None	
Acknowledgement	None	
Test result		yes/no
Acceptance criteria	Message is displayed on control panel	
·	- Wessage is displayed on control pariet	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.4.3 ME 27: BATCH CONTROL: BATCH IS ACTIVE

Test objective	Test whether correct message is displayed on control panel	
Test procedure		
Test prerequisites	 Machine is ready in automatic mode A released recipe is active A new batch has been created 	
Required operations	Start batch	
Consequence	Message is displayed on control panel	
Comments	Test can be done together with test of ME 25 and ME 26	
Acknowledgement	Press "finish" on "Batch control" and confirm indicated message	
Test result		yes/no
Acceptance criteria	Message is displayed on control panel	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

ME 40: OPERATING MODE: AUTOMATIC 9.4.4

Test objective	Test whether correct message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in setup mode	
Required operations	 Activate SWSOPM 1 "Automatic" Press "Reset" Press "Start" 	
Consequence	Message is displayed on control panel	
Comments	Test can be done together with test of SWSOPM 1	
Acknowledgement	• None	
Test result	yes/n	0
Acceptance criteria	Message is displayed on control panel	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

ME 41: OPERATING MODE: SETUP 9.4.5

Test objective	Test whether correct message is displayed on control panel
Test procedure	
Test prerequisites	Machine is ready in automatic mode
Required operations	 Activate SWSOPM 2 "Setup" Press "Reset" Press "Start"
Consequence	Message is displayed on control panel
Comments	Test can be done together with test of SWSOPM 2
Acknowledgement	• None
Test result	yes/no
Acceptance criteria	Message is displayed on control panel
Comments	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.4.6 ME 42: OPERATING MODE: INCHING MODE PANEL

Test objective	Test whether correct message is displayed on control panel
Test procedure	
Test prerequisites	Machine is ready in automatic mode
Required operations	 Activate SWSOPM 3 "Inching panel" Close all guard doors Press "Reset" Press "Start" and hold the button "Start"
Consequence	Message is displayed on control panel as long as the button "Start" is pressed"
Comments	Test can be done together with test of SWSOPM 3
Acknowledgement	• None
Test result	yes/no
Acceptance criteria	Message is displayed on control panel
Comments	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

ME 43: OPERATING MODE: INCHING MODE CABLE FRONT 9.4.7

Test objective	Test whether correct message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in inching mode	
Required operations	 Plug in "inching cable front" Activate SWSOPM 4 "Inching cable" Close all guard doors Press "Reset" Press button on "inching cable front" 	
Consequence	Message is displayed on control panel as long as the button on inching cable is pressed"	!
Comments	Test can be done together with test of SWSOPM 4	
Acknowledgement	• None	
Test result	yes/no)
Acceptance criteria	Message is displayed on control panel	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.4.8 ME 47: OPERATING MODE: INCHING MODE SINGLE AXES

Test objective	Test whether correct message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in inching mode	
Required operations	 Activate SWSOPM 8 "Inching single axis" Activate SWSOPM_DRV 1054 "Robot 1: Drill axis" Press "Start" and hold the button 	
Consequence	Message is displayed on control panel	
Comments	• None	
Acknowledgement	• None	
Test result		yes/no
Acceptance criteria	Message is displayed on control panel	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.4.9 ME 49: OPERATING MODE: INCHING MODE INFEED

Test objective	Test whether correct message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	 Activate SWSOPM 11 "Inching infeed" Press "Start" at the operating panel of the infeed system 	
Consequence	Message is displayed on control panel	
Comments	• None	
Acknowledgement	None	
Test result		yes/no
Acceptance criteria	Message is displayed on control panel	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.4.10 ME 160: MACHINE: OPERATIONAL

Test objective	Test whether correct message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeNo fault messages pending	
Required operations	Press "Reset" in case of pending fault messages	
Consequence	Message is displayed on control panel	
Comments	• None	
Acknowledgement	• None	
Test result		yes/no
Acceptance criteria	Message is displayed on control panel	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.4.11 ME 448: ROBOT 1: REFERENCE RUN

Test objective	Test whether correct message is displayed on control panel
Test procedure	
Test prerequisites	Machine is ready in automatic modeRobot arm is not installed
Required operations	 Activate SWSOPM 9 "Infeed: Reference run" Confirm indicated message Press "Start"
Consequence	 Message is displayed on control panel Reference run of robot is performed
Comments	• None
Acknowledgement	• None
Test result	yes/no
Acceptance criteria	Message is displayed on control panel
Comments	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

ME 548: ROBOT 2: REFERENCE RUN 9.4.12

Test objective	Test whether correct message is displayed on control panel
Test procedure	
Test prerequisites	Machine is ready in automatic modeRobot arm is not installed
Required operations	 Activate SWSOPM 9 "Infeed: Reference run" Confirm indicated message Press "Start"
Consequence	 Message is displayed on control panel Reference run of robot is performed
Comments	• None
Acknowledgement	• None
Test result	yes/no
Acceptance criteria	Message is displayed on control panel
Comments	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

Test protocols - Warnings (WA) 9.5

WA 27: CONTROL-PC: M-GUARD ENABLED 9.5.1

Test objective	Test whether correct warning message is displayed on control panel
Test procedure	
Test prerequisites	Machine is ready in automatic mode
Required operations	Turn key at control cabinet to "I" to enable M-Guard
Consequence	Warning is displayed on control panel
Comments	• None
Acknowledgement	Switch off key
Test result	yes/no
Acceptance criteria	Warning is displayed on control panel
Comments	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.5.2 WA 28: CONTROL-PC: REMOTE MAINTENANCE ENABLED

Results comply	yes/no Date/Initials
Results approved	Date/Initials

WA 29: OPERATOR: NO USER LOGGED IN 9.5.3

Test objective	Test whether correct warning message is displayed on control panel	
Test procedure		
Test prerequisites	User is logged in	
Required operations	Log out userPress "Start"	
Consequence	Warning is displayed on control panel	
Comments	Window "Please log in first!" is displayed on control panel	
Acknowledgement	Log in userPress "Reset"	
Test result	yes/	no
Acceptance criteria	Warning is displayed on control panel	
	Machine cannot be started as long as warning is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

WA 45: CONTROL CABINET: UPS NOT READY 9.5.4

Test objective	Test whether correct fault message is displayed on control panel
Test procedure	
Test prerequisites	Machine is ready in automatic mode
Required operations	 Switch off the machine at the main switch Wait until all LED lamps on the displkay of the USV are off completely Remove USV (XUSB in) Switch on the machine Log in a user
Consequence	Warning is displayed on control panel
Comments	 Together with the warning message FLT 45 "CONTROLCABINET: UPS NOT READY" is displayed on control panel In order to acknowledge the fault message follow: menu "system configuration" > submenu "diagnostic" > Button "device diagnostic" > submenu "diagnostic USV" > Press button "Monitoring USV off" The simulation of the warning WA "45" can be done together with fault message FLT "45"
Acknowledgement	 Close the HMI and shut down the IPC Switch off the main switch Reconnect the USV Switch on the main switch Log in a user
Test result	yes/no
Acceptance criteria	Warning is displayed on control panel
Comments	
Results comply	yes/no Date/Initials
Results approved	Date/Initials

WA 220: CARTON: PREWARNING LOW STOCK - IB32_B01 9.5.5

Test objective	Test whether correct warning message is displayed on control panel
Test procedure	
Test prerequisites	Machine is ready in automatic mode
Required operations	Remove folding cartons from the magazine until sensor "=CAR1.B32-B01" is uncovered
Consequence	Warning is displayed on control panel
Comments	• None
Acknowledgement	Fill cartons into magazine
Test result	yes/no
Acceptance criteria	Warning is displayed on control panel
Comments	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.5.6 WA 380: LEAFLET INSERTER TURRET MAGAZINE: PREWARNIN LOW LEVEL

Test objective	Test whether correct warning message is displayed on control panel
Test procedure	
Test prerequisites	Machine is ready in automatic modeSWS 300 "Leaflet device" is activated
Required operations	Disconnect input: "CAR1.X161-KI02:2"
Consequence	Warning is displayed on control panel
Comments	• None
Acknowledgement	Reconnect input: "CAR1.X161-KI02:2"
Test result	yes/no
Acceptance criteria	Warning is displayed on control panel
Comments	

Results comply	yes/no	Date/Initials
Results approved	Date/Initials	

WA 381: LEAFLET INSERTER TURRET MAGAZINE: NOT READY 9.5.7

Test objective	Test whether correct warning message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 300 "Leaflet device" is activated	
Required operations	Disconnect input: "CAR1.X161-KI02:3"	
Consequence	Warning is displayed on control panel	
Comments	• None	
Acknowledgement	Reconnect input: "CAR1.X161-KI02:3"	
Test result	у	es/no
Acceptance criteria	Warning is displayed on control panel	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.5.8 WA 800: LEAFLET INSERTER : MINIMAL SUPPLY

Test objective	Test whether correct warning message is displayed on control panel
Test procedure	
Test prerequisites	Machine is ready in automatic modeLeaflets are available
Required operations	 Activate SWS 300 "Leaflet device" Call up leaflets until these are below the sensor "IX161" for minimum accumulation
Consequence	Warning is displayed on control panel
Comments	• None
Acknowledgement	None
Test result	yes/no
Acceptance criteria	Warning is displayed on control panel
Comments	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.5.9 WA 900: BROCHURE INFEED: MINIMAL SUPPLY

Test objective	Test whether correct warning message is displayed on control panel
Test procedure	
Test prerequisites	Machine is ready in automatic modeBrochures are available
Required operations	 Activate SWS 900 "Brochure infeed" Call up brochures until this are below the sensor for minimum accumulation
Consequence	Warning is displayed on control panel
Comments	• None
Acknowledgement	• None
Test result	yes/no
Acceptance criteria	Warning is displayed on control panel
Comments	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.5.10 WA 980: BROCHURE INFEED: TURRET MAGAZINE: MINIMAL SUPPLY

Test objective	Test whether correct warning message is displayed on control panel
Test procedure	
Test prerequisites	Machine is ready in automatic modeBrochures are available
Required operations	Activate SWS 900 "Brochure infeed"Disconnect input "CAR1.X261-KI02:2"
Consequence	Warning is displayed on control panel
Comments	• None
Acknowledgement	• None
Test result	yes/no
Acceptance criteria	Warning is displayed on control panel
Comments	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.5.11 WA 981: BROCHURE INFEED: TURRET MAGAZINE: NOT IN OPERATION MODE

Test objective	Test whether correct warning message is displayed on control panel
Test procedure	
Test prerequisites	 Machine is ready in automatic mode SWS 900 "Brochure infeed" is activated
Required operations	Disconnect input "CAR1.X261-KI02:3"
Consequence	Warning is displayed on control panel
Comments	• None
Acknowledgement	Reconnect input
Test result	yes/no
Acceptance criteria	Warning is displayed on control panel
Comments	

Results comply	yes/no	Date/Initials
Results approved	Date/Initials	

WA 2001: ROBOT 1: UPLINE MACHINE IS NOT RUNNING 9.5.12

Test objective	Test whether correct warning message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 400 "Robot 1" is activated	
Required operations	Disconnect input: "CAR1.W150-KI10:4"	
Consequence	Warning is displayed on control panel	
Comments	• None	
Acknowledgement	Reconnect input: "CAR1.W150-KI10:4"	
Test result	у	es/no
Acceptance criteria	Warning is displayed on control panel	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

WA 2002: ROBOT 2: UPLINE MACHINE IS NOT RUNNING 9.5.13

Test objective	Test whether correct warning message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 500 "Robot 2" is activated	
Required operations	Disconnect input: "CAR1.W250-KI10:1"	
Consequence	Warning is displayed on control panel	
Comments	• None	
Acknowledgement	Reconnect input: "CAR1.W250-KI10:1"	
Test result	у	es/no
Acceptance criteria	Warning is displayed on control panel	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.5.14 WA 2003: ROBOT 2: DOWNLINE MACHINE IS NOT READY

Test objective	Test whether correct warning message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 500 "Robot 2" is activated	
Required operations	Disconnect input: "CAR1.W250-KI10:16"	
Consequence	Warning is displayed on control panel	
Comments	• None	
Acknowledgement	Reconnect input: "CAR1.W250-KI10:16"	
Test result	уе	s/no
Acceptance criteria	Warning is displayed on control panel	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6 Test protocols - Faults (FLT)

9.6.1 FLT 3: OPERATOR: NORMAL STOP ACTIVATED

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Press "Stop"	
Consequence	Fault message is displayed on control panel	
Comments	• None	
Acknowledgement	Press "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.2 FLT 4: OPERATOR: INCORRECT OPERATING MODE SELECTED

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in setup mode	
Required operations	Activate SWSOPM 4 "Inching cable"	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Activate SWSOPM 2 "Setup"Press "Reset"	
Test result		yes/no
Test result Acceptance criteria	Fault message is displayed on control panel	yes/no
	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	yes/no
		yes/no

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.3 FLT 5: COMPRESSED AIR: MONITORING INLET PRESSURE

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Close the compressed air supply on cartoning machine	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Open the compressed air supply on cartoning machinePress "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 6: ETHERCAT: FAULT 9.6.4

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Disconnect ethernet plugPress "Start"	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	 Observe without fail! Control Cabinet door is open! Additionally other fault messages are displayed on control panel The variable contains further information about ethercat connection 	
Acknowledgement	 Reconnect ethernet plug Switch off the main switch of the machine Switch on the main switch of the machine 	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 7: COMPRESSED AIR: MONITORING SAFE SHUTDOWN 9.6.5

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Disconnect "=CAR1.G95-Q01:S1"	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Reconnect "=CAR1.G95-Q01:S1"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.6 FLT 11: MOTION DRIVES: OVERLOAD POWER SUPPLY - IK00_Q10

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Disconnect input: "=CAR1.K00-KI04:10"	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Reconnect input: "=CAR1.K00-KI04:10"	
Test result		yes/no
Test result Acceptance criteria	Fault message is displayed on control panel	yes/no
	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	yes/no

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.7 FLT 17: OPERATOR: OPERATING MODE CHANGE

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in setup mode	
Required operations	Change over to automatic mode	
Consequence	Machine is ready in automatic mode Fault message is displayed on control panel Machine cannot be started as long as fault is active	
Comments	None	
Acknowledgement	Press "Reset"	
Test result	yes/no	
Acceptance criteria	Fault message is displayed on control panel	\bigcup
	Machine cannot be started as long as fault is active	\bigcup
Comments		
Acknowledgement Test result Acceptance criteria	None Press "Reset" yes/no Fault message is displayed on control panel	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 22: VACUUM: OVERLOAD PUMP 9.6.8

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Switch off protective motor switch: "=CAR1.G05-Q10"	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	None	
Acknowledgement	 Switch on protective motor switch: "=CAR1.G05-Q10" Press "Reset" 	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.9 FLT 23: VACUUM: MONITORING NEGATIVE PRESSURE

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	 Activate SWS 54 "Machine vacuum" Disconnect sensor "=CAR1.G05-B80" Press "Start" 	
Consequence	 Machine stops Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Reconnect sensor: "=CAR1.G05-B80"Press "Reset"	
Test result		yes/no
Acceptance criteria	Machine stops	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.10 FLT 34: MACHINE: REPAIR SWITCH DRIVES SWITCHED OFF - IW00_Q02

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Switch off the repair switch	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Switch on the repair switchPress "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 35: CONTROL CABINET: MONITORING MAIN 9.6.11 SWITCH - IW00_Q01

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Switch off main switch and switch on immediately	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Switch off the main switch of the machineSwitch on the main switch of the machine	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 42: CONTROL CABINET S1: OVERTEMPERATURE - IW00_B13 9.6.12

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Adjust value of temperature controller "=CAR1.W00-B13" less than the act temperature in the control cabinet	ual
Consequence	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments	• None	
Acknowledgement	 Set temperature controller "=CAR1.W00-B13" back to previous value of 45 °C 	
Test result	yes	/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 79: EMERGENCY STOP ACTIVATED: LINE 9.6.13

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Line is running at reduced speed	
Required operations	Press "EMERGENCY STOP" strike button of downstream/ upstream	machine
Consequence	 Line stops Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	 Unlock "EMERGENCY STOP" strike button of downstream/ upstream machine Press "Reset" 	n
Test result		yes/no
Acceptance criteria	Line stops	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT0079_AFT_CUKCUC(100223_1)

FLT 81: EMERGENCY STOP ACTIVATED: OPERATING 9.6.14 PANEL - IF91_S90

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is running with reduced speed	
Required operations	Press "EMERGENCY STOP" strike button	
Consequence	 Machine stops Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Unlock "EMERGENCY STOP" strike buttonPress "Reset"	
Test result		yes/no
Test result Acceptance criteria	Machine stops	yes/no
	 Machine stops Fault message is displayed on control panel 	yes/no
		yes/no

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.15 FLT 86: GUARD OPEN: BOTTOM [7]

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Open guard door: botton	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Close guard door: bottomPress "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.16 FLT 87: GUARD OPEN: CARTON BELT [5]

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Open guard door: front carton belt	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Close guard door: front carton beltPress "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.17 FLT 88: GUARD OPEN: INSERTION REAR [20]

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Open guard door: insertion rear	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Close guard door: insertion rearPress "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 111: PRODUCT SENSING: CONSECUTIVE FAULT - IB10_B01 9.6.18

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	 Machine is ready in automatic mode Counter 112 "Product: Consecutive fault sensing" is set to 3 Product scanner is switched on Product is available 	
Required operations	 Press "Start" and run the machine until enough product is in front of th "=CAR1.B10-B01" Press "Stop" Remove 3x product from product chain successively in front of sensor "=CAR1.B10-B01" 	
Consequence	 Machine stop after detecting a missing product three times Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Press "Reset"	
Test result		yes/no
Acceptance criteria	Machine stops	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		
Results comply	yes/no Date/Initials	
Results comply Results approved	yes/no Date/Initials Date/Initials	

FLT 112: PRODUCT SENSING: SENSOR MONITORING - IB10-B01 9.6.19

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Cover sensor: "=CAR1.B10-B01" permanentlyPress "Start"	
Consequence	 Machine stops Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Uncover sensor: "=CAR1.B10-B01"Press "Reset"	
Test result		yes/no
Acceptance criteria	Machine stops	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 120: INSERTION: SAFETY SENSOR - IB36_B01 9.6.20

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	 Press "Start" and stop the machine after placing the product and inserts into the carton Simulate incomplete insertion (e.g. pull leaflet out of the carton which has jubeen inserted or insert a piece of paper so that it protrudes out of the carton about 2 cm) Press "Start" 	ıst
Consequence	 Machine stops Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Remove cartonPress "Reset"	
Test result	yes/r	10
Acceptance criteria	Machine stops	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		
Results comply	yes/no Date/Initials	\supset
Results approved	Date/Initials	

9.6.21 FLT 122: PRODUCT CHAIN: OVERFILL COLLECTING BOX - IB80_B02

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Open guard door: at insertionCover sensor "=CAR1.B80-B02"	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Uncover sensorPress "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 126: INSERTION: OVERLOAD PREINSERTION 9.6.22

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	 Open guard door: at insertion Simulate overload at sensor "=CAR1.W40-B02" (or disconnect input "=K00-KI12:1") Press "Start" 	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Press "Reset"	
Test result	у	es/no
Acceptance criteria	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	Iviacrime cannot be started as long as fault is active	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.23 FLT 129: INSERTION: MONITORING SENSOR SAFETY SENSOR - IB36_B01

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Disconnect sensor "=CAR1.B36-B01"	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Reconnect sensor: "=CAR1.B36-B01"Press "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.24 FLT 200: MACHINE: JAM AT DISCHARGE

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Cover sensor: "=CAR1.X85-B01"	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Uncover sensor: "=CAR1.X85-B01"Press "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.25 FLT 201: DOWNSTREAM: STOP FROM MACHINE 1

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	 Machine is ready in automatic mode SWSOPM 60 "Line mode Downstream equipment" is activated 	
Required operations	Press "Start" on cartoning machinePress Stop on machine 1	
Consequence	 Cartoning machine stop from machine 1 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Press "Reset" on cartoning machinePress "Reset" on machine 1	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.26 FLT 203: DOWNSTREAM: STOP FROM MACHINE 2

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	 Machine is ready in automatic mode SWSOPM 60 "Line mode Downstream equipment" is activated 	
Required operations	Press "Start" on cartoning machinePress Stop on Laser unit	
Consequence	 Cartoning machine stop from laser unit Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Press "Reset" on cartoning machinePress "Reset" on laser unit	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.27 FLT 210: GUARD OPEN: CARTON MAGAZINE [6]

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Open guard door: carton magazine	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	Test can only be performed at standstill, because guard is locked while running	
Acknowledgement	Close guard door: carton magazinePress "Reset"	
Test result	yes	s/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no	Date/Initials
Results approved	Date/Initials	

9.6.28 FLT 212: GUARD OPEN: DISCHARGE TOP [2]

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Open guard door: discharge top	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Close guard door: discharge topPress "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.29 FLT 213: CARTON: OVERLOAD SIDE FLAP FOLDER - IB37_B01

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	 Open guard door: of cartoner Simulate overload at sensor "=CAR1.B37-B01" (move side flap folder back) 	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	 Place side flap folder back to current position Press "Reset" 	
Test result	yes/no)
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.30 FLT 220: CARTON: STOCK AT THE END - IA30_B02

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Empty carton magazine nearly to sensor "=CAR1.A30-B02"Press "Start"	
Consequence	 Machine stops Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Fill up carton magazinePress "Reset"	
Test result		yes/no
Acceptance criteria	Machine stops	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no	Date/Initials
Results approved	Date/Initials	

FLT 221: MACHINE: COVER RAIL OPEN - IA30-B03 9.6.31

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeGuard doors are closed	
Required operations	Turn on "=CAR1.A30-S20" and lift up the cover rail	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Lower cover railTurn off "=CAR1.A30-S20"Press "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 222: CARTON: CONSECUTIVE FAULT CROSSCHECK - IA30_B01 9.6.32

est whether correct fault message is displayed on control panel flachine is ready in automatic mode roduct is available counter 210 "Carton: consecutive fault cross check" is set to 3 demove ≥ 3x cartons from carton magazine in front of sensor =CAR1.A30-B01" one after the other	
roduct is available counter 210 "Carton: consecutive fault cross check" is set to 3 demove ≥ 3x cartons from carton magazine in front of sensor = CAR1.A30-B01" one after the other	
roduct is available counter 210 "Carton: consecutive fault cross check" is set to 3 demove ≥ 3x cartons from carton magazine in front of sensor = CAR1.A30-B01" one after the other	
=CAR1.A30-B01" one after the other	
ress "Start"	
fter 3 undetected cartons machine stops ault message is displayed on control panel	
lone	
ress "Reset"	
	yes/no
fter 3 undetected cartons machine stops ault message is displayed on control panel	
f	ter 3 undetected cartons machine stops

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT0222_AFT_CUKCUC(3)

9.6.33 FLT 223: CARTON: PICKUP CHECK - IA30_B04

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Cover sensor: "=CAR1.A30-B04"Press "Start"	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Uncover sensor: "=CAR1.A30-B04"Press "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.34 FLT 224: CARTON BELT: OVERLOAD

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Switch off protective motor switch: "=CAR1.W31-Q10"	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	 Switch on protective motor switch: "=CAR1.W31-Q10" Press "Reset" 	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.35 FLT 225: CARTON: CONSECUTIVE FAULT CODE READER - IB75_B01

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	 Machine is ready in automatic mode with product SWS 200 "Code reader carton" is activated Counter 240 "Carton: fault code reading" is set to 3 	
Required operations	Place 3 cartons without or with wrong code in the carton magazinePress "Start"	
Consequence	 Machine stop after detecting a defective bar code three times Fault message is displayed on control panel 	
Comments	• None	
Acknowledgement	Press "Reset"	
Test result		yes/no
Acceptance criteria	Machine stops	
	Fault message is displayed on control panel	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.36 FLT 226: CARTON:\FAULT CODE READER - IB75_B01

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	 Machine is ready in automatic mode with product SWS 200 "Code reader carton" is activated Counter 240 "Carton: fault code reading" is set to 3 	
Required operations	Disconnect input "=CAR1.K00-KI04:12"Press "Start"	
Consequence	Fault message is displayed on control panel	
Comments	• None	
Acknowledgement	Reconnect input "=CAR1.K00-KI04:12"Press "Reset"	
Test result	У	es/no
Acceptance criteria	Fault message is displayed on control panel	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 227: CARTON: SENSOR MONITORING CROSSCHECK - IA30_B01 9.6.37

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	 Cover sensor: "=CAR1-A30-B01" Press "Start" and run the machine more than 1 cycle 	
Consequence	 Machine stops Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Uncover sensor: "=CAR1-A30-B01"Press "Reset"	
Test result		yes/no
Acceptance criteria	Machine stops	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 237: CARTON: SENSOR MONITORING PICKUP CHECK - IA30_B04 9.6.38

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Put 24 V to input "=CAR1.K00-KI08:7"Press "Start"	
Consequence	 Machine stops Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	• None	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine stops	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.39 FLT 243: LASER PRINTING: DEVICE NOT READY

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	 Machine is ready in automatic mode SWSOPM 60 "Line mode Downstream equipment" is activated SWS 201 "Carton laser printing" is activated 	
Required operations	Disconnect input: "CAR1.P78-KI00:7"	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	 Reconnect input: "CAR1.P78-KI00:7" Press "Reset" on cartoning machine Press "Reset" on laser unit 	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no	Date/Initials
Results approved	Date/Initials	

FLT 250: CARTON EJECTION 1: OVERFILL COLLECTING CONTAINER 9.6.40 - IQ82_B03

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Cover sensor: "=CAR1.Q82-B03"	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	None	
Acknowledgement	Uncover sensor: "=CAR1.Q82-B03"Press "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 251: CARTON EJECTION 1: CROSS CHECK - IQ82 B01 9.6.41

	-	
Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Cover sensor: "=CAR1.Q82-B01"	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Uncover sensor: "=CAR1.Q82-B01"Press "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 254: CARTON REJECT 1: SENSOR MONITORING CROSS CHECK 9.6.42 - IQ82_B01

Test objective	Test whether correct warning message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	_
Required operations	 Open guard door: discharge top Bridge guard door discharge top Press "Start" Disconnect sensor "=CAR1.Q82-B01" 	
Consequence	 Machine stops Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	 Reconnect sensor: "=CAR1.Q82-B01" Remove bridge guard door discharge top Close guard door: discharge top Press "Reset" 	
Test result	yes/no	
Acceptance criteria	Machine stops Fault message is displayed on control panel)
Comments	Machine cannot be started as long as fault is active	
		_
Results comply	yes/no Date/Initials)
Results approved	Date/Initials)

FLT 270: CARTON: OPEN FLAP: SENSOR MONITORING COVER 9.6.43 FLAP FRONT - IB86_B01

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	 Machine is ready in automatic mode SWS 256 "Check carton open cover flap" is activated 	
Required operations	Disconnect sensor "CAR1.B86-B01"	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Reconnect sensor: "CAR1.B86-B01"Press "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT0270_AFT_CUC(100223)

9.6.44 FLT 271: CARTON: OPEN FLAP: SENSOR MONITORING COVER FLAP REAR - IB86_B02

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	 Machine is ready in automatic mode SWS 256 "Check carton open cover flap" is activated 	
Required operations	Disconnect sensor "CAR1.B86-B02"	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Reconnect sensor: "CAR1.B86-B02"Press "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 272: CARTON: OPEN FLAP: SENSOR MONITORING SIDE FLAP 9.6.45 FRONT - IB86_B03

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	 Machine is ready in automatic mode SWS 257 "Check carton open side flap" is activated 	
Required operations	Disconnect sensor "CAR1.B86-B03"	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Reconnect sensor: "CAR1.B86-B03"Press "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 273: CARTON: OPEN FLAP: SENSOR MONITORING SIDE FLAP REAR 9.6.46 - IB86_B04

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	 Machine is ready in automatic mode SWS 257 "Check carton open side flap" is activated 	
Required operations	Disconnect sensor "CAR1.B86-B04"	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Reconnect sensor: "CAR1.B86-B04"Press "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT0273_AFT_CUC(100223)

FLT 289: CARTON: OPEN FLAP: CONSECUTIVE FAULT COVER 9.6.47 FLAP -IB86_B01 / IB86_B02

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	 Machine is ready in automatic mode Product is available SWS 256 "Check carton open cover flap" is activated Counter 256 "Carton: Consecutive fault open flap" is set to 1 	
Required operations	 Press "Start" Press "Stop" until some cartons has been closed Put out one cover flap of the carton before the carton reach the sensors "CAR1.B86-B01"/ "CAR1.B86-B02" Press "Start" 	
Consequence	 Machine stop after detecting a open cover flap Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Remove open cartonPress "Reset"	
Test result	yes/no	
Acceptance criteria	Machine stop after detecting a open cover flap	\bigcup
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 294: CARTON: OPEN FLAP: CONSECUTIVE FAULT SIDE 9.6.48 FLAP -IB86_B03 / IB86_B04

Test objective	Test whether correct fault message is displayed on control panel
Test procedure	
Test prerequisites	Machine is ready in automatic mode
	Product is available
	SWS 257 "Check carton open side flap" is activated
	Counter 257 "Carton: Consecutive fault open flap" is set to 3
Required operations	Press "Start"
	Press "Stop" until some cartons has been closed
	 Put out one side flap of the carton before the cartons reach the sensors "CAR1.B86-B03"/ "CAR1.B86-B04"
	Press "Start" and stop the machine before the next carton is checked
	Remove open carton
	Press "Reset"
	 Put out one side flap of the carton before the cartons reach the sensors "CAR1.B86-B03"/ "CAR1.B86-B04"
	Press "Start" and stop the machine before the next carton is checked
	Remove open carton
	Press "Reset"
	 Put out one side flap of the carton before the cartons reach the sensors "CAR1.B86-B03"/ "CAR1.B86-B04"
	Press "Start"
Consequence	Machine stop after detecting a open side flap three times
	Fault message is displayed on control panel
	Machine cannot be started as long as fault is active
Comments	• None
Acknowledgement	Press "Reset"

Test result		yes/no
Acceptance criteria	Machine stop after detecting a open side flap three times	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 301: LEAFLET: CONSECUTIVE FAULT CROSS CHECK 9.6.49 IN PRODUCT CHAIN - IX161_B10

Test objective	Test whether correct fault message is displayed on control panel
Test procedure	
Test prerequisites	 Machine is ready in automatic mode SWS 300 "Leaflet device" is activated Product is available
Required operations	 Activate SWS 300 "Leaflet device" Deactivate SWS 100 "Insert defective product" if necessary Press "Start" with product until leaflet gripper chain is filled up Press "Stop" Remove leaflet before sensor "=CAR1.X161-B10" and after sensor "=CAR1.X161-B11" Press "Start"
Consequence	 Machine stops Fault message is displayed on control panel Machine cannot be started as long as fault is active
Comments	• None
Acknowledgement	Press "Reset"

Test result		yes/no
Acceptance criteria	Machine stops	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.50 FLT 302: LEAFLET: CROSS CHECK EJECTION - IX161_B12

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 300 "Leaflet device" is activated	
Required operations	Cover sensor: "=CAR1.X161-B12"Press "Start"	
Consequence	 Machine stops Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Uncover sensor: "=CAR1.X161-B12"Press "Reset"	
Test result		yes/no
Acceptance criteria	Machine stops	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no	Date/Initials
Results approved	Date/Initials	

FLT 304: LEAFLET: CONSECUTIVE FAULT CROSS CHECK 9.6.51 IN UNIT - IX161_B11

Test objective	Test whether correct fault message is displayed on control panel
Test procedure	
Test prerequisites	 Machine is ready in automatic mode Counter 310 "Leaflet: Consecutive fault cross check in device" is set to 3 SWS 300 "Leaflet device" is activated Product is available
Required operations	 Deactivate SWS 100 "Insert defective product" if necessary Disconnect sensor "=CAR1.X161-B11" Press "Start"
Consequence	 Machine stops Fault message is displayed on control panel Machine cannot be started as long as fault is active
Comments	• None
Acknowledgement	Reconnect sensor: "=CAR1.X161-B11"Press "Reset"

Test result		yes/no
Acceptance criteria	Machine stops	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no	Date/Initials
Results approved	Date/Initials	

FLT 321: LEAFLET: CONSECUTIVE FAULT CODE READER 9.6.52 REAR - IB75_B05

Test objective	Test whether correct fault message is displayed on control panel
Test procedure	
Test prerequisites	 Machine is ready in automatic mode Counter 331 "Leaflet: consecutive fault code reader rear" is set to 3 leaflets with no or wrong code are prepared
Required operations	 Activate SWS 300 "Leaflet device" Deactivate SWS 100 "Insert defective product" Activate SWS 302 "Leaflet: code reader rear" Insert 3 leaflets with no or wrong code in leaflet stack successively Press "Start" with product
Consequence	 Leaflets with wrong / missing code will be ejected After 3 defective leaflets machine stops Fault message is displayed on control panel Machine cannot be started as long as fault is active
Comments	To manipulate the code of the leaflet use a black felt-tip pen (e.g. Edding)
Acknowledgement	Press "Reset"

Test result		yes/no
Acceptance criteria	Leaflets with wrong / missing code will be ejected	
	After 3 defective leaflets machine stops	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 322: LEAFLET: CONSECUTIVE FAULT CROSS CHECK IN 9.6.53 CARTON - IB75_B04

Test objective	Test whether correct fault message is displayed on control panel
Test procedure	
Test prerequisites	 Machine is ready in automatic mode SWS 300 "Leaflet device" is activated Counter 303 "Leaflet: Consecutive fault cross check in carton" is set to 3 Product is available
Required operations	 Press "Start" and run the machine until enough product is in front of the sensor "=CAR1.B75-B04" Press "Stop" Open guard door: Remove 3x leaflet before sensor "=CAR1.B75-B04" Close guard door: Press "Reset" Press "Start" with product
Consequence	 After 3 missing leaflets machine stops Fault message is displayed on control panel Machine cannot be started as long as fault is active
Comments Acknowledgement	 Cartons without leaflets will be ejected Press "Reset"

Test result		yes/no
Acceptance criteria	After 3 missing leaflets machine stops	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no	Date/Initials
Results approved	Date/Initials	

FLT 324: LEAFLET: FAULT CODE READER REAR - IB75_B05 9.6.54

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modePrepare one leaflet with no or wrong code	
Required operations	 Activate SWS 300 "Leaflet device" Deactivate SWS 100 "Insert defective product" Activate SWS 302 "Leaflet: code reader rear" Disconnect input: "=CAR1.X161-KI00:11" Press "Start" with product 	
Consequence	 Machine stops Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	To manipulate the code of the leaflet use a black felt-tip pen (e.g. Edding)	ng)
Acknowledgement	Reconnect input: "=CAR1.X161-KI00:11"Press "Reset"	
Test result	у	es/no
Acceptance criteria	Machine stops	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		
Results comply	yes/no Date/Initials	

FLT 325: LEAFLET: SENSOR MONITORING CROSS CHECK IN 9.6.55 CARTON - IB75-B04

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Cover sensor: "=CAR1.B75-B04"Press "Start"	
Consequence	 Machine stops Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Uncover sensor: "=CAR1.B75-B04"Press "Reset"	
Test result		yes/no
Acceptance criteria	Machine stops	
Acceptance criteria	Machine stopsFault message is displayed on control panel	
Acceptance criteria		
Acceptance criteria Comments	Fault message is displayed on control panel	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.56 FLT 326: LEAFLET: SENSOR MONITORING CROSS CHECK IN PRODUCT CHAIN - IX161_B10

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 300 "Leaflet device" is activated	
Required operations	Cover sensor: "=CAR1.X161-B10"Press "Reset"Press start	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Uncover sensor: "=CAR1.X161-B10"Press "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 408: ROBOT 1: SENSOR MONITORING PRODUCT 9.6.57 SENSING EXCEED HEIGHT

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 400 "Robot 1" is activated	
Required operations	Disconnect sensor "=CAR1-W150-B02"Press "Start"	
Consequence	 Machine stops Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Reconnect sensor: "=CAR1-W150-B02"Press "Reset"	
Test result		yes/no
Acceptance criteria	Machine stops	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		
Results comply	yes/no Date/Initials	
Results approved	Date/Initials	

FLT 412: ROBOT 1: EMERGENCY STOP FRONT ACTIVATED 9.6.58

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	 SWSOPM 60 "Line mode downstream equipment" is activated SWSOPM 62 "Machine: Automatic restart" is activated SWS 400 "Robot 1" is activated Line is running at reduced speed 	
Required operations	Press "EMERGENCY STOP" strike button on Robot 1	
Consequence	 Line stops Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Required operations	Press "Reset"	
Comments	Test can be done together whit test of SWSOPM 62 "Machine: Autorestart"	omatic
Acknowledgement	 Unlock "EMERGENCY STOP" strike button on Robot 1 Press "Reset" 	
Test result		yes/no
Acceptance criteria	Line stops	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		
Results comply	yes/no Date/Initials	
Results approved	Date/Initials	

9.6.59 FLT 421: ROBOT 1: GUARD OPEN: FRONT LEFT [17]

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 400 "Robot 1" is activated	
Required operations	Open guard door: front left	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Close guard door: front leftPress "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.60 FLT 422: ROBOT 1: GUARD OPEN: FRONT RIGHT [16]

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 400 "Robot 1" is activated	
Required operations	Open guard door: front right	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Close guard door: front rightPress "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.61 FLT 430: ROBOT 1: BELT 1: OVERLOAD

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 400 "Robot 1" is activated	
Required operations	Switch off protective motor switch: "=CAR1.W166-Q20"Press "Start"	
Consequence	 Machine stops Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	None	
Acknowledgement	Switch on protective motor switch: "=CAR1.W166-Q20"Press "Reset"	
Test result		yes/no
Acceptance criteria	Machine stops	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.62 FLT 432: ROBOT 1: BELT 2: OVERLOAD

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 400 "Robot 1" is activated	
Required operations	Switch off protective motor switch: "=CAR1.W166-Q30"Press "Start"	
Consequence	 Machine stops Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Switch on protective motor switch: "=CAR1.W166-Q30"Press "Reset"	
Test result		yes/no
Acceptance criteria	Machine stops	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 450: ROBOT 1: PRODUCT SENSING: EXCEED HEIGHT 9.6.63

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 400 "Robot 1" is activated	
Required operations	Cover sensor: "=CAR1-W150-B02"Press "Start"	
Consequence	 Machine stops Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Uncover sensor: "=CAR1-W150-B02"Press "Reset"	
Test result		yes/no
Acceptance criteria	Machine stops	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no	Date/Initials
Results approved	Date/Initials	

FLT 451: ROBOT 1: BELT 1: PRODUCT JAM 9.6.64

Test objective	Test whether correct fault message is displayed on control panel			
Test procedure				
Test prerequisites	Machine is ready in automatic modeSWS 400 "Robot 1" is activated			
Required operations	Cover sensor: "=CAR1-W166-B05"Press "Start"			
Consequence	 Machine stops Fault message is displayed on control panel Machine cannot be started as long as fault is active 			
Comments	• None			
Acknowledgement	Uncover sensor: "=CAR1-W166-B05"Press "Reset"			
Test result		yes/no		
Acceptance criteria	Machine stops			
	Fault message is displayed on control panel			
	Machine cannot be started as long as fault is active			
Comments				

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 457: ROBOT 1: MONITORING COMPRESSED AIR 9.6.65

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 400 "Robot 1" is activated	
Required operations	Close the compressed air supply on robot 1	
Consequence	Fault message is displayed on control panelMachine cannot be started as long as fault is active	
Comments	• None	
Acknowledgement	Open the compressed air supply on robot 1	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 477: ROBOT 1: CONTROL CABINET: OVERTEMPERATURE 9.6.66

Test objective	Test whether correct fault message is displayed on control panel
Test procedure	
Test prerequisites	Machine is ready in automatic modeSWS 400 "Robot 1" is activated
Required operations	Adjust value of temperature controller "=CAR1.W150-B19" less than the actual temperature in the control cabinet
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active
Comments	• None
Acknowledgement	• Set temperature controller "=CAR1.W150-B19" back to previous value of 45°C
Test result	yes/no
Acceptance criteria	Fault message is displayed on control panel
	Machine cannot be started as long as fault is active
Comments	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT0477_AFT_CUC(100223)

FLT 512: ROBOT 2: EMERGENCY STOP FRONT ACTIVATED 9.6.67

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	 SWSOPM 60 "Line mode downstream equipment" is activated SWSOPM 62 "Machine: Automatic restart" is activated SWS 500 "Robot 2" is activated Line is running at reduced speed 	
Required operations	• Press "EMERGENCY STOP" strike button on Robot 2	
Consequence	 Line stops Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Required operations	Press "Reset"	
Comments	Test can be done together whit test of SWSOPM 62 "Machine: Autorestart"	matic
Acknowledgement	 Unlock "EMERGENCY STOP" strike button on Robot 2 Press "Reset" 	
Test result		yes/no
Acceptance criteria	Line stops	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		
Results comply	yes/no Date/Initials	
Results approved	Date/Initials	

9.6.68 FLT 521: ROBOT 2: GUARD OPEN: FRONT LEFT [19]

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	 Machine is ready in automatic mode SWS 500 "Robot 2" is activated 	
Required operations	Open guard door: front left	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Close guard door: front leftPress "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.69 FLT 522: ROBOT 2: GUARD OPEN: FRONT RIGHT [18]

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 500 "Robot 2" is activated	
Required operations	Open guard door: front right	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Close guard door: front rightPress "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no	Date/Initials
Results approved	Date/Initials	

FLT 547: ROBOT 2: EMPTYING 9.6.70

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 500 "Robot 2" is activated	
Required operations	 Lock the cylinder "=CAR1.W266-M2661" in the middle position (betwee position "=CAR1.W266-B50" and "=CAR1.W266-B51") "=CAR1.W266-B05" Press "Reset" Press start 	en end
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	None	
Acknowledgement	Unlock cylinder "=CAR1.W266-M2661"Press "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.71 FLT 577: ROBOT 2: CONTROL CABINET: OVERTEMPERATURE

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 500 "Robot 2" is activated	
Required operations	Adjust value of temperature controller "=CAR1.W250-B19" less than the act temperature in the control cabinet	ual
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Set temperature controller "=CAR1.W250-B19" back to previous value of 45	;°C
Test result	yes/n	10
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 807: LEAFLET INSERTER: NOT IN OPERATIONAL 9.6.72

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	 Machine is ready in automatic mode Product is available SWS 300 "Leaflet device" is activated 	
Required operations	Disconnect input: "=CAR1.X161-KI00:1"	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Reconnect input "=CAR1.X161-KI00:1"Press "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.73 FLT 808: LEAFLET INSERTER: GUARD OPEN: REAR [23]

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 300 "Leaflet device" is activated	
Required operations	Open guard door: rear	
Consequence	Fault message is displayed on control panelMachine cannot be started as long as fault is active	
Comments	• None	
Acknowledgement	Close guard door: rearPress "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 822: LEAFLET INSERTER: EMERGENCY STOP ACTIVATED 9.6.74

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	 SWSOPM 60 "Line mode downstream equipment" is activated SWSOPM 62 "Machine: Automatic restart" is activated SWS 300 "Leaflet device" is activated Line is running at reduced speed 	
Required operations	Press "EMERGENCY STOP" strike button on leaflet unit	
Consequence	 Line stops Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Required operations	Press "Reset"	
Comments	• None	
Acknowledgement	 Unlock "EMERGENCY STOP" strike button on leaflet unit Press "Reset" 	
Test result		yes/no
Acceptance criteria	Line stops	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		
Results comply	yes/no Date/Initials	
Results approved	Date/Initials	

9.6.75 FLT 843: LEAFLET INSERTER: CONTROL CABINET: VOLTAGE MONITORING

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 300 "Leaflet device" is activated	
Required operations	• Switch off fuse "=CAR1.X161-F3090"	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Switch on fuse "=CAR1.X161-F3090"Press "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.76 FLT 902: BROCHURE INFEED: CONSECUTIVE FAULT CODE READER

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure Test prerequisites	 Machine is ready in automatic mode SWS 100 "Product: Insert defective product" is deactivated if nece SWS 900 "Brochure infeed" is activated SWS 901 "Brochure infeed: Code reader" is activated Counter 901 for "Brochure infeed: Consecutive fault code reading" 	·
Required operations	 Generate 3 brochures with a defective ID code with a black felt-tip Edding 400) Put these brochures one after the other into the brochure stack Press "Start" 	pen (e.g.
Consequence	 Machine stop after detecting a defective ID code three times Fault message is displayed on control panel Machine cannot be started as long as fault is active The defective brochures are transferred out 	
Comments	• None	
Acknowledgement	Press "Reset"	
Test result		yes/no
Test result Acceptance criteria	 Machine stop after detecting a defective ID code three times Fault message is displayed on control panel 	yes/no
		yes/no
	Fault message is displayed on control panel	yes/no
Acceptance criteria Comments	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	yes/no
Acceptance criteria	Fault message is displayed on control panel	yes/no

9.6.77 FLT 903: BROCHURE INFEED: CONSECUTIVE FAULT CROSS CHECK

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	 Machine is ready in automatic mode Product is available SWS 100 "Insert defective product" is deactivated (if necessary) SWS 900 "Brochure infeed" is activated Counter 903 "Brochure infeed: Consecutive fault cross check" is set to 	o "1"
Required operations	 Press "Start" with product until brochure gripper chain is filled up Press "Stop" Remove all brochures (in front of sensor "=CAR1.B75-B05") Cover sensor cross-check sensor 93B1 (rontech unit) Press "Reset" Press "Start" 	
Consequence	 If a brochure is missing the machine stops Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Press "Reset"	
Test result		yes/no
Acceptance criteria	If a brochure is missing the machine stops	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.78 FLT 907: BROCHURE INFEED: NOT IN OPERATIONAL

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	 Machine is ready in automatic mode Product is available SWS 900 "Brochure infeed" is activated 	
Required operations	Disconnect input: "=CAR1.X261-KI02:3"	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Reconnect input "=CAR1.X261-KI02:3"Press "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.79 FLT 916: BROCHURE INFEED: GUARD OPEN: PRODUCT CHAIN INFEED [14]

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 900 "Brochure infeed" is activated	
Required operations	Open guard door:, product chain infeed	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Close guard door:, product chain infeedPress "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.80 FLT 920: BROCHURE INFEED: NOT IN CORRECT POSITION [15]

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	 Machine is ready in automatic mode Product is available SWS 900 "Brochure infeed" is activated 	
Required operations	Pull out the brochure infeed from the CUC line	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Press "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 922: BROCHURE INFEED: EMERGENCY STOP ACTIVATED 9.6.81

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure Test prerequisites	 Machine is ready in automatic mode SWSOPM 60 "Line mode downstream equipment" is activated SWSOPM 62 "Machine: Automatic restart" is activated SWS 900 "Brochure infeed" is activated 	
Required operations Consequence	 Line is running at reduced speed Press "EMERGENCY STOP" strike button on brochure unit Line stops Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Required operations	Press "Reset"	
Comments	• None	
Acknowledgement	Unlock "EMERGENCY STOP" strike button on brochure unitPress "Reset"	
Test result		yes/no
Acceptance criteria	 Line stops Fault message is displayed on control panel 	
	Machine cannot be started as long as fault is active	
Comments		
Comments		
Results comply	yes/no Date/Initials	

FLT 943: BROCHURE INFEED: CONTROL CABINET: VOLTAGE 9.6.82 **MONITORING**

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 900 "Brochure infeed" is activated	
Required operations	• Switch off fuse "=CAR1.X261-F3090"	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Switch on fuse "=CAR1.X261-F3090"Press "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT0943_AFT_CUC(100223)

9.6.83 FLT 1212: GUARD OPEN: COLLECTING CONTAINER 1 [3]

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	 Open guard door: Guard ejector (=CAR1.Q82-B802) Close guard door: Collecting box (=CAR1.Q82-B801) 	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	None	
Acknowledgement	 Close guard door: Guard ejector (=CAR1.Q82-B802) Open guard door: Collecting box (=CAR1.Q82-B801) Press "Reset" 	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 1213: CARTON EJECTION 1: OVERFILL GUARD COLLECTING 9.6.84 **CONTAINER**

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	 Machine is ready in automatic mode Counter 214 for "Carton ejection 1: Overfill guard collecting container" is set to 1 	
Required operations	 Create a defective item (e. g. incorrect code), run the machine until the defective item is rejected Close guard door: Guard ejector (=CAR1.Q82-B802) Press "Start" 	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	 Open guard door: Guard ejector (=CAR1.Q82-B802) Press "Reset" 	
Test result	yes	/no
Acceptance criteria	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 1230: DOWNSTREAM: GUARD OPEN: MACHINE 1 9.6.85

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	 Machine is ready in automatic mode SWSOPM 60 "Line mode downstream equipment" is activated 	
Required operations	Open guard door: Machine 1	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Close guard door: Machine 1Press "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 1413: ROBOT 1: REMOVE PRODUCT BEFORE 9.6.86 REFERENCE/ SYNCHRON RUN

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	 Machine is ready in automatic mode SWS 400 "Robot 1" is activated SWS 500 "Robot 2" is activated Product is available 	
Required operations	 Set product data "2010: product width cycle belt" to 50 mm Press "Start" 	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	 Set product data "2010: product width cycle belt" to originale length Press "Reset" 	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.87 FLT 1416: ROBOT 1: EMERGENCY STOP REAR ACTIVATED

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure Test prerequisites	 Machine is ready in automatic mode SWSOPM 60 "Line mode downstream equipment" is activated SWSOPM 62 "Machine: Automatic restart" is activated SWS 400 "Robot 1" is activated Line is running at reduced speed 	
Required operations Consequence	 Press "EMERGENCY STOP" strike button on Robot 1 rear Line stops Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Required operations	Press "Reset"	
Comments	Test can be done together whit test of SWSOPM 62 "Machine: Automatic restart"	3
Acknowledgement	 Unlock "EMERGENCY STOP" strike button on Robot 1 rear Press "Reset" 	
Test result	уе	s/no
Acceptance criteria	 Line stops Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments		
Results comply	yes/no Date/Initials	
Results approved	Date/Initials	

9.6.88 FLT 1417: ROBOT 1: GUARD OPEN: REAR LEFT [25]

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 400 "Robot 1" is activated	
Required operations	Open guard door: Rear left	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Close guard door: Rear leftPress "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.89 FLT 1418: ROBOT 1: GUARD OPEN: REAR RIGHT [24]

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 400 "Robot 1" is activated	
Required operations	Open guard door: Rear right	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Close guard door: Rear rightPress "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 1419: ROBOT 1: RANGE LIMIT REACHED 9.6.90

Test objective	Test whether correct fault message is displayed on control panel
Test procedure	
Test prerequisites	 Machine is ready in automatic mode SWS 400 "Robot 1" is activated
Required operations	 Open guard door: Turn safety switch "=CAR1.W150-S20" (Open stop break robot) to position "1" Move the robot arm outside the operational limit by hand manually Turn safety switch "=CAR1.W150-S20" (Open stop break robot) to position "0" Close the guard doors from robot Press "Reset"
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active
Comments	• None
Acknowledgement	 Move the robot arm to the initial position by hand manually Press "Reset"
Test result	yes/no
Acceptance criteria	 Fault message is displayed on control panel Machine cannot be started as long as fault is active
Comments	
Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.91 FLT 1516: ROBOT 2: EMERGENCY STOP REAR ACTIVATED

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	 Machine is ready in automatic mode SWSOPM 60 "Line mode downstream equipment" is activated SWSOPM 62 "Machine: Automatic restart" is activated SWS 500 "Robot 2" is activated Line is running at reduced speed 	
Required operations	Press "EMERGENCY STOP" strike button on Robot 2 rear	
Consequence	 Line stops Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Required operations	Press "Reset"	
Comments	Test can be done together whit test of SWSOPM 62 "Machine: Auto restart"	matic
Acknowledgement	Unlock "EMERGENCY STOP" strike button on Robot 2 rearPress "Reset"	
Test result		yes/no
Acceptance criteria	Line stopsFault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments	Machine cannot be started as long as fault is active	
Comments	Machine cannot be started as long as fault is active	
Comments Results comply	Machine cannot be started as long as fault is active yes/no Date/Initials	

9.6.92 FLT 1517: ROBOT 2: GUARD OPEN: REAR LEFT [27]

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 500 "Robot 2" is activated	
Required operations	Open guard door: Rear left	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Close guard door: Rear leftPress "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.93 FLT 1518: ROBOT 2: GUARD OPEN: REAR RIGHT [26]

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	 Machine is ready in automatic mode SWS 500 "Robot 2" is activated 	
Required operations	Open guard door: Rear right	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Close guard door: Rear rightPress "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 1519: ROBOT 2: RANGE LIMIT REACHED 9.6.94

Test objective	Test whether correct fault message is displayed on control panel
Test procedure	
Test prerequisites	 Machine is ready in automatic mode SWS 500 "Robot 2" is activated
Required operations	 Open guard door: Turn safety switch "=CAR1.W250-S20" (Open stop break robot) to position "1" Move the robot arm outside the operational limit by hand manually Turn safety switch "=CAR1.W250-S20" (Open stop break robot) to position "0" Close the guard doors from robot Press "Reset"
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active
Comments	• None
Acknowledgement	 Move the robot arm to the initial position by hand manually Press "Reset"
Test result	yes/no
Acceptance criteria	 Fault message is displayed on control panel Machine cannot be started as long as fault is active
Comments	
Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 2001: GUARD OPEN: BACK SIDE BELOW LOADING [21] 9.6.95

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Open guard door: Back side below loading	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Close guard door: Back side below loadingPress "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.96 FLT 2002: GUARD OPEN: BACK SIDE VACUUM PUMP [22]

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Open guard door: Back side vacuum pump	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Close guard door: Back side vacuum pumpPress "Reset"	
Test result	Y	yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 2004: GUARD OPEN: PRODUCT EJECTION COLLECTING BOX [4] 9.6.97

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Open guard door: Product ejection collecting box	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Close guard door: Product ejection collecting boxPress "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.98 FLT 2005: LASER: SENSOR MONITORING PRODUCT TRIGGER - IP78_B11

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	 Machine is ready in automatic mode SWSOPM 60 "Line mode Downstream equipment" is activated SWS 201 "Carton laser printing" is activated 	
Required operations	Disconnect sensor "=CAR1.IP78-B11"	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Reconnect sensor: "=CAR1.IP78-B11"Press "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 2009: LEAFLET INSERTER: POSITION CHECK LEAFLET INSERTER 9.6.99

Test objective	Test whether correct fault message is displayed on control panel
Test procedure	
Test prerequisites	Machine is ready in automatic modeSWS 300 "Leaflet device" is activated
Required operations	 Disconnect sensor "=CAR1.X161-B820" or "=CAR1.X161-B505" on leaflet unit Press "Start"
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active
Comments	• None
Acknowledgement	 Reconnect sensor: "=CAR1.X161-B820" or "=CAR1.X161-B505" on leaflet unit Press "Reset"
Test result	yes/no
Acceptance criteria	Fault message is displayed on control panel
	Machine cannot be started as long as fault is active
Comments	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT2009_AFT_CUC(100223-2)

FLT 2010: LEAFLET INSERTER: GUARD OPEN: BOTTOM LEFT [12] 9.6.100

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 300 "Leaflet device" is activated	
Required operations	Open guard door: Bottom left on leaflet unit	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	 Close guard door: Bottom left on leaflet unit Press "Reset" 	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 2011: LEAFLET INSERTER: GUARD OPEN: BOTTOM RIGHT [11] 9.6.101

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 300 "Leaflet device" is activated	
Required operations	Open guard door:, bottom right on leaflet unit	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Close guard door: bottom right on leaflet unitPress "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.102 FLT 2012: LEAFLET INSERTER: GUARD OPEN: TOP LEFT [9]

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 300 "Leaflet device" is activated	
Required operations	Open guard door: Top left on leaflet unit	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Close guard door: Top left on leaflet unitPress "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 2013: LEAFLET INSERTER: GUARD OPEN: TOP RIGHT [8] 9.6.103

Test whether correct fault message is displayed on control panel	
Machine is ready in automatic modeSWS 300 "Leaflet device" is activated	
Open guard door: Top right on leaflet unit	
 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
• None	
Close guard door: Top right on leaflet unitPress "Reset"	
	yes/no
Fault message is displayed on control panel	
Machine cannot be started as long as fault is active	
	 Machine is ready in automatic mode SWS 300 "Leaflet device" is activated Open guard door: Top right on leaflet unit Fault message is displayed on control panel Machine cannot be started as long as fault is active None Close guard door: Top right on leaflet unit Press "Reset" Fault message is displayed on control panel

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.104 FLT 2014: LEAFLET INSERTER: POSITION CHECKINFEED BELT [10]

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 300 "Leaflet device" is activated	
Required operations	Disconnect sensor "93W5" on leaflet unitPress "Start"	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Reconnect sensor "93W5" on leaflet unitPress "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 2015: LEAFLET INSERTER: EMERGENCY STOP 9.6.105 **ACTIVATED: TURRET MAGAZINE**

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	 Machine is ready in automatic mode SWS 300 "Leaflet device" is activated Line is running at reduced speed 	
Required operations	Press "EMERGENCY STOP" strike button on leaflet unit	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	 Unlock "EMERGENCY STOP" strike button " on leaflet unit Press "Reset" 	
Test result		yes/no
Acceptance criteria	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 2016: BROCHURE INFEED: EMERGENCY STOP 9.6.106 **ACTIVATED: TURRET MAGAZINE**

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	 Machine is ready in automatic mode SWS 900 "Brochure infeed" is activated Line is running at reduced speed 	
Required operations	Press "EMERGENCY STOP" strike button on brochure unit	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	 Unlock "EMERGENCY STOP" strike button " on brochure unit Press "Reset" 	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.107 FLT 2017: BROCHURE INFEED: COVER OPEN: TRANSFER SECTION [13]

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 900 "Brochure infeed" is activated	
Required operations	Open cover transfer section on brochure unit	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Close cover transfer section on brochure unitPress "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.108 FLT 2019: GUARD OPEN: LASER DOOR [1]

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 201 "Carton laser printing" is activated	
Required operations	Open guard door:, Laser unit	
Consequence	Fault message is displayed on control panelMachine cannot be started as long as fault is active	
Comments	• None	
Acknowledgement	Close guard door: Laser unitPress "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 2021: EMERGENCY STOP ACTIVATED: LASER 9.6.109

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	 Machine is ready in automatic mode SWS 201 "Carton laser printing" is activated Line is running at reduced speed 	
Required operations	Press "EMERGENCY STOP" strike button on laser unit	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	 Unlock "EMERGENCY STOP" strike button " on laser unit Press "Reset" 	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 2022: LEAFLET: MONITORING PROPER SIGNAL CODE 9.6.110 READING LEAFLET REAR - IB75_B05

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	 Machine is ready in automatic mode SWS 300 "Leaflet unit" is activated SWS 302 "Code reader leaflet rear" is activated 	
Required operations	Put 24 V to input "=CAR1.X161-KI00:10"Press "Start"	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	 Remove the 24 V from input "=CAR1.X161-KI00:10" Press "Reset" 	
Test result		yes/no
Acceptance criteria	Machine stops	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		
Results comply	yes/no Date/Initials	

Results approved

Date/Initials

9.6.111 FLT 2023: BROCHURE INSERTER: MONITORING PROPER SIGNAL CODE READING BROCHURE - IB75_B06

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 201 "Carton laser printing" is activated	
Required operations	Put 24 V to input "=CAR1.X261-KI00:10"Press "Start"	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Remove the 24 V from input "=CAR1.X261-KI00:10"Press "Reset"	
Test result		yes/no
Acceptance criteria	Machine stops	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments	Machine cannot be started as long as fault is active	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.112 FLT 2024: BROCHURE: MONITORING PROPER SIGNAL CODE READING BROCHURE IN CARTON - IB75_B07

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 2002 "Code reader glued-in brochure" is activated	
Required operations	Put 24 V to input "=CAR1.K00-KI14:3"Press "Start"	
Consequence	Fault message is displayed on control panelMachine cannot be started as long as fault is active	
Comments	• None	
Acknowledgement	Remove the 24 V from input "=CAR1.K00-KI14:3"Press "Reset"	
- , ,		
Test result		yes/no
Acceptance criteria	Machine stops	yes/no
	 Machine stops Fault message is displayed on control panel 	yes/no
	·	yes/no
	Fault message is displayed on control panel	yes/no

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 2025: CODE READER GLUED-IN BROCHURE NOT READY 9.6.113

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 2002 "Code reader glued-in brochure" is activated	
Required operations	Disconnect input: "=CAR1.K00-KI14:4"Press "Start"	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Reconnect input: "=CAR1.K00-KI14:4"Press "Reset"	
Test result		yes/no
Acceptance criteria	Machine stops	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.114 FLT 2026: GLUED IN BROCHURE: CONSECUTIVE FAULT CODE READER - IB75_B07

Test objective	Test whether correct fault message is displayed on control panel
Test procedure	
Test prerequisites	 Machine is ready in automatic mode SWS 100 "Product: Insert defective product" is deactivated if necessary SWS 2002 "Code reader glued-in brochure" is activated Counter 2001 "Glued-in brochure: Consecutive fault code reader" is set to 3
Required operations	 Generate 3 booklets with a defective ID code Put these booklets with defective ID codes into the insert stack Press "Start"
Consequence	 Machine stop after detecting a defective code three times Fault message is displayed on control panel Machine cannot be started as long as fault is active
Comments	The defective booklets are transferred out
Acknowledgement	Press "Reset"Attach booklets with manipulated ID code to this document
Test result	yes/no
Acceptance criteria	Machine stop after detecting a defective ID code three times
	Fault message is displayed on control panel
	Machine cannot be started as long as fault is active
Comments	
Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.115 FLT 2029: LASER: PRODUCT JAM

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 201 "Carton laser printing" is activated	
Required operations	• Cover sensor: "=CAR1.P78-B01"	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Uncover sensor: "=CAR1.P78-B01"Press "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
Acceptance criteria	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.116 FLT 2036: ROBOT 1: CROSS CHECK EJECTOR

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
	SWS 400 "Robot 1" is activated	
Required operations	• Cover sensor: "=CAR1.W166-B03"	
	Cover sensor: "=CAR1.W166-B01"Cover sensor: "=CAR1.W166-B02"	
	Cover sensor: =CAR1.vv166-B02Press "Start"	
Consequence	Machine stops	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments	• None	
Acknowledgement	Uncover sensor: "=CAR1.W166-B03"	
	• Uncover sensor: "=CAR1.W166-B01"	
	Uncover sensor: "=CAR1.W166-B02"	
	Press "Reset"	
Test result		yes/no
Acceptance criteria	Machine stops	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is setive	
	Machine cannot be started as long as fault is active	
Comments		
)
Results comply	yes/no Date/Initials	
Results approved	Date/Initials	

FLT 2037: ROBOT 1: EXCESS FILL COLLECTING BOX 9.6.117 ACCELERATION BELT

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 400 "Robot 1" is activated	
Required operations	Cover sensor: "=CAR1.W166-B04"Press "Start"	
Consequence	 Machine stops Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Uncover sensor: "=CAR1.W166-B04"Press "Reset"	
Test result		yes/no
Acceptance criteria	 Machine stops Fault message is displayed on control panel 	
Comments	Machine cannot be started as long as fault is active	
Results comply	yes/no Date/Initials	
Results approved	Date/Initials	

FLT 2038: ROBOT 1: EXCESS FILL COLLECTING BOX CYCLE BELT 9.6.118

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeSWS 400 "Robot 1" is activated	
Required operations	Cover sensor: "=CAR1.W166-B10"Press "Start"	
Consequence	 Machine stops Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Uncover sensor: "=CAR1.W166-B10"Press "Reset"	
Test result		yes/no
Acceptance criteria	Machine stops	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 2040: ROBOT 2: EMERGENCY STOP ACTIVATED: PUCK 9.6.119 **OUTFEED BELT**

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	 Machine is ready in automatic mode SWSOPM 60 "Line mode downstream equipment" is activated SWS 500 "Robot 2" is activated Line is running at reduced speed 	
Required operations	Press "EMERGENCY STOP" strike button at puck outfeed belt	
Consequence	 Line stops Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Required operations	Press "Reset"	
Comments	• None	
Acknowledgement	 Unlock "EMERGENCY STOP" strike button at puck outfeed belt Press "Reset" 	
Test result		yes/no
Acceptance criteria	Line stops	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		
Results comply	yes/no Date/Initials	
Results comply Results approved	yes/no Date/Initials Date/Initials	

FLT 2041: ROBOT 2: EMERGENCY STOP ACTIVATED: DISKUS BELT 9.6.120

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	 Machine is ready in automatic mode SWSOPM 60 "Line mode downstream equipment" is activated SWS 500 "Robot 2" is activated Line is running at reduced speed 	
Required operations	Press "EMERGENCY STOP" strike button at diskus belt	
Consequence	 Line stops Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Required operations	Press "Reset"	
Comments	• None	
Acknowledgement	Unlock "EMERGENCY STOP" strike button at diskus beltPress "Reset"	
Test result		yes/no
Acceptance criteria	Line stops	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		
Results comply	yes/no Date/Initials	
Results approved	Date/Initials	

FLT 2043: ROBOT 2: PUCK OUTFEED BELT: OVERLOAD 9.6.121

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeActivate SWS 500 "Robot 2"	
Required operations	Switch off protective motor switch: "=CAR1.W266-Q10"	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	 Switch on protective motor switch: "=CAR1.W266-Q10" Press "Reset" 	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.122 FLT 2044: ROBOT 2: PUCK INFEED BELT: OVERLOAD

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeActivate SWS 500 "Robot 2"	
Required operations	Switch off protective motor switch: "=CAR1.W266-Q20"	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	 Switch on protective motor switch: "=CAR1.W266-Q20" Press "Reset" 	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 2045: ROBOT 2: DISKUS BELT: OVERLOAD 9.6.123

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeActivate SWS 500 "Robot 2"	
Required operations	Switch off protective motor switch: "=CAR1.W266-Q30"	
Consequence	Fault message is displayed on control panelMachine cannot be started as long as fault is active	
Comments	• None	
Acknowledgement	Switch on protective motor switch: "=CAR1.W266-Q30"Press "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 2046: ROBOT 2: PUCK NOT EMPTY AT OUTFEED BELT 9.6.124

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic modeActivate SWS 500 "Robot 2"	
Required operations	Cover sensor: "=CAR1.W266-B04"Press "Start"	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Uncover sensor "=CAR1.W266-B04"Press "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.125 FLT 2055: ROBOT 1: MONITORING SENSOR PRODUCT CYCLE BELT -W166_B06

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	Disconnect sensor "=CAR1.W166-B06"	
Consequence	 Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Reconnect sensor: "=CAR1.W166-B06"Press "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 2056: ROBOT 1: MONITORING SENSOR FLOWPACK CYCLE 9.6.126 BELT -W166_B07

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	 Cover sensor: "=CAR1.W166-B07" Press "Reset" Press start 	
Consequence	 Machine starts but stop after some cycles Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Uncover sensor: "=CAR1.W166-B07"Press "Reset"	
Test result		yes/no
Acceptance criteria	 Machine starts but stop after some cycles Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments		
Results comply	yes/no Date/Initials	
Results approved	Date/Initials	

FLT 2057: CARTON REJECT 1: CROSS CHECK CARTON 9.6.127 REJECT IQ82_B10

Test objective	Test whether correct fault message is displayed on control panel
Test procedure	
Test prerequisites	 Machine is ready in automatic mode Product is available SWS 100 "Product: Insert defective product" is deactivated SWS 200 "Code reader carton" is activated Counter 240 for "Carton: Consecutive fault code reader" is set to 3
Required operations	 Put one carton with manipulated ID Code into the carton stack Press "Start" wait until the carton is scanned by sensor "CAR1.B75-B01" Stop the machine when the out-feed belt conveys the carton out of the machine Remove the carton out of the out-feed belt Press "Reset" Press "Start"
Consequence	 If the sensor "CAR1.Q82-B10 detects that the faulty carton is missing, the machine stops Fault message is displayed on control panel Machine cannot be started as long as fault is active
Comments	• None
Acknowledgement	Press "Reset"

Test result		yes/no
Acceptance criteria	If the sensor "CAR1.Q82-B10 detects that the faulty carton is missing, the machine stops	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

9.6.128 FLT 2058: CARTON REJECT 1: SENSOR MONITORING CROSS CHECKCARTON REJECTED - IQ82_B10

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	 Cover sensor: "=CAR1.Q82-B10" Press "Reset" Press start 	
Consequence	 Machine starts and stop after some cycles Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Uncover sensor: "=CAR1.Q82-B10"Press "Reset"	
Test result		yes/no
Acceptance criteria	Machine starts and stop after some cycles	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		
Results comply	yes/no Date/Initials	

Date/Initials

Results approved

FLT 2059: ROBOT 1: MONITORING SENSOR DISTANCE CONTROL 9.6.129 1 INFEED BELT -W166-B01

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	 Disconnect sensor "=CAR1.W166-B01" Press "Reset" Press start 	
Consequence	 Machine starts and stop after some cycles Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Reconnect sensor: "=CAR1.W166-B01"Press "Reset"	
Test result		yes/no
Acceptance criteria	Machine starts and stop after some cycles	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		
Results comply	yes/no Date/Initials	

Results approved

Date/Initials

9.6.130 FLT 2060: ROBOT 1: MONITORING SENSOR DISTANCE CONTROL 2 INFEED BELT -W166-B02

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	 Disconnect sensor "=CAR1.W166-B02" Press "Reset" Press start 	
Consequence	 Machine starts and stop after some cycles Fault message is displayed on control panel Machine cannot be started as long as fault is active 	
Comments	• None	
Acknowledgement	Reconnect sensor: "=CAR1.W166-B02"Press "Reset"	
Test result		yes/no
Acceptance criteria	Machine starts and stop after some cycles	
	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		
Results comply	yes/no Date/Initials	
Results approved	Date/Initials	

9.6.131 FLT 2061: INSERTION: MONITORING SENSOR: OVERLOAD PREINSERTION -W40-B02

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	 Disconnect sensor "=CAR1.W40-B02" Press "Reset" Press "Start" 	
Consequence	Fault message is displayed on control panelMachine cannot be started as long as fault is active	
Comments	• None	
Acknowledgement	Reconnect sensor: "=CAR1.W40-B02"Press "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT 2062: ROBOT 1: MONITORING SENSOR: EXCESS FILL 9.6.132 COLLECTION BOX CYCLE BELT -W166_B10

Test objective	Test whether correct fault message is displayed on control panel	
Test procedure		
Test prerequisites	Machine is ready in automatic mode	
Required operations	 Disconnect sensor "=CAR1.W166-B10" Press "Reset" Press "Start" 	
Consequence	Fault message is displayed on control panelMachine cannot be started as long as fault is active	
Comments	• None	
Acknowledgement	Reconnect sensor: "=CAR1.W166-B10"Press "Reset"	
Test result		yes/no
Acceptance criteria	Fault message is displayed on control panel	
	Machine cannot be started as long as fault is active	
Comments		

Results comply	yes/no Date/Initials
Results approved	Date/Initials

FLT2062_AFT_CUC(100223)

10 Notes

Alarm and Function Testing Report

CUC 2002 Cartoning machine

100223

Glaxowellcome Production France

Alarm and Function Testing Report approval by ROTZINGER PharmaPack GmbH (after execution)

Function	Name	Position/Company	Signature	Date
Approved by				
Approved by				

The AFT was performed in	from	until
--------------------------	------	-------

Summary result The test execution has shown that,	yes/no
the AFT has been completed successfully without deviations.	
the AFT has been performed with deviations. A deviation report has been created and follow up actions have been defined.	

Final conclusion	yes/no
The next qualification step can be started.	
Before starting the next qualification step, the deviations must be closed.	

Alarm and Function Testing Report approval by Glaxowellcome Production (after execution)

Function	Name	Position/Company	Signature	Date
Approved by				
Approved by				
Approved by				

Document history:

4.0	28.Mar.2023	Version for approval Rev. 4.0, software changes after FAT 3	Lorena Fausten
3.0	17.Mar.2023	Version for approval Rev. 3.0, hardware and software changes after FAT 2	Lorena Fausten
2.0	19.Jan.2023	Version for approval Rev. 2.0, hardware and software changes after FAT 1	Lorena Fausten
1.0	05.Dec.2022	Version for approval	Lorena Fausten
Rev.	Date	Description	Author